

# Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Design](#) / **Blended Bifocals**

## Shamir Insights Inc. – Blended Bifocals

															Sphere	Add	Cyl to	MFH	
	GR	BR	GG	S	AME	AMB	E	R	Available in style mirrors										
Blended Bifocal (digital)																			
<b>1.50 Standard Index Blended Bifocal (digital)</b>																			
Shamir Duo	GR	BR							GR				-9.00 to +6.00	+0.75 to +4.00	-6.00	15			
<b>1.53 Trivex® Blended Bifocal (digital)</b>																			
Shamir Duo	GR	BR							GR				-13.25 to +6.00	+0.75 to +4.00	-6.50	15			
<b>1.59 Polycarbonate Blended Bifocal (digital)</b>																			
Shamir Duo	GR	BR							GR				-14.75 to +6.00	+0.75 to +4.00	-7.00	15			
<b>1.60 High Index Blended Bifocal (digital)</b>																			
Shamir Duo	GR	BR							GR				-15.00 to +6.00	+0.75 to +4.00	-7.00	15			
<b>1.67 Super High Index Blended Bifocal (digital)</b>																			
Shamir Duo	GR	BR							GR				-16.75 to +8.00	+0.75 to +4.00	-8.00	15			
<b>1.74 Ultra High Index Blended Bifocal (digital)</b>																			
Shamir Duo	GR	BR											-18.75 to +12.00	+0.75 to +4.00	-9.00	15			

### TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.