
















Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Index](#) / **1.50 Index**

1.50 Standard Index

	Transitions® Gen#S								Transitions® XTRACTIVE™			Transitions® XTRACTIVE™ POLARIZED*			
Finished Single Vision												Sphere	Cyl to		
1.50 Standard Index FSV															
Essilor FSV Spherical	GR	BR								GR		-4.00 to +4.00	-2.00		
Essilor FSV Aspheric	GR											-4.00 to +4.00	-2.00		
Essilor FSV 1.50 Crizal® Sapphire™ HR	GR									GR		-4.00 to +4.00	-2.00		
Essilor FSV 1.50 Crizal® Rock™	GR									GR		-4.00 to +4.00	-2.00		
Essilor FSV 1.50 Crizal® Easy Pro	GR									GR		-4.00 to +4.00	-2.00		

	Transitions® Gen#S								Transitions® XTRACTIVE™			Transitions® XTRACTIVE™ POLARIZED*			
Semi-Finished Single Vision												Sphere	Cyl to		
1.50 Standard Index SFSV															
Essilor Spherical Single Vision	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.50 to +10.00	-4.00	
Essilor Stylistic® Single Vision	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	-6.00	
Eyezen® +0	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	-6.00	
Eyezen® +1	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	-6.00	
Eyezen® +2	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	-6.00	
Eyezen® +3	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	-6.00	
Eyezen® +4	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	-6.00	
Eyezen® Start	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +6.00	-6.00	

	Transitions® Gen#S								Transitions® XTRACTIVE™			Transitions® XTRACTIVE™ POLARIZED*					
Progressives												Sphere	Add	Cyl to	MFH		
1.50 Standard Index																	
Varilux® XR Track™	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	NVB Measurement	
Varilux® XR Design™	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14	
Varilux® X 4D™	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14	
Varilux® X FIT™	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14	
Varilux® X Design™	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14	
Varilux® Physio® W3+ Eyecode	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14	
Varilux® Physio® W3+ FIT	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14	
Varilux® Physio® W3+	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14	
Varilux Comfort® Max Fit	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	14	
Varilux Comfort® Max	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	14	
Varilux Physio DRx™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50 & +4.00	-6.00	17	
Varilux Physio Short DRx™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.00	-6.00	14	
Varilux Comfort DRx™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-6.00	17	
Varilux Comfort Short DRx™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.00	-6.00	14	
Varilux® Liberty™ 3.0	GR	BR											-10.00 to +7.50	+0.75 to +1.50	-6.00	14	
Varilux Comfort®	GR	BR							GR	BR	GG		-10.00 to +7.00	+0.75 to +3.50	-6.00	17	
Varilux Comfort Short™	GR	BR							GR	BR	GG		-10.00 to +7.00	+0.75 to +3.00	-6.00	14	
Varilux® Stylistic™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-6.00	17	











Essilor Ideal Advanced™ FIT	GR	BR	GG	S	AME	AMB	E	R						-10.00 to +7.50	+0.75 to +3.50	-4.00	17
Essilor Ideal Advanced™ Wrap	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-4.00	17
Essilor Ideal® FIT	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	18
Essilor Ideal® Driving	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	13
Essilor Ideal®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Ideal® Short	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Natural® Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.50	-6.00	15
Essilor Accolade™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Brillenz®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Brillenz® Short	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Ovation® Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Ovation®	GR	BR												-10.00 to +6.00	+0.75 to +3.50	-4.00	17
Essilor SmallFit® Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.00	-6.00	14











HOYA	Transitions Gen#S			Transitions XTRACTIVE™			Transitions XTRACTIVE™ POLARIZED®				
Finished Single Vision										Sphere	Cyl to
1.50 Standard Index FSV											
HOYA FSV	GR									-4.00 to +4.00	-2.00








HOYA	Transitions Gen#S			Transitions XTRACTIVE™			Transitions XTRACTIVE™ POLARIZED®				
Semi-Finished Single Vision										Sphere	Cyl to
1.50 Standard Index SFSV											
HOYA SFSV 1.50	GR	BR					GR	BR		-8.00 to +9.00	-4.00
HOYA MySV 150	GR	BR					GR	BR	GR	-8.00 to +9.00	-4.00
HOYA SV iQ 1.50	GR	BR					GR	BR	GR	-8.00 to +9.00	-4.00
Sync III	GR	BR					GR	BR	GR	-10.00 to +6.00	-4.00





HOYA	Transitions Gen#S			Transitions XTRACTIVE™			Transitions XTRACTIVE™ POLARIZED®						
Progressives										Sphere	Add	Cyl to	MFH
1.50 Standard Index													
HOYALUX iD LifeStyle®4 Indoor 1.50	GR!	BR!								-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Outdoor 1.50	GR!	BR!								-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Urban 1.50	GR!	BR!								-8.00 to +6.00	+0.75 to +4.00	-6.00	14
Amp HD3 1.50	GR!	BR!								-8.00 to +6.00	+0.75 to +4.00	-6.00	18
Amp Mini HD3 1.50	GR!	BR!								-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX ARRAY™ 11, 13, 15, 17, VL 1.50	GR	BR					GR	BR	GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL 1.50	GR	BR					GR	BR	GR	-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ 2 11, 13, 15, 17, VL 1.50	GR	BR					GR	BR	GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ 2 Wrap 11, 13, 15, 17, VL 1.50	GR	BR					GR	BR	GR	-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS 1.50	GR	BR					GR	BR	GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS 1.50	GR	BR					GR	BR	GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	14
Summit cd BKS 1.50	GR	BR					GR	BR	GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	14
Summit ecp BKS 1.50	GR	BR					GR	BR	GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	18
Hoya GP Wide 1.50	GR	BR					GR	BR	GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	18

1.50 Standard Index SFSV													
Signetek Processed SV	GR	BR	GG	S	AME	AMB	E	R	GR			-8.00 to +6.00	-6.00
KODAK Digital SV	GR	BR	GG	S	AME	AMB	E	R	GR			-8.00 to +6.00	-6.00
KODAK PowerUp Lens (.40/.66/1.00D)	GR	BR										-8.00 to +6.00	-6.00
Crossbows Custom SV	GR	BR										-8.00 to +6.00	-6.00

	Transitions gen#S			Transitions XTRACTIVE			Transitions XTRACTIVE <small>POLARIZED</small>					
Bifocal Signature™										Sphere	Add	Cyl to
1.50 Standard Index SFST-28												
Signetek™ Processed 1.50 FT-28									GR	-10.00 to +6.00	+1.00 to +3.00	-6.00

	Transitions gen#S			Transitions XTRACTIVE			Transitions XTRACTIVE <small>POLARIZED</small>						
Progressives										Sphere	Add	Cyl to	MFH
1.50 Standard Index													
KODAK Precise® 1.50 SF	GR									-9.00 to +5.00	+1.00 to +3.00	-4.00	17
KODAK Precise® Digital 1.50	GR									-8.00 to +5.00	+1.00 to +3.00	-6.00	17
KODAK Precise® Plus 1.50	GR	BR	GG	S	AME	AMB	E	R	GR	-8.00 to +4.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus 1.50	GR	BR	GG	S	AME	AMB	E	R	GR	-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Precise® PB 1.50	GR	BR							GR	-8.00 to +4.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB 1.50	GR	BR							GR	-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique™ 1.50	GR	BR	GG	S	AME	AMB	E	R	GR	-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique™ HD 1.50	GR	BR								-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ 1.50	GR	BR	GG	S	AME	AMB	E	R	GR	-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.50	GR	BR	GG	S	AME	AMB	E	R	GR	-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Easy 14 1.50	GR	BR	GG	S	AME	AMB	E	R		-8.00 to +4.00	+0.75 to +3.50	-6.00	14
KODAK Easy 18 1.50	GR	BR	GG	S	AME	AMB	E	R		-8.00 to +4.00	+0.75 to +3.50	-6.00	18
DirecTek 1.50	GR	BR								-8.00 to +4.00	+0.75 to +3.50	-6.00	16.5
DirecTek Short 1.50	GR	BR								-8.00 to +4.00	+0.75 to +3.50	-6.00	13.5
Navigator FBS 1.50	GR	BR								-8.00 to +4.00	+0.75 to +3.50	-6.00	18
Navigator Short FBS 1.50	GR	BR								-8.00 to +4.00	+0.75 to +3.50	-6.00	14
Crossbows Custom U 1.50	GR	BR								-8.00 to +4.00	+0.75 to +3.50	-6.00	13
Crossbows Custom V 1.50	GR	BR								-8.00 to +4.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 13 Progressive 1.50	GR	BR								-8.00 to +4.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 17 Progressive 1.50	GR	BR								-8.00 to +4.00	+0.75 to +3.50	-6.00	17
Crossbows Junior	GR	BR								-8.00 to +4.00	+0.75 to +3.50	-6.00	13

	Transitions gen#S			Transitions XTRACTIVE			Transitions XTRACTIVE <small>POLARIZED</small>					
Semi-Finished Single Vision										Sphere	Cyl to	
1.50 Standard Index SFSV												
TRN 1.50 SV (72mm) UC/HC DIA.72									GR BR	Base 0.50, 2.25, 4.25, 5.25, 6.25, 8.25		
TRN 1.50 DS (76mm) UC/HC DIA.72	GR	BR							GR BR	GR	Base 0.50, 8.75, 10.50	
TRN 1.50 DS (76mm) UC/HC DIA.76	GR	BR							GR BR	GR	Base 2.00, 3.50, 4.25, 5.25, 6.25, 7.25	

	Transitions gen#S			Transitions XTRACTIVE			Transitions XTRACTIVE <small>POLARIZED</small>					
Semi-Finished Single Vision										Sphere	Cyl to	
1.50 Standard Index SFSV												
Unity SVx	GR	BR							GR	GR	-12.25 to +6.00	-6.00
Unity SVxtra	GR	BR							GR	GR	-12.25 to +6.00	-6.00

Unity SVxtreme	GR	BR							GR			GR	-6.00 to +6.00	-6.00
----------------	----	----	--	--	--	--	--	--	----	--	--	----	----------------	-------

	Transitions gen#S								Transitions XTRACTIVE™			Transitions XTRACTIVE™ POLARIZED™						
Progressives														Sphere	Add	Cyl to	MFH	
1.50 Standard Index																		
Unity Via II	GR	BR								GR			GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12	
Unity Via Plus II	GR	BR								GR			GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12	
Unity Via Elite II	GR	BR								GR			GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12	
Unity Via Mobile II	GR	BR								GR			GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12	
Unity Via Wrap II	GR	BR								GR			GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12	
Ethos	GR	BR								GR			GR	-9.00 to +4.50	0.75 to +3.50	-4.00	14	
Ethos Plus	GR	BR								GR			GR	-9.00 to +4.50	0.75 to +3.50	-4.00	14	

	Transitions gen#S								Transitions XTRACTIVE™			Transitions XTRACTIVE™ POLARIZED™						
Finished Single Vision														Sphere	Cyl to			
1.50 Standard Index FSV																		
Younger Optics FSV AR Super Hydrophobic	GR	BR												-4.00 to +4.00	-2.00			

	Transitions gen#S								Transitions XTRACTIVE™			Transitions XTRACTIVE™ POLARIZED™						
Semi-Finished Single Vision														Sphere	Cyl to			
1.50 Standard Index SFSV																		
Younger Optics SFSV (76mm)	GR	BR	GG	SA	AME	AMB	E	R		GR	BR			-9.00 to +6.00	-4.00			
Younger Optics SFSV DIGITAL													GR	-10.00 to +8.50	-4.00			

	Transitions gen#S								Transitions XTRACTIVE™			Transitions XTRACTIVE™ POLARIZED™						
Digital														Sphere	Add			
1.50 Standard Index																		
Camber™										GR	BR			-11.25 to +6.25	-4.00			

	Transitions gen#S								Transitions XTRACTIVE™			Transitions XTRACTIVE™ POLARIZED™						
Bifocal Signature™														Sphere	Add	Cyl to		
1.50 Standard Index SFST-28																		
Younger Optics SFFT 28										GR				-11.00 to +6.00	+1.00 to +3.00	-4.00		

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRACTIVE are registered trademarks, and XTRACTIVE Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

Unless indicated otherwise, all registered trademarks and trademarks are the property of Essilor International and/or its subsidiaries in the United States and in other countries.

All other trademarks are the property of their respective owners.