






# Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Index](#) / [1.60 Index](#)

## 1.60 High Index

													Sphere	Cyl to	
	Semi-Finished Single Vision	GR	BR	GG	S	AME	AMB	E	R						
<b>1.60 High Index Single Vision</b>															
Essilor Spherical Single Vision Thin&Lite® 1.60	GR	BR	GG	S	AME	AMB	E	R						-13.00 to +8.00	-4.00
Eyezen® +0 Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	-6.00
Eyezen® +1 Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	-6.00
Eyezen® +2 Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	-6.00
Eyezen® +3 Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	-6.00
Eyezen® +4 Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	-6.00
Eyezen® Start Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	-6.00

													Sphere	Add	Cyl to	MFH	
	Progressives	GR	BR	GG	S	AME	AMB	E	R								
<b>1.60 High Index</b>																	
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.60	GR	BR	GG											-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Thin&Lite® 1.60	GR	BR	GG											-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Thin&Lite® 1.60	GR	BR	GG											-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit Thin&Lite® 1.60	GR	BR	GG											-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max Thin&Lite® 1.60	GR	BR	GG											-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™ Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux Comfort® Thin&Lite® 1.60	GR													-12.00 to +8.00	+0.75 to +3.50	-4.00	17
Varilux Comfort Short™ Thin&Lite® 1.60	GR													-12.00 to +8.00	+0.75 to +3.00	-4.00	14
Essilor Ideal Advanced™ FIT Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Sport Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Driving Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	+0.75 to +3.50	-6.00	13
Essilor Ideal® Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Short Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Essilor Natural® Digital Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Ovation® Digital Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	+0.75 to +3.50	-6.00	17
SmallFit® Digital Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	+0.75 to +3.00	-6.00	14

																
--	---	--	--	--	--	--	--	--	---	--	--	---	--	--	--	--



<b>1.60 SFSV</b>													
Shamir Autograph® III Single Vision	GR	BR						GR				-14.00 to +8.75	-7.00
Shamir Relax™	GR	BR						GR				-14.00 to +8.75	-7.00

Blended Bifocal (digital)												<b>Sphere</b>	<b>Add</b>	<b>Cyl to</b>	<b>MFH</b>
<b>1.60 High Index Blended Bifocal (digital)</b>															
Shamir Duo	GR	BR						GR				-15.00 to +6.00	+0.75 to +4.00	-7.00	15

Progressives												<b>Sphere</b>	<b>Add</b>	<b>Cyl to</b>	<b>MFH</b>
<b>1.60 High Index</b>															
Shamir Autograph Intelligence™	GR	BR						GR				-14.75 to +8.75	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR						GR				-15.00 to +8.00	+0.75 to +3.50	-7.00	15, 18
Shamir Autograph® II+	GR	BR						GR				-15.00 to +6.00	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Golf™ SuperLite™ 1.6	GR	BR						GR				-6.00 to +5.25	+1.00 to +3.50	-7.00	19
Shamir Spectrum+™	GR	BR						GR				-15.00 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
Shamir Element™	GR	BR										-10.00 to +6.00	+0.75 to +3.50	-7.00	16, 19
Shamir GENESIS HD	GR	BR						GR				-6.00 to +6.00	+0.75 to +3.50	-7.00	16
Shamir Computer™	GR	BR						GR				-6.00 to +8.00	+0.75 to +2.00	-7.00	14 Lower, 10 Upper
Shamir Workspace™	GR	BR						GR				-6.00 to +8.00	+0.75 to +2.50	-7.00	14 Lower, 10 Upper
Shamir InTouch™	GR	BR						GR				-15.00 to +6.00	+0.75 to +3.50	-7.00	15, 18
Shamir Attitude III Sport	GR	BR						GR				-6.50 to +6.00	+0.75 to +3.50	-7.00	18
Shamir Attitude III Fashion	GR	BR						GR				-6.50 to +6.00	+0.75 to +3.50	-7.00	15, 18

Semi-Finished Single Vision												<b>Sphere</b>	<b>Cyl to</b>		
<b>1.60 High Index SFSV</b>															
Signetek Processed SV	GR	BR										-10.00 to +7.50	-6.00		
KODAK Digital SV	GR	BR										-10.00 to +7.50	-6.00		
Crossbows Custom SV	GR	BR										-10.00 to +7.50	-6.00		

Progressives												<b>Sphere</b>	<b>Add</b>	<b>Cyl to</b>	<b>MFH</b>
<b>1.60 High Index</b>															
KODAK Unique	GR	BR										-10.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique HD	GR	BR										-10.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ 1.60	GR	BR										-10.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.60	GR	BR										-10.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Custom U 1.60	GR	BR										-10.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Custom V 1.60	GR	BR										-10.00 to +6.00	+0.75 to +3.50	-6.00	13

Semi-Finished Single Vision												<b>Sphere</b>	<b>Cyl to</b>		
<b>1.60 High Index Single Vision</b>															
TOG 1.60 SV (75mm) UC/HC DIA.65	GR	BR						GR	BR			Base 10.00, 11.00			

TOG 1.60 SV (75mm) UC/HC DIA.70	GR	BR							GR	BR			Base 9.00
TOG 1.60 SV (75mm) UC/HC DIA.75	GR	BR							GR	BR			Base 1.00, 2.50, 4.00, 5.00, 6.00, 7.00, 8.00
TOG 1.60 SV (80mm) UC/HC DIA.80									GR	BR			Base 1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.50
TRN 1.60 SV (80mm) UC/HC DIA.76									GR				Base 8.50
TRN 1.60 SV (80mm) UC/HC DIA.80									GR				Base 1.75, 2.75, 3.75, 4.75, 6.00, 6.75

	Transitions <b>gen#S</b>							Transitions <b>XTRACTIVE</b>			Transitions <b>XTRACTIVE</b> POLARIZED			
Semi-Finished Single Vision													<b>Sphere</b>	<b>Cyl to</b>
<b>1.60 High Index</b>														
Unity SVx	GR	BR											-14.00 to +6.00	-6.00
Unity SVxtra	GR	BR											-14.00 to +6.00	-6.00
Unity SVxtreme	GR	BR											-6.50 to +6.00	-6.00

	Transitions <b>gen#S</b>							Transitions <b>XTRACTIVE</b>			Transitions <b>XTRACTIVE</b> POLARIZED					
Progressives													<b>Sphere</b>	<b>Add</b>	<b>Cyl to</b>	<b>MFH</b>
<b>1.60 High Index</b>																
Unity Via II	GR	BR											-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR											-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR											-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR	BR											-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR											-6.50 to +5.25	+0.50 to +4.50	-6.00	12
Ethos	GR	BR											-9.00 to +5.00	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR											-9.00 to +5.00	0.75 to +3.50	-4.00	14

	Transitions <b>gen#S</b>							Transitions <b>XTRACTIVE</b>			Transitions <b>XTRACTIVE</b> POLARIZED			
Semi-Finished Single Vision													<b>Sphere</b>	<b>Cyl to</b>
<b>1.60 High Index Single Vision</b>														
Younger Optics Single Vision 1.60	GR	BR	GG	SA	AME	AMB	E	R					-11.00 to +8.00	-4.00

	Transitions <b>gen#S</b>							Transitions <b>XTRACTIVE</b>			Transitions <b>XTRACTIVE</b> POLARIZED			
Digital													<b>Sphere</b>	<b>Add</b>
<b>1.60 High Index</b>														
Camber™	GR	BR											-11.50 to +7.25	-4.00

**TRANSITIONS® LIGHT INTELLIGENT LENSES™**



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

Unless indicated otherwise, all registered trademarks and trademarks are the property of Essilor International and/or its subsidiaries in the United States and in other countries.

All other trademarks are the property of their respective owners.