



Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Index](#) / **1.67 Index**

1.67 Super High Index

	Transitions Gen [®] S							Transitions XTRACTIVE [®]			Transitions XTRACTIVE POLARIZED [®]		Sphere	Cyl to	
	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR			BR
Semi-Finished Single Vision															
1.67 Super High Index Single Vision															
Essilor Spherical Single Vision Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-4.00
Essilor Aspheric Single Vision Thin&Lite [®] 1.67	GR	BR	GG											-16.00 to +8.00	-4.00
Essilor Eyecode™ SV Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00
Essilor FIT™ SV Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00
Essilor Azio™ SV Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00
Essilor SV 360°™ Spherical Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00
Essilor Stylistic [®] Single Vision Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00
Eyezen [®] +0 Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00
Eyezen [®] +1 Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00
Eyezen [®] +2 Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00
Eyezen [®] +3 Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00
Eyezen [®] +4 Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00
Eyezen [®] Start Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00

	Transitions Gen [®] S							Transitions XTRACTIVE [®]			Transitions XTRACTIVE POLARIZED [®]		Sphere	Add	Cyl to	MFH	
	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR					BR
Progressives																	
1.67 Super High Index																	
Varilux [®] XR Track™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-12.00 to +9.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux [®] XR Design™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux [®] X 4D™ Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux [®] X FIT™ Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux [®] X Design™ Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux [®] Physio [®] W3+ Eyecode Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux [®] Physio [®] W3+ FIT Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux [®] Physio [®] W3+ Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort [®] Max Fit Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Comfort [®] Max Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™ Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Varilux [®] Liberty™ 3.0	GR	BR												-14.00 to +9.00	+0.75 to +1.50	-6.00	14
Varilux Comfort [®] Thin&Lite [®] 1.67	GR	BR												-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Varilux Comfort [®] Short™ Thin&Lite [®] 1.67	GR	BR												-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Varilux [®] Stylistic™ Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-5.00 to +5.00	+0.75 to +3.50	-3.00	17
Essilor Ideal Advanced™ FIT Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R						-14.00 to +9.00	+0.75 to +3.00	-6.00	18





Essilor Ideal Advanced Wrap™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Sport Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Driving Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Essilor Ideal® Short	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Essilor Natural® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-14.00 to +9.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Brillenz® Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Ovation® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Ovation® Thin&Lite® 1.67	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-4.00	17
Essilor SmallFit® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-14.00 to +9.00	+0.75 to +3.00	-6.00	14

HOYA	Transitions Gen#S			Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED					
Finished Single Vision										Sphere	Cyl to	
1.67 Super High Index FSV												
HOYA 1.67 Aspheric FSV	GR									-8.00 to +2.00	-2.00	
HOYA 1.67 Aspheric FSV AR	GR									-8.00 to +2.00	-2.00	





HOYA	Transitions Gen#S			Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED					
Semi-Finished Single Vision										Sphere	Cyl to	
1.67 Super High Index SFSV												
HOYA SFSV 1.67	GR	BR					GR	BR		-10.00 to +8.00	-4.00	
HOYA MySV 1.67	GR	BR					GR	BR	GR	-10.00 to +8.00	-4.00	
HOYA SV IQ 1.67	GR	BR					GR	BR	GR	-10.00 to +8.00	-4.00	
Sync III	GR	BR					GR	BR	GR	-13.00 to +7.50	-4.00	





HOYA	Transitions Gen#S			Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED						
Progressives										Sphere	Add	Cyl to	MFH
1.67 Super High Index													
HOYALUX iD LifeStyle®4 Indoor 1.67	GR!	BR!								-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Outdoor 1.67	GR!	BR!								-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Urban 1.67	GR!	BR!								-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle 3 Indoor 1.67	GR!	BR!								-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle 3 Outdoor 1.67	GR!	BR!								-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle 3 Urban 1.67	GR!	BR!								-13.00 to +8.00	+0.75 to +4.00	-6.00	14
Amp HD3 1.67	GR!	BR!								-13.00 to +8.00	+0.75 to +4.00	-6.00	18
Amp Mini HD3 1.67	GR!	BR!								-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX ARRAY 11, 13,15,17,VL 1.67	GR	BR					GR	BR	GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY Wrap 11,13,15,17,VL 1.67	GR	BR					GR	BR	GR	-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13,15,17,VL 1.67	GR	BR					GR	BR	GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11,13,15,17,VL 1.67	GR	BR					GR	BR	GR	-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21

1.67																	
Amplitude BKS 1.67	GR	BR						GR	BR		GR		-13.00 to +7.50	+0.75 to +4.00	-6.00	19	
Amplitude Mini BKS 1.67	GR	BR						GR	BR		GR		-13.00 to +7.50	+0.75 to +4.00	-6.00	14	
Summit cd BKS 1.67	GR	BR						GR	BR		GR		-13.00 to +7.50	+0.75 to +4.00	-6.00	14	
Summit ecp BKS 1.67	GR	BR						GR	BR		GR		-13.00 to +7.50	+0.75 to +4.00	-6.00	18	
Hoya GP Wide 1.67	GR	BR						GR	BR		GR		-13.00 to +7.50	+0.75 to +4.00	-6.00	18	
Tact 40 BKS 1.67	GR	BR						GR	BR		GR		-13.00 to +7.50	+0.75 to +4.00	-6.00	18	
Tact 60 BKS 1.67	GR	BR						GR	BR		GR		-13.00 to +7.50	+0.75 to +4.00	-6.00	18	




	Transitions Gen#S									Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED				
Semi-Finished Single Vision														Sphere	Cyl to		
1.67 Super High Index SFSV																	
Nikon SeeMax Infinite	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-14.00 to +9.00	-6.00
Nikon Infinite Single Vision	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-14.00 to +9.00	-6.00
Nikon Smart Boost	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-14.00 to +9.00	-6.00

	Transitions Gen#S									Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED						
Progressives														Sphere	Add	Cyl to	MFH		
1.67 Super High Index																			
Nikon SeeMax Ultimate	GR	BR	GG	S	AME	AMB	E			GR	BR	GG		GR		-14.00 to +9.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	BR	GG	S	AME	AMB	E			GR	BR	GG		GR		-14.00 to +9.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR	GG	S	AME	AMB	E			GR	BR	GG		GR		-14.00 to +9.00	+0.50 to +4.00	-6.00	15, 17
Nikon Move	GR	BR	GG	S	AME	AMB	E			GR	BR	GG		GR		-14.00 to +9.00	+0.50 to +4.00	-6.00	15, 17

	Transitions Signature GEN8									Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED				
Semi-Finished Single Vision														Sphere	Cyl to		
1.67 Super High Index SFSV																	
Seiko PL 1.67 Transitions®										GR						-12.00 to +10.00	-5.00





	Transitions Gen#S									Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED				
Semi-Finished Single Vision														Sphere	Cyl to		
1.67 SFSV																	
Shamir Autograph® III Single Vision	GR	BR								GR	BR			GR		-16.00 to +16.00	-8.00
Shamir Relax™	GR	BR								GR	BR			GR		-16.75 to +16.00	-8.00

	Transitions Gen#S									Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED						
Blended Bifocal (digital)														Sphere	Add	Cyl to	MFH		
1.67 Super High Index Blended Bifocal (digital)																			
Shamir Duo	GR	BR								GR						-16.75 to +8.00	+0.75 to +4.00	-8.00	15

	Transitions Gen#S									Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED						
Progressives														Sphere	Add	Cyl to	MFH		
1.67 Super High Index																			
Shamir Autograph Intelligence™	GR	BR								GR				GR		-16.75 to +10.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR								GR				GR		-16.75 to +9.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR								GR				GR		-16.75 to +8.25	+0.75 to +4.00	-8.00	11, 13, 15, 18

Shamir Spectrum+™	GR	BR							GR		GR	-16.75 to +8.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Element™ SuperLite™ 1.67	GR	BR							GR		GR	-12.50 to +8.00	+0.75 to +3.50	-8.00	16, 19
Shamir GENESIS HD	GR	BR							GR		GR	-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir InTouch™	GR	BR							GR		GR	-16.75 to +8.00	+0.75 to +3.50	-8.00	15, 18
Shamir Attitude III Sport	GR	BR							GR		GR	-6.50 to +6.00	+0.75 to +3.50	-8.00	18
Shamir Attitude III Fashion	GR	BR							GR		GR	-6.50 to +6.00	+0.75 to +3.50	-8.00	15, 18

															
Semi-Finished Single Vision														Sphere	Cyl to
1.67 Super High Index SFSV															
Signetek Processed SV	GR	BR	GG	S	AME	AMB	E	R	GR					-11.00 to +7.50	-6.00
KODAK Digital SV	GR	BR	GG	S	AME	AMB	E	R	GR					-11.00 to +7.50	-6.00
KODAK PowerUp Lens (.40/.66/1.00D)	GR	BR												-11.00 to +6.00	-6.00
Crossbows Custom SV	GR	BR												-11.00 to +7.50	-6.00

																	
Progressives														Sphere	Add	Cyl to	MFH
1.67 Super High Index																	
KODAK Precise® PB 1.67	GR	BR	GG						GR					-11.00 to +6.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB 1.67	GR	BR	GG						GR					-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Precise® Plus Poly	GR	BR	GG	S	AME	AMB	E	R	GR					-11.00 to +6.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus Poly	GR	BR	GG	S	AME	AMB	E	R	GR					-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique	GR	BR	GG	S	AME	AMB	E	R	GR					-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique HD	GR	BR												-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ 1.67	GR	BR	GG	S	AME	AMB	E	R	GR					-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.67	GR	BR	GG	S	AME	AMB	E	R	GR					-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Easy 14 Poly	GR	BR	GG	S	AME	AMB	E	R						-11.00 to +6.00	+0.75 to +3.50	-6.00	14
KODAK Easy 18 Poly	GR	BR	GG	S	AME	AMB	E	R						-11.00 to +6.00	+0.75 to +3.50	-6.00	18
DirecTek 1.67	GR	BR												-11.00 to +6.00	+0.75 to +3.50	-6.00	16.5
DirecTek Short 1.67	GR	BR												-11.00 to +6.00	+0.75 to +3.50	-6.00	13.5
Navigator FBS 1.67	GR	BR												-11.00 to +6.00	+0.75 to +3.50	-6.00	18
Navigator Short FBS 1.67	GR	BR												-11.00 to +6.00	+0.75 to +3.50	-6.00	14
Crossbows Custom U 1.67	GR	BR												-11.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Custom V 1.67	GR	BR												-11.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 13 Progressive 1.67	GR	BR												-11.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 17 Progressive 1.67	GR	BR												-11.00 to +6.00	+0.75 to +3.50	-6.00	17
Crossbows Junior 1.67	GR	BR												-11.00 to +6.00	+0.75 to +3.50	-6.00	13
1.74 Ultra High Index																	
KODAK Unique	GR	BR												-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique HD	GR	BR												-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™	GR	BR												-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.74	GR	BR												-13.00 to +9.00	+0.75 to +3.50	-6.00	13

															
Semi-Finished Single Vision														Sphere	Cyl to
1.67 Super High Index Single Vision															
TOG 1.67 SV (65mm) UC/HC DIA.65	GR	BR							GR	BR				Base 10.00, 11.00	
TOG 1.67 SV (70mm) UC/HC DIA.70	GR	BR							GR	BR				Base 9.00	
TOG 1.67 SV (75mm) UC/HC DIA.75	GR	BR							GR	BR				Base 1.00, 2.50, 4.00, 5.00, 6.00, 7.00, 8.00	
TOG 1.67 SV (80mm) UC/HC DIA.80									GR	BR				Base 1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.50	

TRN 1.67 SV (80mm) UC/HC DIA.80											GR	Base 1.75, 2.75, 3.75, 4.75, 6.00
TRN 1.67 SV (80mm) UC/HC DIA.81											GR	Base 6.75, 8.00

	Transitions Gen#S			Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED						
Semi-Finished Single Vision												Sphere	Cyl to
1.67 Super High Index SFSV													
Unity SVx	GR	BR							GR			-16.00 to +8.00	-6.00
Unity SVxtra	GR	BR							GR			-16.00 to +8.00	-6.00
Unity SVxtreme	GR	BR							GR			-6.00 to +6.00	-6.00

	Transitions Gen#S			Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED								
Progressives												Sphere	Add	Cyl to	MFH
1.67 Super High Index															
Unity Via II	GR	BR							GR			-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR							GR			-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR							GR			-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR	BR							GR			-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR							GR			-6.00 to +5.00	+0.50 to +4.50	-6.00	12
Ethos	GR	BR							GR			-12.00 to +8.00	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR							GR			-12.00 to +8.00	0.75 to +3.50	-4.00	14

	Transitions Gen#S			Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED						
Semi-Finished Single Vision												Sphere	Cyl to
1.67 Super High Index Single Vision													
Younger Optics Single Vision 1.67 MR-10™	GR	BR	GG	SA	AME	AMB	E	R	GR	BR		-11.00 to +7.00	-4.00

	Transitions Gen#S			Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED						
Digital												Sphere	Add
1.67 Super High Index													
Camber™	GR	BR							GR	BR		-14.00 to +8.00	-4.00

	Transitions Gen#S			Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED								
Progressives												Sphere	Add	Cyl to	MFH
1.67 Super High Index															
Younger Image® 1.67 MR-10™	GR	BR										-10.00 to +8.50	+1.00 to +3.50	-4.00	18

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited.

©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

Unless indicated otherwise, all registered trademarks and trademarks are the property of Essilor International and/or its subsidiaries in the United States and in other countries.

All other trademarks are the property of their respective owners.