







Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Transitions® Signature® GEN S™](#) / [Signet Armorlite](#)

Signet Armorlite – Transitions® GEN S™

											
Semi-Finished Single Vision										Sphere	Cyl to
1.50 Standard Index											
Signetek Processed SV	GR BR	GG	S	AME	AMB	E	R			-8.00 to +6.00	-6.00
KODAK Digital SV	GR BR	GG	S	AME	AMB	E	R			-8.00 to +6.00	-6.00
KODAK PowerUp Lens (.40/.66/1.00D)	GR BR									-8.00 to +6.00	-6.00
Crossbows Custom SV	GR BR									-8.00 to +6.00	-6.00
1.53 TRIVEX® Single Vision											
Signetek Processed SV	GR BR	GG	S	AME	AMB	E	R			-10.00 to +6.00	-6.00
KODAK Digital SV	GR BR	GG	S	AME	AMB	E	R			-10.00 to +6.00	-6.00
Crossbows Custom SV	GR BR									-10.00 to +6.00	-6.00
1.59 Polycarbonate Single Vision											
Signetek Processed SV	GR BR	GG	S	AME	AMB	E	R			-10.00 to +6.50	-6.00
KODAK Digital SV	GR BR	GG	S	AME	AMB	E	R			-10.00 to +6.50	-6.00
KODAK PowerUp Lens (.40/.66/1.00D)	GR BR									-10.00 to +6.50	-6.00
Crossbows Custom SV	GR BR									-10.00 to +6.50	-6.00
1.60 High Index Single Vision											
Signetek Processed SV	GR BR									-10.00 to +7.50	-6.00
KODAK Digital SV	GR BR									-10.00 to +7.50	-6.00
Crossbows Custom SV	GR BR									-10.00 to +7.50	-6.00
1.67 Super High Index Single Vision											
Signetek Processed SV	GR BR	GG	S	AME	AMB	E	R			-11.00 to +7.50	-6.00
KODAK Digital SV	GR BR	GG	S	AME	AMB	E	R			-11.00 to +7.50	-6.00
KODAK PowerUp Lens (.40/.66/1.00D)	GR BR									-11.00 to +6.00	-6.00
Crossbows Custom SV	GR BR									-11.00 to +7.50	-6.00
1.74 Ultra High Index Single Vision											
Signetek™ Processed SV	GR BR									-13.00 to +9.00	-6.00
KODAK Digital SV	GR BR									-13.00 to +9.00	-6.00

													
Progressives										Sphere	Add	Cyl to	MFH
1.50 Standard Index													
KODAK Precise® 1.50 SF	GR									-9.00 to +5.00	+1.00 to +3.00	-4.00	17
KODAK Precise® Digital 1.50	GR									-8.00 to +5.00	+1.00 to +3.00	-6.00	17
KODAK Precise® Plus 1.50	GR BR	GG	S	AME	AMB	E	R			-8.00 to +4.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus 1.50	GR BR	GG	S	AME	AMB	E	R			-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Precise® PB 1.50	GR BR									-8.00 to +4.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB 1.50	GR BR									-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique™ 1.50	GR BR	GG	S	AME	AMB	E	R			-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique™ HD 1.50	GR BR									-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ 1.50	GR BR	GG	S	AME	AMB	E	R			-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.50	GR BR	GG	S	AME	AMB	E	R			-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Easy 14 1.50	GR BR	GG	S	AME	AMB	E	R			-8.00 to +4.00	+0.75 to +3.50	-6.00	14
KODAK Easy 18 1.50	GR BR	GG	S	AME	AMB	E	R			-8.00 to +4.00	+0.75 to +3.50	-6.00	18
DirecTek 1.50	GR BR									-8.00 to +4.00	+0.75 to +3.50	-6.00	16.5

KODAK Unique DRO™ HD 1.67	GR	BR	GG	S	AME	AMB	E	R	-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Easy 14 Poly	GR	BR	GG	S	AME	AMB	E	R	-11.00 to +6.00	+0.75 to +3.50	-6.00	14
KODAK Easy 18 Poly	GR	BR	GG	S	AME	AMB	E	R	-11.00 to +6.00	+0.75 to +3.50	-6.00	18
DirecTek 1.67	GR	BR							-11.00 to +6.00	+0.75 to +3.50	-6.00	16.5
DirecTek Short 1.67	GR	BR							-11.00 to +6.00	+0.75 to +3.50	-6.00	13.5
Navigator FBS 1.67	GR	BR							-11.00 to +6.00	+0.75 to +3.50	-6.00	18
Navigator Short FBS 1.67	GR	BR							-11.00 to +6.00	+0.75 to +3.50	-6.00	14
Crossbows Custom U 1.67	GR	BR							-11.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Custom V 1.67	GR	BR							-11.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 13 Progressive 1.67	GR	BR							-11.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 17 Progressive 1.67	GR	BR							-11.00 to +6.00	+0.75 to +3.50	-6.00	17
Crossbows Junior 1.67	GR	BR							-11.00 to +6.00	+0.75 to +3.50	-6.00	13
1.74 Ultra High Index												
KODAK Unique	GR	BR							-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique HD	GR	BR							-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™	GR	BR							-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.74	GR	BR							-13.00 to +9.00	+0.75 to +3.50	-6.00	13

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.