







Transitions® Lenses Product Availability Guide




< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Transitions® GEN S™](#) / [Younger Optics](#)

Younger Optics – Transitions® GEN S™

											
Finished Single Vision										Sphere	Cyl to
1.50 Standard Index											
Younger Optics FSV AR Super Hydrophobic	GR BR									-4.00 to +4.00	-2.00
1.59 Polycarbonate Finished Single Vision											
Younger FSV Poly AR Super Hydrophobic	GR BR									-6.00 to +4.00	-2.00

											
Semi-Finished Single Vision										Sphere	Cyl to
1.50 Standard Index											
Essilor Spherical Single Vision	GR BR GG SA AME AMB E R									-9.00 to +6.00	-4.00
1.53 TRIVEX® Single Vision											
Younger Optics Trilogy® Spherical SV	GR BR									-8.00 to +7.50	-4.00
1.59 Polycarbonate Single Vision											
Younger Optics SV	GR BR GG SA AME AMB E R									-9.00 to +6.00	-4.00
1.60 High Index Single Vision											
Younger Optics Single Vision 1.60	GR BR GG SA AME AMB E R									-11.00 to +8.00	-4.00
1.67 Super High Index Single Vision											
Younger Optics Single Vision 1.67 MR-10™	GR BR GG SA AME AMB E R									-11.00 to +7.00	-4.00
1.74 Ultra High Index Single Vision											
Younger Optics Single Vision 1.74	GR BR GG SA AME AMB E R									-18.00 to +12.00	-4.00

											
Digital										Sphere	Cyl to
1.53 TRIVEX®											
Camber™	GR BR									-11.00 to +6.25	-4.00
1.59 Polycarbonate											
Camber™	GR BR									-10.00 to +6.25	-4.00
1.60 High Index											
Camber™	GR BR									-11.50 to +7.25	-4.00
1.67 Super High Index											
Camber™	GR BR									-14.00 to +8.00	-4.00
1.74 Ultra High Index											
Camber™	GR BR									-19.00 to +9.00	-4.00

													
Progressives										Sphere	Add	Cyl to	MFH
1.53 TRIVEX®													
Younger Trilogy® Image®	GR									-8.00 to +8.00	+1.00 to +3.00	-4.00	18
1.59 index polycarbonate													
Younger Image®	GR BR									-9.00 to +7.00	+1.00 to +3.00	-4.00	18
1.67 Super High Index													

TRANSITIONS® LIGHT INTELLIGENT LENSES™**GR**
Gray**BR**
Brown**GG**
Graphite
Green**S**
Sapphire**AME**
Amethyst**AMB**
Amber**E**
Emerald**R**
RubyAvailable
in style
mirrors

Transitions, *Transitions Signature*, and *XTRActive* are registered trademarks, and *XTRActive Polarized*, the *Transitions* logo and *Transitions Light Intelligent Lenses* are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.