

# Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Suppliers](#) / [Hoya Vision Care](#)

<

## Hoya Vision Care — Entire Product Line

HOYA	Transitions® Gen#S								Transitions® XTRACTIVE™			Transitions® XTRACTIVE™ POLARIZED™		Sphere	Cyl to	
	GR	BR	GR!	BR!	GR	BR	GR	BR	GR	BR	GR	BR				
Finished Single Vision																
<b>1.50 Standard Index FSV</b>																
HOYA FSV	GR														-4.00 to +4.00	-2.00
<b>1.53 TRIVEX® Finished Single Vision</b>																
HOYA Phoenix™ FSV	GR!														-6.00 to +2.00	-2.00
<b>1.59 Polycarbonate Finished Single Vision</b>																
HOYA 1.59 Aspheric FSV	GR														-6.00 to +4.00	-2.00
HOYA 1.59 FSV AR	GR!														-6.00 to +4.00	-2.00
<b>1.67 Finished Single Vision</b>																
HOYA 1.67 Aspheric FSV	GR														-8.00 to +2.00	-2.00
HOYA 1.67 Aspheric FSV AR	GR														-8.00 to +2.00	-2.00
Semi-Finished Single Vision																
<b>1.50 Standard Index SFSV</b>																
HOYA SFSV 1.50	GR	BR							GR	BR					-8.00 to +9.00	-4.00
HOYA MySV 150	GR	BR							GR	BR		GR			-8.00 to +9.00	-4.00
HOYA SV iQ 1.50	GR	BR							GR	BR		GR			-8.00 to +9.00	-4.00
Sync III	GR	BR							GR	BR		GR			-10.00 to +6.00	-4.00
<b>1.53 TRIVEX® Single Vision</b>																
HOYA SFSV Phoenix™	GR!	BR!							GR	BR					-13.00 to +9.00	-4.00
HOYA MySV Phoenix™	GR!	BR!							GR	BR		GR			-13.00 to +9.00	-4.00
HOYA SV iQ Phoenix™	GR!	BR!							GR	BR					-13.00 to +9.00	-4.00
Sync III	GR!	BR!													-10.00 to +6.00	-4.00
<b>1.59 Polycarbonate Single Vision</b>																
HOYA SFSV Polycarbonate	GR	BR							GR	BR					-6.00 to +6.00	-4.00
HOYA MySV Polycarbonate	GR	BR							GR	BR		GR			-10.00 to +6.00	-4.00
HOYA SV iQ Polycarbonate	GR	BR							GR	BR		GR			-10.00 to +6.00	-4.00
Sync III	GR	BR							GR	BR		GR			-13.00 to +5.25	-4.00
<b>1.60 High Index Single Vision</b>																
HOYA SFSV 1.60	GR!	BR!													-10.00 to +8.00	-4.00
HOYA SV iQ 1.60	GR!	BR!													-10.00 to +8.00	-4.00
Sync III	GR!	BR!													-10.00 to +6.00	-4.00
<b>1.67 Super High Index Single Vision</b>																
HOYA SFSV 1.67	GR	BR							GR	BR					-10.00 to +8.00	-4.00
HOYA MySV 1.67	GR	BR							GR	BR		GR			-10.00 to +8.00	-4.00
HOYA SV iQ 1.67	GR	BR							GR	BR		GR			-10.00 to +8.00	-4.00
Sync III	GR	BR							GR	BR		GR			-13.00 to +7.50	-4.00
<b>1.74 ultra High Index Single Vision</b>																
HOYA SFSV 1.74	GR!	BR!							GR	BR					-13.00 to +10.00	-4.00
HOYA SV iQ 1.74	GR!	BR!													-13.00 to +10.00	-4.00
Sync III	GR!	BR!													-13.00 to +8.00	-4.00



Progressives											Sphere	Add	Cyl to	MFH	
<b>1.50 Standard Index</b>															
HOYALUX iD LifeStyle®4 Indoor 1.50	GR! BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	14	
HOYALUX iD LifeStyle®4 Outdoor 1.50	GR! BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	14	
HOYALUX iD LifeStyle®4 Urban 1.50	GR! BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	14	
Amp HD3 1.50	GR! BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	18	
Amp Mini HD3 1.50	GR! BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	14	
HOYALUX ARRAY™ 11, 13, 15, 17, VL 1.50	GR BR							GR BR		GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL 1.50	GR BR							GR BR		GR	-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	
HOYALUX ARRAY™ 2 11, 13, 15, 17, VL 1.50	GR BR							GR BR		GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	
HOYALUX ARRAY™ 2 Wrap 11, 13, 15, 17, VL 1.50	GR BR							GR BR		GR	-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	
Amplitude BKS 1.50	GR BR							GR BR		GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	19	
Amplitude Mini BKS 1.50	GR BR							GR BR		GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	14	
Summit cd BKS 1.50	GR BR							GR BR		GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	14	
Summit ecp BKS 1.50	GR BR							GR BR		GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	18	
Hoya GP Wide 1.50	GR BR							GR BR		GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	18	
Tact 40 BKS 1.50	GR BR							GR BR		GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	18	
Tact 60 BKS 1.50	GR BR							GR BR		GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	18	
<b>1.53 TRIVEX®</b>															
HOYALUX iD LifeStyle®4 Indoor Phoenix™	GR! BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	14	
HOYALUX iD LifeStyle®4 Outdoor Phoenix™	GR! BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	14	
HOYALUX iD LifeStyle®4 Urban Phoenix™	GR! BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	14	
HOYALUX iD LifeStyle®3 Indoor Phoenix™	GR! BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	14	
HOYALUX iD LifeStyle®3 Outdoor Phoenix™	GR! BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	14	
HOYALUX iD LifeStyle®3 Urban Phoenix™	GR! BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	14	
Amp HD3 Phoenix™	GR! BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	18	
Amp Mini HD3 Phoenix™	GR! BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	14	
HOYALUX ARRAY™ 11, 13, 15, 17, VL Phoenix™	GR BR							GR BR		GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL Phoenix™	GR BR							GR BR		GR	-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	
HOYALUX ARRAY 2 11, 13, 15, 17, VL Phoenix™	GR BR							GR BR		GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	
HOYALUX ARRAY 2 Wrap 11, 13, 15, 17, VL Phoenix™	GR BR							GR BR		GR	-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	
Amplitude BKS Phoenix™	GR BR							GR BR		GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	19	
Amplitude Mini BKS Phoenix™	GR BR							GR BR		GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	14	
Summit cd BKS Phoenix™	GR BR							GR BR		GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	18	
Summit ecp BKS Phoenix™	GR BR							GR BR		GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	18	
Hoya GP Wide Phoenix™	GR BR							GR BR		GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	18	
Tact 40 BKS Phoenix™	GR BR							GR BR		GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	18	
Tact 60 BKS Phoenix™	GR BR							GR BR		GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	18	
<b>1.59 Polycarbonate</b>															
HOYALUX iD LifeStyle®4 Indoor Poly	GR! BR!										-6.00 to +6.00	+0.75 to +4.00	-6.00	14	
HOYALUX iD LifeStyle®4 Outdoor Poly	GR! BR!										-6.00 to +6.00	+0.75 to +4.00	-6.00	14	
HOYALUX iD LifeStyle®4 Urban Poly	GR! BR!										-6.00 to +6.00	+0.75 to +4.00	-6.00	14	
HOYALUX iD LifeStyle®3 Indoor Poly	GR! BR!										-6.00 to +6.00	+0.75 to +4.00	-6.00	14	
HOYALUX iD LifeStyle®3 Outdoor Poly	GR! BR!										-6.00 to +6.00	+0.75 to +4.00	-6.00	14	
HOYALUX iD LifeStyle®3 Urban Poly	GR! BR!										-6.00 to +6.00	+0.75 to +4.00	-6.00	14	
Amp HD3 Poly	GR! BR!										-6.00 to +6.00	+0.75 to +4.00	-6.00	18	
Amp Mini HD3 Poly	GR! BR!										-6.00 to +6.00	+0.75 to +4.00	-6.00	14	
HOYALUX ARRAY™ 11, 13, 15, 17, VL Poly	GR BR							GR BR		GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL Poly	GR BR							GR BR		GR	-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	
HOYALUX ARRAY 2 11, 13, 15, 17, VL Poly	GR BR							GR BR		GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	

HOYALUX ARRAY 2 Wrap 11, 13, 15, 17, VL Poly	GR	BR						GR	BR			GR	-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS Poly	GR	BR						GR	BR			GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS Poly	GR	BR						GR	BR			GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	14
Summit cd BKS Poly	GR	BR						GR	BR			GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	14
Summit ecp BKS Poly	GR	BR						GR	BR			GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	18
Hoya GP Wide Poly	GR	BR						GR	BR			GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	18
Tact 40 BKS Poly	GR	BR						GR	BR			GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	18
Tact 60 BKS Poly	GR	BR						GR	BR			GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	18
<b>1.60 High Index</b>																
HOYALUX iD LifeStyle®4 Indoor 1.60	GR!	BR!											-11.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Outdoor 1.60	GR!	BR!											-11.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Urban 1.60	GR!	BR!											-11.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Indoor 1.60	GR!	BR!											-11.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Outdoor 1.60	GR!	BR!											-11.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Urban 1.60	GR!	BR!											-11.00 to +8.00	+0.75 to +4.00	-6.00	14
Amp HD3 1.60	GR!	BR!											-11.00 to +8.00	+0.75 to +4.00	-6.00	14
Amp Mini HD3 1.60	GR!	BR!											-11.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX ARRAY 11, 13, 15, 17, VL 1.60	GR	BR											-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY Wrap 11,13,15,17,VL 1.60	GR	BR											-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13, 15, 17, VL 1.60	GR	BR											-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11,13,15,17,VL 1.60	GR	BR											-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS 1.60	GR	BR											-11.00 to +8.00	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS 1.60	GR	BR											-11.00 to +8.00	+0.75 to +4.00	-6.00	14
Summit cd BKS 1.60	GR	BR											-11.00 to +8.00	+0.75 to +4.00	-6.00	14
Summit ecp BKS 1.60	GR	BR											-11.00 to +8.00	+0.75 to +4.00	-6.00	18
Hoya GP Wide 1.60	GR	BR											-11.00 to +8.00	+0.75 to +4.00	-6.00	18
Tact 40 BKS 1.60	GR	BR											-11.00 to +8.00	+0.75 to +4.00	-6.00	18
Tact 60 BKS 1.60	GR	BR											-11.00 to +8.00	+0.75 to +4.00	-6.00	18
<b>1.67 Super High Index</b>																
HOYALUX iD LifeStyle®4 Indoor 1.67	GR!	BR!											-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Outdoor 1.67	GR!	BR!											-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Urban 1.67	GR!	BR!											-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle 3 Indoor 1.67	GR!	BR!											-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle 3 Outdoor 1.67	GR!	BR!											-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle 3 Urban 1.67	GR!	BR!											-13.00 to +8.00	+0.75 to +4.00	-6.00	14
Amp HD3 1.67	GR!	BR!											-13.00 to +8.00	+0.75 to +4.00	-6.00	18
Amp Mini HD3 1.67	GR!	BR!											-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX ARRAY 11, 13,15,17,VL 1.67	GR	BR						GR	BR			GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY Wrap 11,13,15,17,VL 1.67	GR	BR						GR	BR			GR	-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13,15,17,VL 1.67	GR	BR						GR	BR			GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11,13,15,17,VL 1.67	GR	BR						GR	BR			GR	-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS 1.67	GR	BR						GR	BR			GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS 1.67	GR	BR						GR	BR			GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	14
Summit cd BKS 1.67	GR	BR						GR	BR			GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	14
Summit ecp BKS 1.67	GR	BR						GR	BR			GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	18
Hoya GP Wide 1.67	GR	BR						GR	BR			GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	18
Tact 40 BKS 1.67	GR	BR						GR	BR			GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	18
Tact 60 BKS 1.67	GR	BR						GR	BR			GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	18
<b>1.74 ultra High Index</b>																
HOYALUX iD LifeStyle®4 Indoor 1.74	GR!	BR!											-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX iD LifeStyle®4 Outdoor 1.74	GR!	BR!											-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX iD LifeStyle®4 Urban 1.74	GR!	BR!											-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX iD LifeStyle 3 Indoor 1.74	GR!	BR!											-13.00 to +8.00	+0.75 to +3.50	-6.00	14

HOYALUX iD LifeStyle 3 Outdoor 1.74	GR!	BR!												-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX iD LifeStyle 3 Urban 1.74	GR!	BR!												-13.00 to +8.00	+0.75 to +3.50	-6.00	14
Amp HD3 1.74	GR!	BR!												-13.00 to +8.00	+0.75 to +3.50	-6.00	14
Amp Mini HD3 1.74	GR!	BR!												-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX ARRAY 11, 13, 15, 17, VL 1.74	GR	BR												-13.00 to +8.00	+0.75 to +3.50	-6.00	14, 17, 19, 21
HOYALUX ARRAY Wrap 11, 13, 15, 17, VL 1.74	GR	BR												-6.00 to +4.00	+0.75 to +3.50	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13, 15, 17, VL 1.74	GR	BR												-13.00 to +8.00	+0.75 to +3.50	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 WRAP 11, 13, 15, 17, VL 1.74	GR	BR												-6.00 to +4.00	+0.75 to +3.50	-6.00	14, 17, 19, 21
Amplitude BKS	GR	BR												-13.00 to +8.00	+0.75 to +3.50	-6.00	19
Amplitude Mini BKS	GR	BR												-13.00 to +8.00	+0.75 to +3.50	-6.00	14
Summit cd BKS 1.74	GR	BR												-13.00 to +8.00	+0.75 to +3.50	-6.00	14
Summit ecp BKS 1.74	GR	BR												-13.00 to +8.00	+0.75 to +3.50	-6.00	19
Hoya GP Wide 1.74	GR	BR												-13.00 to +8.00	+0.75 to +3.50	-6.00	18
Tact 40 BKS 1.74	GR	BR												-13.00 to +8.00	+0.75 to +3.50	-6.00	18
Tact 60 BKS 1.74	GR	BR												-13.00 to +8.00	+0.75 to +3.50	-6.00	18

TRANSITIONS® LIGHT INTELLIGENT LENSES™



GR  
Gray



BR  
Brown



GG  
Graphite  
Green



S  
Sapphire



AME  
Amethyst



AMB  
Amber



E  
Emerald



R  
Ruby



Available  
in style  
mirrors

Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.