








Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Suppliers](#) / [Nikon Optical USA](#)

Nikon Optical USA – Entire Product Line

 Semi-Finished Single Vision	 Transitions Gen ^S								 Transitions XTRACTIVE			 Transitions XTRACTIVE POLARIZED	Sphere	Cyl to	
	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR			
1.50 Standard Index Single Vision															
Nikon SeeMax Infinite	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	-6.00
Nikon Infinite Single Vision	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	-6.00
Nikon Smart Boost	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	-6.00
1.53 TRIVEX® Single Vision															
Nikon SeeMax Infinite	GR	BR	GG											-10.00 to +7.00	-6.00
Nikon Infinite Single Vision	GR	BR	GG											-10.00 to +7.00	-6.00
Nikon Smart Boost	GR	BR	GG											-10.00 to +7.00	-6.00
1.59 Polycarbonate Single Vision															
Nikon SeeMax Infinite	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	-6.00
Nikon Infinite Single Vision	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	-6.00
Nikon Smart Boost	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	-6.00
1.60 High Index Single Vision															
Nikon SeeMax Infinite	GR	BR												-12.00 to +8.00	-6.00
Nikon Infinite Single Vision	GR	BR												-12.00 to +8.00	-6.00
Nikon Smart Boost	GR	BR												-12.00 to +8.00	-6.00
1.67 Super High Index SFSV															
Nikon SeeMax Infinite	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00
Nikon Infinite Single Vision	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00
Nikon Smart Boost	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00
1.74 Ultra High Index SFSV															
Nikon SeeMax Infinite	GR	BR												-16.00 to +10.00	-6.00
Nikon Infinite Single Vision	GR	BR												-16.00 to +10.00	-6.00
Nikon Smart Boost	GR	BR												-16.00 to +10.00	-6.00

 Progressives	 Transitions Gen ^S								 Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED	Sphere	Add	Cyl to	MFH	
	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR					
1.50 Standard Index																	
Nikon SeeMax Ultimate	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
Nikon Move	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
1.53 TRIVEX®																	
Nikon SeeMax Ultimate	GR	BR	GG											-10.00 to +7.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	BR	GG											-10.00 to +7.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR	GG											-10.00 to +7.00	+0.50 to +4.00	-6.00	15, 17
Nikon Move	GR	BR	GG											-10.00 to +7.00	+0.50 to +4.00	-6.00	15, 17
1.59 Polycarbonate																	
Nikon SeeMax Ultimate	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-11.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-11.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-11.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
Nikon Move	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-11.00 to +8.00	+0.50 to +4.00	-6.00	15, 17

1.60 High Index																	
Nikon SeeMax Ultimate	GR	BR												-12.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	BR												-12.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR												-12.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
Nikon Move	GR	BR												-12.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
1.67 Super High Index																	
Nikon SeeMax Ultimate	GR	BR	GG	S	AME	AMB	E		GR	BR	GG		GR	-14.00 to +9.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	BR	GG	S	AME	AMB	E		GR	BR	GG		GR	-14.00 to +9.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR	GG	S	AME	AMB	E		GR	BR	GG		GR	-14.00 to +9.00	+0.50 to +4.00	-6.00	15, 17
Nikon Move	GR	BR	GG	S	AME	AMB	E		GR	BR	GG		GR	-14.00 to +9.00	+0.50 to +4.00	-6.00	15, 17
1.74 Ultra High Index																	
Nikon SeeMax Ultimate	GR	BR												-16.00 to +10.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	BR												-16.00 to +10.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR												-16.00 to +10.00	+0.50 to +4.00	-6.00	15, 17

TRANSITIONS® LIGHT INTELLIGENT LENSES™



GR
Gray



BR
Brown



GG
Graphite
Green



S
Sapphire



AME
Amethyst



AMB
Amber



E
Emerald



R
Ruby



Available
in style
mirrors

Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.