


Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Suppliers](#) / [Shamir Insights Inc.](#)

Shamir Insights Inc. — Entire Product Line

	Transitions Gen#S									Transitions XTRACTIVE™			Transitions XTRACTIVE™ POLARIZED™			Sphere	Cyl to	
Finished Single Vision																		
1.59 Polycarbonate FSV																		
Shamir FSV with Glacier PLUS™ AR	GR																-6.00 to +4.00	-2.00

	Transitions Gen#S									Transitions XTRACTIVE™			Transitions XTRACTIVE™ POLARIZED™			Sphere	Cyl to	
Semi-Finished Single Vision																		
1.50 SFSV																		
Shamir Autograph® III Single Vision	GR	BR								GR	BR			GR			-12.00 to +7.25	-6.00
Shamir Relax™	GR	BR								GR	BR			GR			-14.00 to +8.00	-6.00
1.53 SFSV																		
Shamir Autograph® III Single Vision	GR	BR								GR	BR			GR			-14.00 to +8.00	-6.50
Shamir Relax™	GR	BR								GR	BR			GR			-14.00 to +8.00	-6.50
1.59 SFSV																		
Shamir Autograph® III Single Vision	GR	BR								GR	BR			GR			-14.00 to +8.75	-7.00
Shamir Relax™	GR	BR								GR	BR			GR			-14.75 to +6.00	-7.00
1.60 SFSV																		
Shamir Autograph® III Single Vision	GR	BR								GR							-14.00 to +8.75	-7.00
Shamir Relax™	GR	BR								GR							-14.00 to +8.75	-7.00
1.67 SFSV																		
Shamir Autograph® III Single Vision	GR	BR								GR	BR			GR			-16.00 to +16.00	-8.00
Shamir Relax™	GR	BR								GR	BR			GR			-16.75 to +16.00	-8.00
1.74 SFSV																		
Shamir Autograph® III Single Vision	GR	GR															-20.00 to +13.50	-9.00

	Transitions Gen#S									Transitions XTRACTIVE™			Transitions XTRACTIVE™ POLARIZED™			Sphere	Add	Cyl to	MFH	
Blended Bifocal (digital)																				
1.50 Standard Index Blended Bifocal (digital)																				
Shamir Duo	GR	BR								GR							-9.00 to +6.00	+0.75 to +4.00	-6.00	15
1.53 Trivex® Blended Bifocal (digital)																				
Shamir Duo	GR	BR								GR							-13.25 to +6.00	+0.75 to +4.00	-6.50	15
1.59 Polycarbonate Blended Bifocal (digital)																				
Shamir Duo	GR	BR								GR							-14.75 to +6.00	+0.75 to +4.00	-7.00	15
1.60 High Index Blended Bifocal (digital)																				
Shamir Duo	GR	BR								GR							-15.00 to +6.00	+0.75 to +4.00	-7.00	15
1.67 Super High Index Blended Bifocal (digital)																				
Shamir Duo	GR	BR								GR							-16.75 to +8.00	+0.75 to +4.00	-8.00	15
1.74 Ultra High Index Blended Bifocal (digital)																				
Shamir Duo	GR	BR															-18.75 to +12.00	+0.75 to +4.00	-9.00	15

Transitions XTRACTIVE™

Transitions XTRACTIVE™ POLARIZED™



Progressives										Sphere	Add	Cyl to	MFH		
1.50 Standard Index															
Shamir Autograph Intelligence™	GR	BR								GR	GR	-12.25 to +6.50	+0.75 to +4.00	-6.00	11, 13,15,18
Shamir Autograph® III	GR	BR								GR	GR	-12.25 to +6.00	+0.75 to +4.00	-6.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR								GR	GR	-12.25 to +6.00	+0.75 to +4.00	-6.00	11, 13, 15, 18
Shamir Intouch™	GR	BR								GR	GR	-12.25 to +6.00	+0.75 to +3.50	-6.00	14, 16, 18
Shamir Spectrum+™	GR	BR								GR	GR	-12.25 to +6.00	+0.75 to +3.50	-6.00	14, 16, 18
Shamir Element™	GR	BR								GR	GR	-7.00 to +4.00	+0.75 to +3.50	-6.00	16, 19
Shamir GENESIS HD	GR	BR								GR	GR	-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir Computer™	GR	BR								GR	GR	-7.00 to +4.00	+0.75 to +2.00	-6.00	14 Lower, 10 Upper
Shamir Workspace™	GR	BR								GR	GR	-6.00 to +6.00	+0.75 to +2.50	-6.00	14 Lower, 10 Upper
Shamir Attitude III Sport	GR	BR								GR	GR	-6.50 to +5.25	+0.75 to +3.50	-6.00	18
Shamir Attitude III Fashion	GR	BR								GR	GR	-6.50 to +5.25	+0.75 to +3.50	-6.00	15, 18
1.53 TRIVEX®															
Shamir Autograph Intelligence™	GR	BR								GR	GR	-13.25 to +7.00	+0.75 to +4.00	-6.50	11, 13, 15, 18
Shamir Autograph® III	GR	BR								GR	GR	-13.25 to +7.00	+0.75 to +4.00	-6.50	11, 13, 15, 18
Shamir Autograph® II+	GR	BR								GR	GR	-13.25 to +6.00	+0.75 to +4.00	-6.50	11, 13, 15, 18
Shamir Intouch™	GR	BR								GR	GR	-13.25 to +6.00	+0.75 to +3.50	-6.50	14, 16, 18
Shamir Golf™	GR	BR								GR	GR	-5.00 to +4.00	+1.00 to +3.50	-6.50	19
Shamir Spectrum+™	GR	BR								GR	GR	-13.25 to +6.00	+0.75 to +3.50	-6.50	14, 16, 18
Shamir Element™	GR	BR								GR	GR	-13.25 to +6.00	+0.75 to +3.50	-6.50	16, 19
Shamir GENESIS HD	GR	BR								GR	GR	-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir Computer™	GR	BR								GR	GR	-6.00 to +8.00	+0.75 to +2.00	-6.50	14 Lower, 10 Upper
Shamir Workspace™	GR	BR								GR	GR	-6.00 to +8.00	+0.75 to +2.50	-6.50	14 Lower, 10 Upper
Shamir Attitude III Sport	GR	BR								GR	GR	-6.50 to +6.00	+0.75 to +3.50	-6.50	18
Shamir Attitude III Fashion	GR	BR								GR	GR	-6.50 to +6.00	+0.75 to +3.50	-6.50	15, 18
1.59 Polycarbonate															
Shamir Autograph Intelligence™	GR	BR								GR	GR	-14.75 to +7.50	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR								GR	GR	-14.75 to +8.00	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR								GR	GR	-14.75 to +6.50	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Intouch™	GR	BR								GR	GR	-14.75 to +6.00	+0.75 to +3.50	-7.00	15, 18
Shamir Golf™	GR	BR								GR	GR	-6.00 to +5.25	+1.00 to +3.50	-7.00	19
Shamir Spectrum+™	GR	BR								GR	GR	-14.75 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
Shamir Element™	GR	BR								GR	GR	-10.00 to +6.00	+0.75 to +3.50	-7.00	16, 19
Shamir GENESIS HD	GR	BR								GR	GR	-10.00 to +6.00	+0.75 to +3.50	-7.00	16
Shamir Computer™	GR	BR								GR	GR	-6.00 to +8.00	+0.75 to +2.00	-7.00	14 Lower, 10 Upper
Shamir Workspace™	GR	BR								GR	GR	-6.00 to +8.00	+0.75 to +2.50	-7.00	14 Lower, 10 Upper
Shamir Attitude III Sport	GR	BR								GR	GR	-6.50 to +6.00	+0.75 to +3.50	-6.50	18
Shamir Attitude III Fashion	GR	BR								GR	GR	-6.50 to +6.00	+0.75 to +3.50	-6.50	15, 18
1.60 High Index															
Shamir Autograph Intelligence™	GR	BR								GR	GR	-14.75 to +8.75	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR								GR	GR	-15.00 to +8.00	+0.75 to +3.50	-7.00	15,18
Shamir Autograph® II+	GR	BR								GR	GR	-15.00 to +6.00	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Golf™ SuperLite™ 1.6	GR	BR								GR	GR	-6.00 to +5.25	+1.00 to +3.50	-7.00	19
Shamir Spectrum+™	GR	BR								GR	GR	-15.00 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
Shamir Element™	GR	BR								GR	GR	-10.00 to +6.00	+0.75 to +3.50	-7.00	16, 19
Shamir GENESIS HD	GR	BR								GR	GR	-6.00 to +6.00	+0.75 to +3.50	-7.00	16
Shamir Computer™	GR	BR								GR	GR	-6.00 to +8.00	+0.75 to +2.00	-7.00	14 Lower, 10 Upper
Shamir Workspace™	GR	BR								GR	GR	-6.00 to +8.00	+0.75 to +2.50	-7.00	14 Lower, 10 Upper
Shamir InTouch™	GR	BR								GR	GR	-15.00 to +6.00	+0.75 to +3.50	-7.00	15, 18
Shamir Attitude III Sport	GR	BR								GR	GR	-6.50 to +6.00	+0.75 to +3.50	-7.00	18

Shamir Attitude III Fashion	GR	BR					GR						-6.50 to +6.00	+0.75 to +3.50	-7.00	15, 18
1.67 Super High Index																
Shamir Autograph Intelligence™	GR	BR					GR		GR				-16.75 to +10.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR					GR		GR				-16.75 to +9.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR					GR		GR				-16.75 to +8.25	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Spectrum+™	GR	BR					GR		GR				-16.75 to +8.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Element™ SuperLite™ 1.67	GR	BR					GR		GR				-12.50 to +8.00	+0.75 to +3.50	-8.00	16, 19
Shamir GENESIS HD	GR	BR					GR		GR				-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir InTouch™	GR	BR					GR		GR				-16.75 to +8.00	+0.75 to +3.50	-8.00	15, 18
Shamir Attitude III Sport	GR	BR					GR		GR				-6.50 to +6.00	+0.75 to +3.50	-8.00	18
Shamir Attitude III Fashion	GR	BR					GR		GR				-6.50 to +6.00	+0.75 to +3.50	-8.00	15, 18
1.74 Ultra High Index																
Shamir Autograph Intelligence™	GR	BR											-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR											-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR											-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir GENESIS HD	GR	BR											-6.00 to +6.00	+0.75 to +4.00	-4.00	16
Shamir Spectrum+™	GR	BR											-18.75 to +11.50	+0.75 to +3.50	-9.00	14, 16, 18
Shamir InTouch™	GR	BR											-18.75 to +15.50	+0.75 to +3.50	-9.00	15, 18

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, *Transitions Signature*, and *XTRActive* are registered trademarks, and *XTRActive Polarized*, the *Transitions* logo and *Transitions Light Intelligent Lenses* are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.