





# Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Design](#) / [Progressives](#) / [Signet Armorlite](#)

## Signet Armorlite – Progressives

													Sphere	Add	Cyl to	MFH	
	Progressives																
<b>1.50 Standard Index</b>																	
KODAK Precise® 1.50 SF	GR													-9.00 to +5.00	+1.00 to +3.00	-4.00	17
KODAK Precise® Digital 1.50	GR													-8.00 to +5.00	+1.00 to +3.00	-6.00	17
KODAK Precise® Plus 1.50	GR	BR	GG	S	AME	AMB	E	R	GR					-8.00 to +4.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus 1.50	GR	BR	GG	S	AME	AMB	E	R	GR					-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Precise® PB 1.50	GR	BR							GR					-8.00 to +4.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB 1.50	GR	BR							GR					-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique™ 1.50	GR	BR	GG	S	AME	AMB	E	R	GR					-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique™ HD 1.50	GR	BR												-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ 1.50	GR	BR	GG	S	AME	AMB	E	R	GR					-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.50	GR	BR	GG	S	AME	AMB	E	R	GR					-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Easy 14 1.50	GR	BR	GG	S	AME	AMB	E	R						-8.00 to +4.00	+0.75 to +3.50	-6.00	14
KODAK Easy 18 1.50	GR	BR	GG	S	AME	AMB	E	R						-8.00 to +4.00	+0.75 to +3.50	-6.00	18
DirecTek 1.50	GR	BR												-8.00 to +4.00	+0.75 to +3.50	-6.00	16.5
DirecTek Short 1.50	GR	BR												-8.00 to +4.00	+0.75 to +3.50	-6.00	13.5
Navigator FBS 1.50	GR	BR												-8.00 to +4.00	+0.75 to +3.50	-6.00	18
Navigator Short FBS 1.50	GR	BR												-8.00 to +4.00	+0.75 to +3.50	-6.00	14
Crossbows Custom U 1.50	GR	BR												-8.00 to +4.00	+0.75 to +3.50	-6.00	13
Crossbows Custom V 1.50	GR	BR												-8.00 to +4.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 13 Progressive 1.50	GR	BR												-8.00 to +4.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 17 Progressive 1.50	GR	BR												-8.00 to +4.00	+0.75 to +3.50	-6.00	17
Crossbows Junior	GR	BR												-8.00 to +4.00	+0.75 to +3.50	-6.00	13
<b>1.53 TRIVEX®</b>																	
KODAK Precise® PB Trivex®	GR	BR							GR					-10.00 to +4.50	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB Trivex®	GR	BR							GR					-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Precise® Plus Trivex®	GR	BR							GR					-10.00 to +4.50	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus Trivex®	GR	BR							GR					-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Unique Trivex®	GR	BR							GR					-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Unique HD Trivex®	GR	BR												-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ Trivex®	GR	BR							GR					-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD Trivex®	GR	BR							GR					-10.00 to +4.50	+0.75 to +3.50	-6.00	13
Crossbows Custom U Trivex®	GR	BR												-10.00 to +4.50	+0.75 to +3.50	-6.00	13
Crossbows Custom V Trivex®	GR	BR												-10.00 to +4.50	+0.75 to +3.50	-6.00	13
<b>1.59 Polycarbonate</b>																	
KODAK Precise® Poly SF	GR													-10.00 to +5.00	+1.00 to +3.00	-4.00	17
KODAK Precise® Digital Poly	GR													-10.00 to +6.00	+1.00 to +3.00	-6.00	17
KODAK Precise® Plus Poly	GR	BR	GG	S	AME	AMB	E	R	GR					-10.00 to +5.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus Poly	GR	BR	GG	S	AME	AMB	E	R	GR					-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ Poly	GR	BR	GG	S	AME	AMB	E	R	GR			GR		-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Precise® PB Poly	GR	BR	GG						GR			GR		-10.00 to +5.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB Poly	GR	BR	GG						GR			GR		-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique Poly	GR	BR	GG	S	AME	AMB	E	R	GR			GR		-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique HD Poly	GR	BR												-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD Poly	GR	BR	GG	S	AME	AMB	E	R	GR			GR		-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Easy 14 Poly	GR	BR	GG	S	AME	AMB	E	R						-10.00 to +5.00	+0.75 to +3.50	-6.00	14

KODAK Easy 18 Poly	GR	BR	GG	S	AME	AMB	E	R										-10.00 to +5.00	+0.75 to +3.50	-6.00	18
DirectTek Poly	GR	BR																-10.00 to +5.00	+0.75 to +3.50	-6.00	16.5
DirectTek Short Poly	GR	BR																-10.00 to +5.00	+0.75 to +3.50	-6.00	13.5
Navigator FBS Poly	GR	BR																-10.00 to +5.00	+0.75 to +3.50	-6.00	18
Navigator Short FBS Poly	GR	BR																-10.00 to +5.00	+0.75 to +3.50	-6.00	14
Crossbows Custom U Poly	GR	BR																-10.00 to +5.00	+0.75 to +3.50	-6.00	13
Crossbows Custom V Poly	GR	BR																-10.00 to +5.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 13 Progressive Poly	GR	BR																-10.00 to +5.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 17 Progressive Poly	GR	BR																-10.00 to +5.00	+0.75 to +3.50	-6.00	17
Crossbows Junior Poly	GR	BR																-10.00 to +5.00	+0.75 to +3.50	-6.00	13
<b>1.60 High Index</b>																					
KODAK Unique	GR	BR																-10.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique HD	GR	BR																-10.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ 1.60	GR	BR																-10.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.60	GR	BR																-10.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Custom U 1.60	GR	BR																-10.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Custom V 1.60	GR	BR																-10.00 to +6.00	+0.75 to +3.50	-6.00	13
<b>1.67 Super High Index</b>																					
KODAK Precise® PB 1.67	GR	BR	GG								GR							-11.00 to +6.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB 1.67	GR	BR	GG								GR							-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Precise® Plus Poly	GR	BR	GG	S	AME	AMB	E	R	GR									-11.00 to +6.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus Poly	GR	BR	GG	S	AME	AMB	E	R	GR									-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique	GR	BR	GG	S	AME	AMB	E	R	GR									-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique HD	GR	BR																-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ 1.67	GR	BR	GG	S	AME	AMB	E	R	GR									-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.67	GR	BR	GG	S	AME	AMB	E	R	GR									-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Easy 14 Poly	GR	BR	GG	S	AME	AMB	E	R										-11.00 to +6.00	+0.75 to +3.50	-6.00	14
KODAK Easy 18 Poly	GR	BR	GG	S	AME	AMB	E	R										-11.00 to +6.00	+0.75 to +3.50	-6.00	18
DirectTek 1.67	GR	BR																-11.00 to +6.00	+0.75 to +3.50	-6.00	16.5
DirectTek Short 1.67	GR	BR																-11.00 to +6.00	+0.75 to +3.50	-6.00	13.5
Navigator FBS 1.67	GR	BR																-11.00 to +6.00	+0.75 to +3.50	-6.00	18
Navigator Short FBS 1.67	GR	BR																-11.00 to +6.00	+0.75 to +3.50	-6.00	14
Crossbows Custom U 1.67	GR	BR																-11.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Custom V 1.67	GR	BR																-11.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 13 Progressive 1.67	GR	BR																-11.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 17 Progressive 1.67	GR	BR																-11.00 to +6.00	+0.75 to +3.50	-6.00	17
Crossbows Junior 1.67	GR	BR																-11.00 to +6.00	+0.75 to +3.50	-6.00	13
<b>1.74 Ultra High Index</b>																					
KODAK Unique	GR	BR																-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique HD	GR	BR																-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™	GR	BR																-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.74	GR	BR																-13.00 to +9.00	+0.75 to +3.50	-6.00	13

**TRANSITIONS® LIGHT INTELLIGENT LENSES™**



**GR**  
Gray



**BR**  
Brown



**GG**  
Graphite  
Green



**S**  
Sapphire



**AME**  
Amethyst



**AMB**  
Amber



**E**  
Emerald



**R**  
Ruby



Available  
in style  
mirrors

All other trademarks are the property of their respective owners.