



Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Design](#) / [Progressive](#) / [Unity Lenses](#)

Unity Lenses – Progressives

 Progressives	 Transitions Gen#S								 Transitions XTRACTIVE			 Transitions XTRACTIVE POL-MILKIE	Sphere	Add	Cyl to	MFH	
	GR	BR															
1.50 Standard Index																	
Unity Via II	GR	BR								GR			GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR								GR			GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR								GR			GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR	BR								GR			GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR								GR			GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Ethos	GR	BR								GR			GR	-9.00 to +4.50	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR								GR			GR	-9.00 to +4.50	0.75 to +3.50	-4.00	14
1.53 TRIVEX®																	
Unity Via II	GR	BR								GR				-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR								GR				-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR								GR				-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR	BR								GR				-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR								GR				-5.00 to +4.00	+0.50 to +4.50	-6.00	12
Ethos	GR	BR								GR				-9.00 to +6.00	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR								GR				-9.00 to +6.00	0.75 to +3.50	-4.00	14
1.59 Polycarbonate																	
Unity Via II	GR	BR								GR			GR	-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR								GR			GR	-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR								GR			GR	-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR	BR								GR			GR	-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR								GR			GR	-6.50 to +6.00	+0.50 to +4.50	-6.00	12
Ethos	GR	BR								GR			GR	-9.00 to +6.00	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR								GR			GR	-9.00 to +6.00	0.75 to +3.50	-4.00	14
1.60 High Index																	
Unity Via II	GR	BR												-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR												-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR												-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR	BR												-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR												-6.50 to +5.25	+0.50 to +4.50	-6.00	12
Ethos	GR	BR												-9.00 to +5.00	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR												-9.00 to +5.00	0.75 to +3.50	-4.00	14
1.67 Super High Index																	
Unity Via II	GR	BR								GR			GR	-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR								GR			GR	-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR								GR			GR	-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR	BR								GR			GR	-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR								GR			GR	-6.00 to +5.00	+0.50 to +4.50	-6.00	12
Ethos	GR	BR								GR			GR	-12.00 to +8.00	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR								GR			GR	-12.00 to +8.00	0.75 to +3.50	-4.00	14

TRANSITIONS® LIGHT INTELLIGENT LENSES™



GR Gray	BR Brown	GG Graphite Green	S Sapphire	AME Amethyst	AMB Amber	E Emerald	R Ruby	Available in style mirrors
-------------------	--------------------	--------------------------------	----------------------	------------------------	---------------------	---------------------	------------------	----------------------------------

Transitions, *Transitions Signature*, and *XTRActive* are registered trademarks, and *XTRActive Polarized*, the *Transitions* logo and *Transitions Light Intelligent Lenses* are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.