

















Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Design](#) / [Progressives](#)

All Suppliers – Progressives







 Progressives	 Transitions Gen#S								 Transitions XTRACTIVE			 Transitions XTRACTIVE PG-ARIZEP	Sphere	Add	Cyl to	MFH	
	GR	BR	GG														
1.50 Standard Index																	
Varilux® XR Track™	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Varilux® Liberty™ 3.0	GR	BR												-10.00 to +7.50	+0.75 to +1.50	-6.00	14
Varilux Comfort®	GR	BR							GR	BR	GG			-10.00 to +7.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short™	GR	BR							GR	BR	GG			-10.00 to +7.00	+0.75 to +3.00	-6.00	14
Varilux® Stylistic™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ FIT	GR	BR	GG	S	AME	AMB	E	R						-10.00 to +7.50	+0.75 to +3.50	-4.00	17
Essilor Ideal Advanced™ Wrap	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-4.00	17
Essilor Ideal® FIT	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	18
Essilor Ideal® Driving	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	13
Essilor Ideal®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Ideal® Short	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Natural® Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.50	-6.00	15
Essilor Accolade™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Brillenz®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Brillenz® Short	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14

Ovation® Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Ovation®	GR	BR										-10.00 to +6.00	+0.75 to +3.50	-4.00	17
Essilor SmallFit® Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.00	-6.00	14














	Transitions Gen#S								Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED					
Progressives												Sphere	Add	Cyl to	MFH		
1.53 TRIVEX®																	
Varilux® XR Track™	GR	BR										-10.00 to +6.00	+0.75 to +4.00	-6.00		NVB Measurement	
Varilux® XR Design™	GR	BR										-10.00 to +6.00	+0.75 to +4.00	-6.00		14	
Varilux® X 4D™ TREXA™	GR	BR										-10.00 to +6.00	+0.75 to +4.00	-6.00		14	
Varilux® X FIT™ TREXA™	GR	BR										-10.00 to +6.00	+0.75 to +4.00	-6.00		14	
Varilux® X Design™ TREXA™	GR	BR										-10.00 to +6.00	+0.75 to +4.00	-6.00		14	
Varilux® Physio® W3+ Eyecode TREXA™	GR	BR							GR	BR		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00		14	
Varilux® Physio® W3+ FIT TREXA™	GR	BR							GR	BR		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00		14	
Varilux® Physio® W3+ TREXA™	GR	BR							GR	BR		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00		14	
Varilux Comfort® Max Fit TREXA™	GR	BR							GR	BR		-10.00 to +6.00	+0.75 to +3.50	-6.00		14	
Varilux Comfort® Max TREXA™	GR	BR							GR	BR		-10.00 to +6.00	+0.75 to +3.50	-6.00		14	
Varilux Physio DRx™ TREXA™	GR	BR	GG						GR	BR		-10.00 to +7.50	+0.75 to +3.50 & +4.00	-6.00		17	
Varilux Physio Short DRx™ TREXA™	GR	BR	GG						GR	BR		-10.00 to +7.50	+0.75 to +3.00	-6.00		14	
Varilux Comfort DRx™ TREXA™	GR	BR	GG						GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00		17	
Varilux Comfort Short DRx™ TREXA™	GR	BR	GG						GR	BR		-10.00 to +7.50	+0.75 to +3.00	-6.00		14	
Essilor Ideal Advanced™ FIT TREXA™	GR	BR	GG									-10.00 to +7.50	+0.75 to +3.50	-6.00		17	
Essilor Ideal Advanced™ TREXA™	GR	BR	GG						GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00		17	
Essilor Ideal® FIT TREXA™	GR	BR	GG						GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00		17	
Essilor Ideal® Sport TREXA™	GR	BR	GG						GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00		18	
Essilor Ideal® Computer TREXA™	GR	BR	GG						GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00		13	
Essilor Ideal® Driving TREXA™	GR	BR	GG						GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00		14	
Essilor Ideal® TREXA™	GR	BR	GG						GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00		14	
Essilor Ideal® Short TREXA™	GR	BR	GG						GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00		14	
Essilor Natural® Digital TREXA™	GR	BR	GG						GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00		17	
Essilor Accolade Freedom™ TREXA™	GR	BR	GG						GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00		15	
Essilor Accolade® TREXA™	GR	BR	GG						GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00		17	
Adaptar® Digital TREXA™	GR	BR	GG						GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00		17	
Adaptar® Short Digital TREXA™	GR	BR	GG						GR	BR		-10.00 to +7.50	+0.75 to +3.00	-6.00		14	
Brillenz® TREXA™	GR	BR	GG						GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00		14	
Brillenz® Short TREXA™	GR	BR	GG						GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00		14	
Ovation® Digital TREXA™	GR	BR	GG						GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00		17	
Essilor SmallFit® Digital TREXA™	GR	BR	GG						GR	BR		-10.00 to +7.50	+0.75 to +3.00	-6.00		14	

	Transitions Gen#S								Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED					
Progressives												Sphere	Add	Cyl to	MFH		
1.59 Polycarbonate																	
Varilux® XR Track™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00		NVB Measurement	








Varilux® XR Design™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Varilux Comfort DRx™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux® Liberty™ 3.0	GR	BR											-10.00 to +8.00	+0.75 to +1.50	-6.00	14
Varilux Comfort® Airwear®	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short™ Airwear®	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +3.00	-6.00	14
Varilux® Stylistic™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ FIT	GR	BR	GG	S	AME	AMB	E	R					-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Wrap	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Essilor Ideal® FIT Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Sport Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Driving Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor Ideal® Short Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor Natural® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +8.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Ovation® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Ovation® Airwear®	GR	BR											-10.00 to +8.00	+0.75 to +3.50	-4.00	17
Essilor SmallFit® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +8.00	+0.75 to +3.00	-6.00	14

															
Progressives											Sphere	Add	Cyl to	MFH	
1.60 High Index															
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.60	GR	BR	GG									-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Thin&Lite® 1.60	GR	BR	GG									-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Thin&Lite® 1.60	GR	BR	GG									-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14

















Varilux Comfort® Max Fit Thin&Lite® 1.60	GR	BR	GG														-12.00 to +8.00	+0.75 to +3.50	-6.00	14	
Varilux Comfort® Max Thin&Lite® 1.60	GR	BR	GG															-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Thin&Lite® 1.60	GR	BR																-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Thin&Lite® 1.60	GR	BR																-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™ Thin&Lite® 1.60	GR	BR																-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Thin&Lite® 1.60	GR	BR																-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux Comfort® Thin&Lite® 1.60	GR																	-12.00 to +8.00	+0.75 to +3.50	-4.00	17
Varilux Comfort Short™ Thin&Lite® 1.60	GR																	-12.00 to +8.00	+0.75 to +3.00	-4.00	14
Essilor Ideal Advanced™ FIT Thin&Lite® 1.60	GR	BR																-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Thin&Lite® 1.60	GR	BR																-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.60	GR	BR																-12.00 to +8.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Sport Thin&Lite® 1.60	GR	BR																-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Thin&Lite® 1.60	GR	BR																-12.00 to +8.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Driving Thin&Lite® 1.60	GR	BR																-12.00 to +8.00	+0.75 to +3.50	-6.00	13
Essilor Ideal® Thin&Lite® 1.60	GR	BR																-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Short Thin&Lite® 1.60	GR	BR																-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Essilor Natural® Digital Thin&Lite® 1.60	GR	BR																-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Thin&Lite® 1.60	GR	BR																-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Thin&Lite® 1.60	GR	BR																-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Ovation® Digital Thin&Lite® 1.60	GR	BR																-12.00 to +8.00	+0.75 to +3.50	-6.00	17
SmallFit® Digital Thin&Lite® 1.60	GR	BR																-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Thin&Lite® 1.60	GR	BR																-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Thin&Lite® 1.60	GR	BR																-12.00 to +8.00	+0.75 to +3.00	-6.00	14

	Transitions Gen#S								Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED														
Progressives																										
1.67 Super High Index																										
Varilux® XR Track™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG												-12.00 to +9.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG												-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG												-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG												-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG												-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG												-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG												-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG												-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG												-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG												-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR											-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR											-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR											-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR											-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Varilux® Liberty™ 3.0	GR	BR																					-14.00 to +9.00	+0.75 to +1.50	-6.00	14
Varilux Comfort® Thin&Lite® 1.67	GR	BR																					-14.00 to +9.00	+0.75 to +3.00	-6.00	14

Varilux Comfort® Short™ Thin&Lite® 1.67	GR	BR												-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Varilux® Stylistic™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-5.00 to +5.00	+0.75 to +3.50	-3.00	17
Essilor Ideal Advanced™ FIT Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R						-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal Advanced Wrap™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Sport Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Driving Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Essilor Ideal® Short	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Essilor Natural® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Brillenz® Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Ovation® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Ovation® Thin&Lite® 1.67	GR	BR												-12.00 to +8.00	+0.75 to +3.50	-4.00	17
Essilor SmallFit® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.00	-6.00	14
















																	
Progressives										Sphere		Add		Cyl to		MFH	
1.74 Ultra High Index																	
Varilux® XR Track™	GR	BR										-14.00 to +13.00	+0.75 to +4.00	-6.00		NVB Measurement	
Varilux® XR Design™	GR	BR										-14.00 to +13.00	+0.75 to +4.00	-6.00		14	
Varilux® X 4D™ Thin&Lite® 1.74	GR	BR										-14.00 to +13.00	+0.75 to +4.00	-6.00		14	
Varilux® X FIT™ Thin&Lite® 1.74	GR	BR										-14.00 to +13.00	+0.75 to +4.00	-6.00		14	
Varilux® X Design™ Thin&Lite® 1.74	GR	BR										-14.00 to +13.00	+0.75 to +4.00	-6.00		14	
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.74	GR	BR										-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00		14	
Varilux® Physio® W3+ FIT Thin&Lite® 1.74	GR	BR										-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00		14	
Varilux® Physio® W3+ Thin&Lite® 1.74	GR	BR										-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00		14	
Varilux Physio DRx™ Thin&Lite® 1.74	GR	BR										-14.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00		17	
Varilux Physio Short DRx™ Thin&Lite® 1.74	GR	BR										-14.00 to +13.00	+0.75 to +3.50	-6.00		14	
Varilux Comfort DRx™ Thin&Lite® 1.74	GR	BR										-20.00 to +13.00	+0.75 to +3.50	-6.00		17	
Varilux Comfort Short DRx™ Thin&Lite® 1.74	GR	BR										-20.00 to +13.00	+0.75 to +3.50	-6.00		17	
Essilor Ideal Advanced™ FIT Thin&Lite® 1.74	GR	BR										-20.00 to +13.00	+0.75 to +3.50	-6.00		17	
Essilor Ideal Advanced™ Thin&Lite® 1.74	GR	BR										-20.00 to +13.00	+0.75 to +3.50	-6.00		17	
Essilor Ideal® FIT Thin&Lite® 1.74	GR	BR										-20.00 to +13.00	+0.75 to +3.50	-6.00		17	
Essilor Ideal® Sport Thin&Lite® 1.74	GR	BR										-20.00 to +13.00	+0.75 to +3.50	-6.00		18	
Essilor Ideal® Computer Thin&Lite® 1.74	GR	BR										-20.00 to +13.00	+0.75 to +3.50	-6.00		18	
Essilor Ideal® Driving Thin&Lite® 1.74	GR	BR										-20.00 to +13.00	+0.75 to +3.50	-6.00		18	

Essilor Ideal® Thin&Lite® 1.74	GR	BR													-20.00 to +13.00	+0.75 to +3.50	-6.00	13
Ovation® Digital Thin&Lite® 1.74	GR	BR													-20.00 to +13.00	+0.75 to +3.50	-6.00	14
Essilor Accolade Freedom™ Thin&Lite® 1.74	GR	BR													-20.00 to +13.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Thin&Lite® 1.74	GR	BR													-14.00 to +9.00	+0.75 to +3.50	-6.00	17

																		
Progressives														Sphere	Add	Cyl to	MFH	
1.50 Standard Index																		
HOYALUX iD LifeStyle®4 Indoor 1.50	GR!	BR!												-8.00 to +6.00	+0.75 to +4.00	-6.00	14	
HOYALUX iD LifeStyle®4 Outdoor 1.50	GR!	BR!												-8.00 to +6.00	+0.75 to +4.00	-6.00	14	
HOYALUX iD LifeStyle®4 Urban 1.50	GR!	BR!												-8.00 to +6.00	+0.75 to +4.00	-6.00	14	
Amp HD3 1.50	GR!	BR!												-8.00 to +6.00	+0.75 to +4.00	-6.00	18	
Amp Mini HD3 1.50	GR!	BR!												-8.00 to +6.00	+0.75 to +4.00	-6.00	14	
HOYALUX ARRAY™ 11, 13, 15, 17, VL 1.50	GR	BR							GR	BR		GR		-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL 1.50	GR	BR							GR	BR		GR		-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	
HOYALUX ARRAY™ 2 11, 13, 15, 17, VL 1.50	GR	BR							GR	BR		GR		-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	
HOYALUX ARRAY™ 2 Wrap 11, 13, 15, 17, VL 1.50	GR	BR							GR	BR		GR		-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	
Amplitude BKS 1.50	GR	BR							GR	BR		GR		-10.00 to +3.75	+0.75 to +4.00	-6.00	19	
Amplitude Mini BKS 1.50	GR	BR							GR	BR		GR		-10.00 to +3.75	+0.75 to +4.00	-6.00	14	
Summit cd BKS 1.50	GR	BR							GR	BR		GR		-10.00 to +3.75	+0.75 to +4.00	-6.00	14	
Summit ecp BKS 1.50	GR	BR							GR	BR		GR		-10.00 to +3.75	+0.75 to +4.00	-6.00	18	
Hoya GP Wide 1.50	GR	BR							GR	BR		GR		-10.00 to +3.75	+0.75 to +4.00	-6.00	18	
Tact 40 BKS 1.50	GR	BR							GR	BR		GR		-10.00 to +3.75	+0.75 to +4.00	-6.00	18	
Tact 60 BKS 1.50	GR	BR							GR	BR		GR		-10.00 to +3.75	+0.75 to +4.00	-6.00	18	
1.53 TRIVEX®																		
HOYALUX iD LifeStyle®4 Indoor Phoenix™	GR!	BR!												-8.00 to +6.00	+0.75 to +4.00	-6.00	14	
HOYALUX iD LifeStyle®4 Outdoor Phoenix™	GR!	BR!												-8.00 to +6.00	+0.75 to +4.00	-6.00	14	
HOYALUX iD LifeStyle®4 Urban Phoenix™	GR!	BR!												-8.00 to +6.00	+0.75 to +4.00	-6.00	14	
HOYALUX iD LifeStyle®3 Indoor Phoenix™	GR!	BR!												-8.00 to +6.00	+0.75 to +4.00	-6.00	14	
HOYALUX iD LifeStyle®3 Outdoor Phoenix™	GR!	BR!												-8.00 to +6.00	+0.75 to +4.00	-6.00	14	
HOYALUX iD LifeStyle®3 Urban Phoenix™	GR!	BR!												-8.00 to +6.00	+0.75 to +4.00	-6.00	14	
Amp HD3 Phoenix™	GR!	BR!												-8.00 to +6.00	+0.75 to +4.00	-6.00	18	
Amp Mini HD3 Phoenix™	GR!	BR!												-8.00 to +6.00	+0.75 to +4.00	-6.00	14	
HOYALUX ARRAY™ 11, 13, 15, 17, VL Phoenix™	GR	BR							GR	BR		GR		-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL Phoenix™	GR	BR							GR	BR		GR		-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	
HOYALUX ARRAY 2 11, 13, 15, 17, VL Phoenix™	GR	BR							GR	BR		GR		-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	
HOYALUX ARRAY 2 Wrap 11, 13, 15, 17, VL Phoenix™	GR	BR							GR	BR		GR		-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	
Amplitude BKS Phoenix™	GR	BR							GR	BR		GR		-10.00 to +6.00	+0.75 to +4.00	-6.00	19	
Amplitude Mini BKS Phoenix™	GR	BR							GR	BR		GR		-10.00 to +6.00	+0.75 to +4.00	-6.00	14	
Summit cd BKS Phoenix™	GR	BR							GR	BR		GR		-10.00 to +6.00	+0.75 to +4.00	-6.00	18	
Summit ecp BKS Phoenix™	GR	BR							GR	BR		GR		-10.00 to +6.00	+0.75 to +4.00	-6.00	18	
Hoya GP Wide Phoenix™	GR	BR							GR	BR		GR		-10.00 to +6.00	+0.75 to +4.00	-6.00	18	
Tact 40 BKS Phoenix™	GR	BR							GR	BR		GR		-10.00 to +6.00	+0.75 to +4.00	-6.00	18	
Tact 60 BKS Phoenix™	GR	BR							GR	BR		GR		-10.00 to +6.00	+0.75 to +4.00	-6.00	18	







											+4.00		
1.59 Polycarbonate													
HOYALUX iD LifeStyle®4 Indoor Poly	GR!	BR!									-6.00 to +6.00	+0.75 to +4.00	-6.00 14
HOYALUX iD LifeStyle®4 Outdoor Poly	GR!	BR!									-6.00 to +6.00	+0.75 to +4.00	-6.00 14
HOYALUX iD LifeStyle®4 Urban Poly	GR!	BR!									-6.00 to +6.00	+0.75 to +4.00	-6.00 14
HOYALUX iD LifeStyle®3 Indoor Poly	GR!	BR!									-6.00 to +6.00	+0.75 to +4.00	-6.00 14
HOYALUX iD LifeStyle®3 Outdoor Poly	GR!	BR!									-6.00 to +6.00	+0.75 to +4.00	-6.00 14
HOYALUX iD LifeStyle®3 Urban Poly	GR!	BR!									-6.00 to +6.00	+0.75 to +4.00	-6.00 14
Amp HD3 Poly	GR!	BR!									-6.00 to +6.00	+0.75 to +4.00	-6.00 18
Amp Mini HD3 Poly	GR!	BR!									-6.00 to +6.00	+0.75 to +4.00	-6.00 14
HOYALUX ARRAY™ 11, 13, 15, 17, VL Poly	GR	BR					GR	BR		GR	-13.00 to +6.00	+0.75 to +4.00	-6.00 14, 17, 19, 21
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL Poly	GR	BR					GR	BR		GR	-6.00 to +4.00	+0.75 to +4.00	-6.00 14, 17, 19, 21
HOYALUX ARRAY 2 11, 13, 15, 17, VL Poly	GR	BR					GR	BR		GR	-13.00 to +6.00	+0.75 to +4.00	-6.00 14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11, 13, 15, 17, VL Poly	GR	BR					GR	BR		GR	-6.00 to +4.00	+0.75 to +4.00	-6.00 14, 17, 19, 21
Amplitude BKS Poly	GR	BR					GR	BR		GR	-13.00 to +6.00	+0.75 to +4.00	-6.00 19
Amplitude Mini BKS Poly	GR	BR					GR	BR		GR	-13.00 to +6.00	+0.75 to +4.00	-6.00 14
Summit cd BKS Poly	GR	BR					GR	BR		GR	-13.00 to +6.00	+0.75 to +4.00	-6.00 14
Summit ecp BKS Poly	GR	BR					GR	BR		GR	-13.00 to +6.00	+0.75 to +4.00	-6.00 18
Hoya GP Wide Poly	GR	BR					GR	BR		GR	-13.00 to +6.00	+0.75 to +4.00	-6.00 18
Tact 40 BKS Poly	GR	BR					GR	BR		GR	-13.00 to +6.00	+0.75 to +4.00	-6.00 18
Tact 60 BKS Poly	GR	BR					GR	BR		GR	-13.00 to +6.00	+0.75 to +4.00	-6.00 18
1.60 High Index													
HOYALUX iD LifeStyle®4 Indoor 1.60	GR!	BR!									-11.00 to +8.00	+0.75 to +4.00	-6.00 14
HOYALUX iD LifeStyle®4 Outdoor 1.60	GR!	BR!									-11.00 to +8.00	+0.75 to +4.00	-6.00 14
HOYALUX iD LifeStyle®4 Urban 1.60	GR!	BR!									-11.00 to +8.00	+0.75 to +4.00	-6.00 14
HOYALUX iD LifeStyle®3 Indoor 1.60	GR!	BR!									-11.00 to +8.00	+0.75 to +4.00	-6.00 14
HOYALUX iD LifeStyle®3 Outdoor 1.60	GR!	BR!									-11.00 to +8.00	+0.75 to +4.00	-6.00 14
HOYALUX iD LifeStyle®3 Urban 1.60	GR!	BR!									-11.00 to +8.00	+0.75 to +4.00	-6.00 14
Amp HD3 1.60	GR!	BR!									-11.00 to +8.00	+0.75 to +4.00	-6.00 14
Amp Mini HD3 1.60	GR!	BR!									-11.00 to +8.00	+0.75 to +4.00	-6.00 14
HOYALUX ARRAY 11, 13, 15, 17, VL 1.60	GR	BR									-10.00 to +6.00	+0.75 to +4.00	-6.00 14, 17, 19, 21
HOYALUX ARRAY Wrap 11,13,15,17,VL 1.60	GR	BR									-6.00 to +4.00	+0.75 to +4.00	-6.00 14, 17, 19, 21
HOYALUX ARRAY 2 11, 13, 15, 17, VL 1.60	GR	BR									-10.00 to +6.00	+0.75 to +4.00	-6.00 14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11,13,15,17,VL 1.60	GR	BR									-6.00 to +4.00	+0.75 to +4.00	-6.00 14, 17, 19, 21
Amplitude BKS 1.60	GR	BR									-11.00 to +8.00	+0.75 to +4.00	-6.00 19
Amplitude Mini BKS 1.60	GR	BR									-11.00 to +8.00	+0.75 to +4.00	-6.00 14
Summit cd BKS 1.60	GR	BR									-11.00 to +8.00	+0.75 to +4.00	-6.00 14
Summit ecp BKS 1.60	GR	BR									-11.00 to +8.00	+0.75 to +4.00	-6.00 18
Hoya GP Wide 1.60	GR	BR									-11.00 to +8.00	+0.75 to +4.00	-6.00 18
Tact 40 BKS 1.60	GR	BR									-11.00 to +8.00	+0.75 to +4.00	-6.00 18
Tact 60 BKS 1.60	GR	BR									-11.00 to +8.00	+0.75 to +4.00	-6.00 18
1.67 Super High Index													
HOYALUX iD LifeStyle®4 Indoor 1.67	GR!	BR!									-13.00 to +8.00	+0.75 to +4.00	-6.00 14
HOYALUX iD LifeStyle®4 Outdoor 1.67	GR!	BR!									-13.00 to +8.00	+0.75 to +4.00	-6.00 14
HOYALUX iD LifeStyle®4 Urban 1.67	GR!	BR!									-13.00 to +8.00	+0.75 to +4.00	-6.00 14

HOYALUX iD LifeStyle 3 Indoor 1.67	GR!	BR!																		-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle 3 Outdoor 1.67	GR!	BR!																		-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle 3 Urban 1.67	GR!	BR!																		-13.00 to +8.00	+0.75 to +4.00	-6.00	14
Amp HD3 1.67	GR!	BR!																		-13.00 to +8.00	+0.75 to +4.00	-6.00	18
Amp Mini HD3 1.67	GR!	BR!																		-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX ARRAY 11, 13, 15, 17, VL 1.67	GR	BR							GR	BR				GR						-13.00 to +7.50	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY Wrap 11, 13, 15, 17, VL 1.67	GR	BR							GR	BR				GR						-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13, 15, 17, VL 1.67	GR	BR							GR	BR				GR						-13.00 to +7.50	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11, 13, 15, 17, VL 1.67	GR	BR							GR	BR				GR						-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS 1.67	GR	BR							GR	BR				GR						-13.00 to +7.50	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS 1.67	GR	BR							GR	BR				GR						-13.00 to +7.50	+0.75 to +4.00	-6.00	14
Summit cd BKS 1.67	GR	BR							GR	BR				GR						-13.00 to +7.50	+0.75 to +4.00	-6.00	14
Summit ecp BKS 1.67	GR	BR							GR	BR				GR						-13.00 to +7.50	+0.75 to +4.00	-6.00	18
Hoya GP Wide 1.67	GR	BR							GR	BR				GR						-13.00 to +7.50	+0.75 to +4.00	-6.00	18
Tact 40 BKS 1.67	GR	BR							GR	BR				GR						-13.00 to +7.50	+0.75 to +4.00	-6.00	18
Tact 60 BKS 1.67	GR	BR							GR	BR				GR						-13.00 to +7.50	+0.75 to +4.00	-6.00	18
1.74 ultra High Index																							
HOYALUX iD LifeStyle®4 Indoor 1.74	GR!	BR!																		-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX iD LifeStyle®4 Outdoor 1.74	GR!	BR!																		-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX iD LifeStyle®4 Urban 1.74	GR!	BR!																		-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX iD LifeStyle 3 Indoor 1.74	GR!	BR!																		-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX iD LifeStyle 3 Outdoor 1.74	GR!	BR!																		-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX iD LifeStyle 3 Urban 1.74	GR!	BR!																		-13.00 to +8.00	+0.75 to +3.50	-6.00	14
Amp HD3 1.74	GR!	BR!																		-13.00 to +8.00	+0.75 to +3.50	-6.00	14
Amp Mini HD3 1.74	GR!	BR!																		-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX ARRAY 11, 13, 15, 17, VL 1.74	GR	BR																		-13.00 to +8.00	+0.75 to +3.50	-6.00	14, 17, 19, 21
HOYALUX ARRAY Wrap 11, 13, 15, 17, VL 1.74	GR	BR																		-6.00 to +4.00	+0.75 to +3.50	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13, 15, 17, VL 1.74	GR	BR																		-13.00 to +8.00	+0.75 to +3.50	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 WRAP 11, 13, 15, 17, VL 1.74	GR	BR																		-6.00 to +4.00	+0.75 to +3.50	-6.00	14, 17, 19, 21
Amplitude BKS 1.74	GR	BR																		-13.00 to +8.00	+0.75 to +3.50	-6.00	19
Amplitude Mini BKS 1.74	GR	BR																		-13.00 to +8.00	+0.75 to +3.50	-6.00	14
Summit cd BKS 1.74	GR	BR																		-13.00 to +8.00	+0.75 to +3.50	-6.00	14
Summit ecp BKS 1.74	GR	BR																		-13.00 to +8.00	+0.75 to +3.50	-6.00	19
Hoya GP Wide 1.74	GR	BR																		-13.00 to +8.00	+0.75 to +3.50	-6.00	18
Tact 40 BKS 1.74	GR	BR																		-13.00 to +8.00	+0.75 to +3.50	-6.00	18
Tact 60 BKS 1.74	GR	BR																		-13.00 to +8.00	+0.75 to +3.50	-6.00	18

																
Progressives												Sphere	Add	Cyl to	MFH	
1.50 Standard Index																
Nikon SeeMax Ultimate	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
Nikon Move	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.50 to	-6.00	15, 17

Shamir Intouch™	GR	BR							GR			GR	-13.25 to +6.00	+0.75 to +3.50	-6.50	14, 16, 18
Shamir Golf™	GR	BR							GR			GR	-5.00 to +4.00	+1.00 to +3.50	-6.50	19
Shamir Spectrum+™	GR	BR							GR			GR	-13.25 to +6.00	+0.75 to +3.50	-6.50	14, 16, 18
Shamir Element™	GR	BR							GR			GR	-13.25 to +6.00	+0.75 to +3.50	-6.50	16, 19
Shamir GENESIS HD	GR	BR							GR			GR	-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir Computer™	GR	BR							GR			GR	-6.00 to +8.00	+0.75 to +2.00	-6.50	14 Lower, 10 Upper
Shamir Workspace™	GR	BR							GR			GR	-6.00 to +8.00	+0.75 to +2.50	-6.50	14 Lower, 10 Upper
Shamir Attitude III Sport	GR	BR							GR			GR	-6.50 to +6.00	+0.75 to +3.50	-6.50	18
Shamir Attitude III Fashion	GR	BR							GR			GR	-6.50 to +6.00	+0.75 to +3.50	-6.50	15, 18
1.59 Polycarbonate																
Shamir Autograph Intelligence™	GR	BR							GR			GR	-14.75 to +7.50	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR							GR			GR	-14.75 to +8.00	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR							GR			GR	-14.75 to +6.50	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Intouch™	GR	BR							GR			GR	-14.75 to +6.00	+0.75 to +3.50	-7.00	15, 18
Shamir Golf™	GR	BR							GR			GR	-6.00 to +5.25	+1.00 to +3.50	-7.00	19
Shamir Spectrum+™	GR	BR							GR			GR	-14.75 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
Shamir Element™	GR	BR							GR			GR	-10.00 to +6.00	+0.75 to +3.50	-7.00	16, 19
Shamir GENESIS HD	GR	BR							GR			GR	-10.00 to +6.00	+0.75 to +3.50	-7.00	16
Shamir Computer™	GR	BR							GR			GR	-6.00 to +8.00	+0.75 to +2.00	-7.00	14 Lower, 10 Upper
Shamir Workspace™	GR	BR							GR			GR	-6.00 to +8.00	+0.75 to +2.50	-7.00	14 Lower, 10 Upper
Shamir Attitude III Sport	GR	BR							GR			GR	-6.50 to +6.00	+0.75 to +3.50	-6.50	18
Shamir Attitude III Fashion	GR	BR							GR			GR	-6.50 to +6.00	+0.75 to +3.50	-6.50	15, 18
1.60 High Index																
Shamir Autograph Intelligence™	GR	BR							GR				-14.75 to +8.75	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR							GR				-15.00 to +8.00	+0.75 to +3.50	-7.00	15, 18
Shamir Autograph® II+	GR	BR							GR				-15.00 to +6.00	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Golf™ SuperLite™ 1.6	GR	BR							GR				-6.00 to +5.25	+1.00 to +3.50	-7.00	19
Shamir Spectrum+™	GR	BR							GR				-15.00 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
Shamir Element™	GR	BR											-10.00 to +6.00	+0.75 to +3.50	-7.00	16, 19
Shamir GENESIS HD	GR	BR							GR				-6.00 to +6.00	+0.75 to +3.50	-7.00	16
Shamir Computer™	GR	BR							GR				-6.00 to +8.00	+0.75 to +2.00	-7.00	14 Lower, 10 Upper
Shamir Workspace™	GR	BR							GR				-6.00 to +8.00	+0.75 to +2.50	-7.00	14 Lower, 10 Upper
Shamir InTouch™	GR	BR							GR				-15.00 to +6.00	+0.75 to +3.50	-7.00	15, 18
Shamir Attitude III Sport	GR	BR							GR				-6.50 to +6.00	+0.75 to +3.50	-7.00	18
Shamir Attitude III Fashion	GR	BR							GR				-6.50 to +6.00	+0.75 to +3.50	-7.00	15, 18
1.67 Super High Index																
Shamir Autograph Intelligence™	GR	BR							GR			GR	-16.75 to +10.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR							GR				-16.75 to +9.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR							GR				-16.75 to +8.25	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Spectrum+™	GR	BR							GR				-16.75 to +8.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Element™ SuperLite™ 1.67	GR	BR							GR				-12.50 to +8.00	+0.75 to +3.50	-8.00	16, 19
Shamir GENESIS HD	GR	BR							GR				-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir InTouch™	GR	BR							GR				-16.75 to +8.00	+0.75 to +3.50	-8.00	15, 18
Shamir Attitude III Sport	GR	BR							GR				-6.50 to +6.00	+0.75 to +3.50	-8.00	18

Shamir Attitude III Fashion	GR	BR							GR		GR	-6.50 to +6.00	+0.75 to +3.50	-8.00	15, 18
1.74 Ultra High Index															
Shamir Autograph Intelligence™	GR	BR										-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR										-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR										-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir GENESIS HD	GR	BR										-6.00 to +6.00	+0.75 to +4.00	-4.00	16
Shamir Spectrum+™	GR	BR										-18.75 to +11.50	+0.75 to +3.50	-9.00	14, 16, 18
Shamir InTouch™	GR	BR										-18.75 to +15.50	+0.75 to +3.50	-9.00	15, 18

																			
Progressives													Sphere	Add	Cyl to	MFH			
1.50 Standard Index																			
KODAK Precise® 1.50 SF	GR												-9.00 to +5.00	+1.00 to +3.00	-4.00	17			
KODAK Precise® Digital 1.50	GR												-8.00 to +5.00	+1.00 to +3.00	-6.00	17			
KODAK Precise® Plus 1.50	GR	BR	GG	S	AME	AMB	E	R	GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	17			
KODAK Precise® Short Plus 1.50	GR	BR	GG	S	AME	AMB	E	R	GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	13			
KODAK Precise® PB 1.50	GR	BR							GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	17			
KODAK Precise® Short PB 1.50	GR	BR							GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	13			
KODAK Unique™ 1.50	GR	BR	GG	S	AME	AMB	E	R	GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	13			
KODAK Unique™ HD 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	13			
KODAK Unique DRO™ 1.50	GR	BR	GG	S	AME	AMB	E	R	GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	13			
KODAK Unique DRO™ HD 1.50	GR	BR	GG	S	AME	AMB	E	R	GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	13			
KODAK Easy 14 1.50	GR	BR	GG	S	AME	AMB	E	R					-8.00 to +4.00	+0.75 to +3.50	-6.00	14			
KODAK Easy 18 1.50	GR	BR	GG	S	AME	AMB	E	R					-8.00 to +4.00	+0.75 to +3.50	-6.00	18			
DirectTek 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	16.5			
DirectTek Short 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	13.5			
Navigator FBS 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	18			
Navigator Short FBS 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	14			
Crossbows Custom U 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	13			
Crossbows Custom V 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	13			
Crossbows Basic 13 Progressive 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	13			
Crossbows Basic 17 Progressive 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	17			
Crossbows Junior	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	13			
1.53 TRIVEX®																			
KODAK Precise® PB Trivex®	GR	BR							GR				-10.00 to +4.50	+0.75 to +3.50	-6.00	17			
KODAK Precise® Short PB Trivex®	GR	BR							GR				-10.00 to +4.50	+0.75 to +3.50	-6.00	13			
KODAK Precise® Plus Trivex®	GR	BR							GR				-10.00 to +4.50	+0.75 to +3.50	-6.00	17			
KODAK Precise® Short Plus Trivex®	GR	BR							GR				-10.00 to +4.50	+0.75 to +3.50	-6.00	13			
KODAK Unique Trivex®	GR	BR							GR				-10.00 to +4.50	+0.75 to +3.50	-6.00	13			
KODAK Unique HD Trivex®	GR	BR											-10.00 to +4.50	+0.75 to +3.50	-6.00	13			
KODAK Unique DRO™ Trivex®	GR	BR							GR				-10.00 to +4.50	+0.75 to +3.50	-6.00	13			
KODAK Unique DRO™ HD Trivex®	GR	BR							GR				-10.00 to +4.50	+0.75 to +3.50	-6.00	13			
Crossbows Custom U Trivex®	GR	BR											-10.00 to +4.50	+0.75 to +3.50	-6.00	13			
Crossbows Custom V Trivex®	GR	BR											-10.00 to +4.50	+0.75 to +3.50	-6.00	13			

Crossbows Custom U 1.67	GR	BR														-11.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Custom V 1.67	GR	BR														-11.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 13 Progressive 1.67	GR	BR														-11.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 17 Progressive 1.67	GR	BR														-11.00 to +6.00	+0.75 to +3.50	-6.00	17
Crossbows Junior 1.67	GR	BR														-11.00 to +6.00	+0.75 to +3.50	-6.00	13
1.74 Ultra High Index																			
KODAK Unique	GR	BR														-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique HD	GR	BR														-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™	GR	BR														-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.74	GR	BR														-13.00 to +9.00	+0.75 to +3.50	-6.00	13

Progressives													Sphere	Add	Cyl to	MFH			
1.50 Standard Index																			
Unity Via II	GR	BR							GR		GR					-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR							GR		GR					-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR							GR		GR					-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR	BR							GR		GR					-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR							GR		GR					-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Ethos	GR	BR							GR		GR					-9.00 to +4.50	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR							GR		GR					-9.00 to +4.50	0.75 to +3.50	-4.00	14
1.53 TRIVEX®																			
Unity Via II	GR	BR							GR							-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR							GR							-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR							GR							-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR	BR							GR							-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR							GR							-5.00 to +4.00	+0.50 to +4.50	-6.00	12
Ethos	GR	BR							GR							-9.00 to +6.00	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR							GR							-9.00 to +6.00	0.75 to +3.50	-4.00	14
1.59 Polycarbonate																			
Unity Via II	GR	BR							GR		GR					-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR							GR		GR					-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR							GR		GR					-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR	BR							GR		GR					-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR							GR		GR					-6.50 to +6.00	+0.50 to +4.50	-6.00	12
Ethos	GR	BR							GR		GR					-9.00 to +6.00	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR							GR		GR					-9.00 to +6.00	0.75 to +3.50	-4.00	14
1.60 High Index																			
Unity Via II	GR	BR														-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR														-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR														-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR	BR														-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR														-6.50 to +5.25	+0.50 to +4.50	-6.00	12
Ethos	GR	BR														-9.00 to +5.00	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR														-9.00 to +5.00	0.75 to +3.50	-4.00	14
1.67 Super High Index																			

Unity Via II	GR	BR								GR			GR	-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR								GR			GR	-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR								GR			GR	-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR	BR								GR			GR	-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR								GR			GR	-6.00 to +5.00	+0.50 to +4.50	-6.00	12
Ethos	GR	BR								GR			GR	-12.00 to +8.00	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR								GR			GR	-12.00 to +8.00	0.75 to +3.50	-4.00	14

YOUNGER OPTICS The Optical Lens Innovators			Transitions Gen#S					Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED					
Progressives													Sphere	Add	Cyl to	MFH
1.53 TRIVEX®																
Younger Trilogy® Image®	GR												-8.00 to +8.00	+1.00 to +3.00	-4.00	18
1.59 index polycarbonate																
Younger Image®	GR	BR								GR	BR		-9.00 to +7.00	+1.00 to +3.00	-4.00	18
1.67 Super High Index																
Younger Image® 1.67 MR-10™	GR	BR											-10.00 to +8.50	+1.00 to +3.50	-4.00	18

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

Unless indicated otherwise, all registered trademarks and trademarks are the property of Essilor International and/or its subsidiaries in the United States and in other countries.

All other trademarks are the property of their respective owners.