


# Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Design](#) / [SFSV](#) / [Nikon Optical USA](#)

## Nikon Optical USA – Semi-Finished Single Vision

 Semi-Finished Single Vision	Transitions Gen#S								Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED	Sphere	Cyl to	
	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR			
<b>1.50 Standard Index SFSV</b>															
Nikon SeeMax Infinite	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
Nikon Infinite Single Vision	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
Nikon Smart Boost	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
<b>1.53 TRIVEX® Single Vision</b>															
Nikon SeeMax Infinite	GR	BR	GG										-10.00 to +7.00	-6.00	
Nikon Infinite Single Vision	GR	BR	GG										-10.00 to +7.00	-6.00	
Nikon Smart Boost	GR	BR	GG										-10.00 to +7.00	-6.00	
<b>1.59 Polycarbonate Single Vision</b>															
Nikon SeeMax Infinite	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
Nikon Infinite Single Vision	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
Nikon Smart Boost	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
<b>1.60 High Index Single Vision</b>															
Nikon SeeMax Infinite	GR	BR											-12.00 to +8.00	-6.00	
Nikon Infinite Single Vision	GR	BR											-12.00 to +8.00	-6.00	
Nikon Smart Boost	GR	BR											-12.00 to +8.00	-6.00	
<b>1.67 Super High Index Single Vision</b>															
Nikon SeeMax Infinite	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	-6.00	
Nikon Infinite Single Vision	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	-6.00	
Nikon Smart Boost	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	-6.00	
<b>1.74 ultra High Index Single Vision</b>															
Nikon SeeMax Infinite	GR	BR											-16.00 to +10.00	-6.00	
Nikon Infinite Single Vision	GR	BR											-16.00 to +10.00	-6.00	
Nikon Smart Boost	GR	BR											-16.00 to +10.00	-6.00	

### TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.