

# Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Design](#) / [SFSV](#) / [Signet Armorlite](#)

## Signet Armorlite – Semi-Finished Single Vision

 Semi-Finished Single Vision	  									Sphere	Cyl to		
	GR	BR	GG	S	AME	AMB	E	R	GR				
<b>1.50 Standard Index SFSV</b>													
Signetek Processed SV	GR	BR	GG	S	AME	AMB	E	R	GR			-8.00 to +6.00	-6.00
KODAK Digital SV	GR	BR	GG	S	AME	AMB	E	R	GR			-8.00 to +6.00	-6.00
KODAK PowerUp Lens (.40/.66/1.00D)	GR	BR										-8.00 to +6.00	-6.00
Crossbows Custom SV	GR	BR										-8.00 to +6.00	-6.00
<b>1.53 TRIVEX® Single Vision</b>													
Signetek Processed SV	GR	BR	GG	S	AME	AMB	E	R	GR			-10.00 to +6.00	-6.00
KODAK Digital SV	GR	BR	GG	S	AME	AMB	E	R	GR			-10.00 to +6.00	-6.00
Crossbows Custom SV	GR	BR										-10.00 to +6.00	-6.00
<b>1.59 Polycarbonate Single Vision</b>													
Signetek Processed SV	GR	BR	GG	S	AME	AMB	E	R	GR	GR		-10.00 to +6.50	-6.00
KODAK Digital SV	GR	BR	GG	S	AME	AMB	E	R	GR	GR		-10.00 to +6.50	-6.00
KODAK PowerUp Lens (.40/.66/1.00D)	GR	BR										-10.00 to +6.50	-6.00
Crossbows Custom SV	GR	BR										-10.00 to +6.50	-6.00
<b>1.60 High Index Single Vision</b>													
Signetek Processed SV	GR	BR										-10.00 to +7.50	-6.00
KODAK Digital SV	GR	BR										-10.00 to +7.50	-6.00
Crossbows Custom SV	GR	BR										-10.00 to +7.50	-6.00
<b>1.67 Super High Index Single Vision</b>													
Signetek Processed SV	GR	BR	GG	S	AME	AMB	E	R	GR			-11.00 to +7.50	-6.00
KODAK Digital SV	GR	BR	GG	S	AME	AMB	E	R	GR			-11.00 to +7.50	-6.00
KODAK PowerUp Lens (.40/.66/1.00D)	GR	BR										-11.00 to +6.00	-6.00
Crossbows Custom SV	GR	BR										-11.00 to +7.50	-6.00
<b>1.74 Ultra High Index Single Vision</b>													
Signetek™ Processed SV	GR	BR										-13.00 to +9.00	-6.00
KODAK Digital SV	GR	BR										-13.00 to +9.00	-6.00

### TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.