





Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Transitions® XTRActive®](#) / **Signet Armorlite**

Signet Armorlite – Transitions® XTRActive®

		Transitions® XTRACTIVE		
Semi-Finished Single Vision			Sphere	Cyl to
1.50 Standard Index				
Signetek Processed SV	GR		-8.00 to +6.00	-6.00
KODAK Digital SV	GR		-8.00 to +6.00	-6.00
1.53 TRIVEX® Single Vision				
Signetek Processed SV	GR		-10.00 to +6.00	-6.00
KODAK Digital SV	GR		-10.00 to +6.00	-6.00
1.59 Polycarbonate Single Vision				
Signetek Processed SV	GR		-10.00 to +6.50	-6.00
KODAK Digital SV	GR		-10.00 to +6.50	-6.00
1.67 Super High Index Single Vision				
Signetek Processed SV	GR		-11.00 to +7.50	-6.00
KODAK Digital SV	GR		-11.00 to +7.50	-6.00

		Transitions® XTRACTIVE			
Bifocal Signature™			Sphere	Add	Cyl to
1.50 Standard Index SFST-28					
Signetek™ Processed 1.50 FT-28	GR		-10.00 to +6.00	+1.00 to +3.00	-6.00

		Transitions® XTRACTIVE				
Progressives			Sphere	Add	Cyl to	MFH
1.50 Standard Index						
KODAK Precise® Plus 1.50	GR		-8.00 to +4.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus 1.50	GR		-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Precise® PB 1.50	GR		-8.00 to +4.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB 1.50	GR		-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique™ 1.50	GR		-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ 1.50	GR		-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.50	GR		-8.00 to +4.00	+0.75 to +3.50	-6.00	13
1.53 TRIVEX®						
KODAK Precise® PB Trivex®	GR		-10.00 to +4.50	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB Trivex®	GR		-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Precise® Plus Trivex®	GR		-10.00 to +4.50	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus Trivex®	GR		-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Unique Trivex®	GR		-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ Trivex®	GR		-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD Trivex®	GR		-10.00 to +4.50	+0.75 to +3.50	-6.00	13
1.59 Polycarbonate						
KODAK Precise® Plus Poly	GR		-10.00 to +5.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus Poly	GR		-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ Poly	GR		-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Precise® PB Poly	GR		-10.00 to +5.00	+0.75 to +3.50	-6.00	17

KODAK Precise® Short PB Poly	GR			-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique Poly	GR			-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD Poly	GR			-10.00 to +5.00	+0.75 to +3.50	-6.00	13
1.67 Super High Index							
KODAK Precise® PB 1.67	GR			-11.00 to +6.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB 1.67	GR			-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Precise® Plus Poly	GR			-11.00 to +6.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus Poly	GR			-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique	GR			-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ 1.67	GR			-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.67	GR			-11.00 to +6.00	+0.75 to +3.50	-6.00	13

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Style Mirrors are available where gray and brown *Transitions® XTRActive®* are available. Specify *Transitions* lenses in style mirrors (**no substitutions**) with your lab to ensure authenticity.

Transitions, *Transitions Signature*, and *XTRActive* are registered trademarks, and *XTRActive Polarized*, the *Transitions* logo and *Transitions Light Intelligent Lenses* are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.