


Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Transitions® XTRActive® Polarized™](#) / [Hoya Vision Care](#)

Hoya Vision Care – Transitions® XTRActive® Polarized™

HOYA		Transitions® XTRACTIVE® POLARIZED™		
Semi-Finished Single Vision			Sphere	Cyl to
1.50 Standard Index				
	HOYA SFSV 1.50		-8.00 to +9.00	-4.00
	HOYA MySV 150	GR	-8.00 to +9.00	-4.00
	HOYA SV iQ 1.50	GR	-8.00 to +9.00	-4.00
	Sync III	GR	-10.00 to +6.00	-4.00
1.53 TRIVEX® Single Vision				
	HOYA MySV Phoenix™	GR	-13.00 to +9.00	-4.00
1.59 Polycarbonate Single Vision				
	HOYA SFSV Polycarbonate		-6.00 to +6.00	-4.00
	HOYA MySV Polycarbonate	GR	-10.00 to +6.00	-4.00
	HOYA SV iQ Polycarbonate	GR	-10.00 to +6.00	-4.00
	Sync III	GR	-13.00 to +5.25	-4.00
1.67 Super High Index Single Vision				
	HOYA MySV 1.67	GR	-10.00 to +8.00	-4.00
	HOYA SV iQ 1.67	GR	-10.00 to +8.00	-4.00
	Sync III	GR	-13.00 to +7.50	-4.00

HOYA		Transitions® XTRACTIVE® POLARIZED™				
Progressives			Sphere	Add	Cyl to	MFH
1.50 Standard Index						
	HOYALUX ARRAY™ 11, 13, 15, 17, VL 1.50	GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
	HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL 1.50	GR	-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
	HOYALUX ARRAY™ 2 11, 13, 15, 17, VL 1.50	GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
	HOYALUX ARRAY™ 2 Wrap 11, 13, 15, 17, VL 1.50	GR	-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
	Amplitude BKS 1.50	GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	19
	Amplitude Mini BKS 1.50	GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	14
	Summit cd BKS 1.50	GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	14
	Summit ecp BKS 1.50	GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	18
	Hoya GP Wide 1.50	GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	18
	Tact 40 BKS 1.50	GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	18
	Tact 60 BKS 1.50	GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	18
1.53 TRIVEX®						
	HOYALUX ARRAY™ 11, 13, 15, 17, VL Phoenix™	GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
	HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL Phoenix™	GR	-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
	HOYALUX ARRAY 2 11, 13, 15, 17, VL Phoenix™	GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
	HOYALUX ARRAY 2 Wrap 11, 13, 15, 17, VL Phoenix™	GR	-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
	Amplitude BKS Phoenix™	GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	19
	Amplitude Mini BKS Phoenix™	GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
	Summit cd BKS Phoenix™	GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	18
	Summit ecp BKS Phoenix™	GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	18
	Hoya GP Wide Phoenix™	GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	18
	Tact 40 BKS Phoenix™	GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	18
	Tact 60 BKS Phoenix™	GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	18

1.59 Polycarbonate						
HOYALUX ARRAY™ 11, 13, 15, 17, VL Poly	GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL Poly	GR	-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	
HOYALUX ARRAY 2 11, 13, 15, 17, VL Poly	GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	
HOYALUX ARRAY 2 Wrap 11, 13, 15, 17, VL Poly	GR	-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	
Amplitude BKS Poly	GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	19	
Amplitude Mini BKS Poly	GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	14	
Summit cd BKS Poly	GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	14	
Summit ecp BKS Poly	GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	18	
Hoya GP Wide Poly	GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	18	
Tact 40 BKS Poly	GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	18	
Tact 60 BKS Poly	GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	18	
1.67 Super High Index						
HOYALUX ARRAY 11, 13,15,17,VL 1.67	GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	14, 17, 19, 21	
HOYALUX ARRAY Wrap 11,13,15,17,VL 1.67	GR	-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	
HOYALUX ARRAY 2 11, 13,15,17,VL 1.67	GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	14, 17, 19, 21	
HOYALUX ARRAY 2 Wrap 11,13,15,17,VL 1.67	GR	-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21	
Amplitude BKS 1.67	GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	19	
Amplitude Mini BKS 1.67	GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	14	
Summit cd BKS 1.67	GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	14	
Summit ecp BKS 1.67	GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	18	
Hoya GP Wide 1.67	GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	18	
Tact 40 BKS 1.67	GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	18	
Tact 60 BKS 1.67	GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	18	

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.