




Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Transitions® XTRActive® Polarized™](#) / [Signet Armorlite](#)

Signet Armorlite – Transitions® XTRActive® Polarized™

						
Progressives			Sphere	Add	Cyl to	MFH
1.50 Standard Index						
	Unity Via II	GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity Via Plus II	GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity Via Elite II	GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity Via Mobile II	GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity Via Wrap II	GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12
	Ethos	GR	-9.00 to +4.50	0.75 to +3.50	-4.00	14
	Ethos Plus	GR	-9.00 to +4.50	0.75 to +3.50	-4.00	14
1.59 Polycarbonate						
	Unity Via II	GR	-10.00 to +6.50	+0.50 to +4.50	-6.00	12
	Unity Via Plus II	GR	-10.00 to +6.50	+0.50 to +4.50	-6.00	12
	Unity Via Elite II	GR	-10.00 to +6.50	+0.50 to +4.50	-6.00	12
	Unity Via Mobile II	GR	-10.00 to +6.50	+0.50 to +4.50	-6.00	12
	Unity Via Wrap II	GR	-6.50 to +6.00	+0.50 to +4.50	-6.00	12
	Ethos	GR	-9.00 to +6.00	0.75 to +3.50	-4.00	14
	Ethos Plus	GR	-9.00 to +6.00	0.75 to +3.50	-4.00	14
1.67 Super High Index						
	Unity Via II	GR	-12.00 to +8.25	+0.50 to +4.50	-6.00	12
	Unity Via Plus II	GR	-12.00 to +8.25	+0.50 to +4.50	-6.00	12
	Unity Via Elite II	GR	-12.00 to +8.25	+0.50 to +4.50	-6.00	12
	Unity Via Mobile II	GR	-12.00 to +8.25	+0.50 to +4.50	-6.00	12
	Unity Via Wrap II	GR	-6.00 to +5.00	+0.50 to +4.50	-6.00	12
	Ethos	GR	-12.00 to +8.00	0.75 to +3.50	-4.00	14
	Ethos Plus	GR	-12.00 to +8.00	0.75 to +3.50	-4.00	14

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.