

# Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Transitions® GEN S™](#) / **Unity Lenses**

## Unity Lenses – Transitions® GEN S™

unity		Transitions Gen#S									
Semi-Finished Single Vision										Sphere	Cyl to
<b>1.50 Standard Index</b>											
	Unity SVx	GR	BR							-12.25 to +6.00	-6.00
	Unity SVxtra	GR	BR							-12.25 to +6.00	-6.00
	Unity SVxtreme	GR	BR							-6.00 to +6.00	-6.00
<b>1.53 TRIVEX® Single Vision</b>											
	Unity SVx	GR	BR							-14.00 to +6.00	-6.00
	Unity SVxtra	GR	BR							-14.00 to +6.00	-6.00
	Unity SVxtreme	GR	BR							-6.50 to +6.00	-6.00
<b>1.59 Polycarbonate Single Vision</b>											
	Unity SVx	GR	BR							-14.00 to +6.00	-6.00
	Unity SVxtra	GR	BR							-14.00 to +6.00	-6.00
	Unity SVxtreme	GR	BR							-6.50 to +6.00	-6.00
<b>1.60 High Index Single Vision</b>											
	Unity SVx	GR	BR							-14.00 to +6.00	-6.00
	Unity SVxtra	GR	BR							-14.00 to +6.00	-6.00
	Unity SVxtreme	GR	BR							-6.50 to +6.00	-6.00
<b>1.67 Super High Index Single Vision</b>											
	Unity SVx	GR	BR							-16.00 to +8.00	-6.00
	Unity SVxtra	GR	BR							-16.00 to +8.00	-6.00
	Unity SVxtreme	GR	BR							-6.00 to +6.00	-6.00

unity		Transitions Gen#S											
Progressives										Sphere	Add	Cyl to	MFH
<b>1.50 Standard Index</b>													
	Unity Via II	GR	BR							-7.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity Via Plus II	GR	BR							-7.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity Via Elite II	GR	BR							-7.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity Via Mobile II	GR	BR							-7.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity Via Wrap II	GR	BR							-7.00 to +6.00	+0.50 to +4.50	-6.00	12
	Ethos	GR	BR							-9.00 to +4.50	0.75 to +3.50	-4.00	14
	Ethos Plus	GR	BR							-9.00 to +4.50	0.75 to +3.50	-4.00	14
<b>1.53 TRIVEX®</b>													
	Unity Via II	GR	BR							-10.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity Via Plus II	GR	BR							-10.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity Via Elite II	GR	BR							-10.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity Via Mobile II	GR	BR							-10.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity Via Wrap II	GR	BR							-5.00 to +4.00	+0.50 to +4.50	-6.00	12
	Ethos	GR	BR							-9.00 to +6.00	0.75 to +3.50	-4.00	14
	Ethos Plus	GR	BR							-9.00 to +6.00	0.75 to +3.50	-4.00	14
<b>1.59 Polycarbonate</b>													
	Unity Via II	GR	BR							-10.00 to +6.50	+0.50 to +4.50	-6.00	12
	Unity Via Plus II	GR	BR							-10.00 to +6.50	+0.50 to +4.50	-6.00	12
	Unity Via Elite II	GR	BR							-10.00 to +6.50	+0.50 to +4.50	-6.00	12

