















Transitions® Lenses Product Availability Guide














< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Suppliers](#) / **Essilor Lenses**

Essilor Lenses – Entire Product Line

















	Transitions Gen#S								Transitions XTRACTIVE			Transitions POLARIZED		Sphere	Cyl to
	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	BR		
Finished Single Vision															
1.50 Standard Index FSV															
Essilor FSV Spherical	GR	BR							GR					-4.00 to +4.00	-2.00
Essilor FSV Aspheric	GR													-4.00 to +4.00	-2.00
Essilor FSV 1.50 Crizal® Sapphire™ HR	GR								GR					-4.00 to +4.00	-2.00
Essilor FSV 1.50 Crizal® Rock™	GR								GR					-4.00 to +4.00	-2.00
Essilor FSV 1.50 Crizal® Easy Pro	GR								GR					-4.00 to +4.00	-2.00
1.59 Polycarbonate Finished Single Vision															
Essilor FSV Spherical Airwear®	GR	BR												-6.00 to +4.00	-2.00
Essilor FSV Aspheric Airwear®	GR	BR							GR					-6.00 to +4.00	-2.00
Essilor FSV Airwear® Crizal® Sapphire™ HR	GR								GR					-6.00 to +4.00	-2.00
Essilor FSV Airwear® Crizal® Rock™	GR								GR					-6.00 to +4.00	-2.00
Essilor FSV Airwear® Crizal® Easy Pro	GR								GR					-6.00 to +4.00	-2.00

	Transitions Gen#S								Transitions XTRACTIVE			Transitions POLARIZED		Sphere	Cyl to
	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	BR		
Semi-Finished Single Vision															
1.50 Standard Index SFSV															
Essilor Spherical Single Vision	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.50 to +10.00	-4.00
Essilor Stylistic® Single Vision	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	-6.00
Eyezen® +0	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	-6.00
Eyezen® +1	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	-6.00
Eyezen® +2	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	-6.00
Eyezen® +3	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	-6.00
Eyezen® +4	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	-6.00
Eyezen® Start	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +6.00	-6.00
1.53 TRIVEX® Single Vision															
Essilor Spherical Single Vision TREXA™	GR	BR	GG						GR	BR				-10.00 to +6.00	-4.00
Essilor Eyecode™ Single Vision TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	-6.00
Essilor FIT™ Single Vision TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	-6.00
Essilor Azio™ Single Vision TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	-6.00
Essilor 360°™ Single Vision TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	-6.00
Eyezen® +0 TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	-6.00
Eyezen® +1 TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	-6.00
Eyezen® +2 TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	-6.00
Eyezen® +3 TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	-6.00
Eyezen® +4 TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	-6.00
Eyezen® Start	GR	BR	GG						GR	BR				-10.00 to +6.00	-6.00
1.59 Polycarbonate Single Vision															
Essilor Spherical Single Vision Airwear®	GR	BR	GG	S	AME	AMB	E	R			GG	GR		-10.00 to	-4.00






	Transitions Gen#S							Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED					
Progressives													Sphere	Add	Cyl to	MFH
1.53 TRIVEX®																
Varilux® XR Track™	GR	BR											-10.00 to +6.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR											-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™ TREXA™	GR	BR											-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ TREXA™	GR	BR											-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ TREXA™	GR	BR											-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode TREXA™	GR	BR						GR	BR				-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT TREXA™	GR	BR						GR	BR				-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ TREXA™	GR	BR						GR	BR				-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit TREXA™	GR	BR						GR	BR				-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max TREXA™	GR	BR						GR	BR				-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ TREXA™	GR	BR	GG					GR	BR				-10.00 to +7.50	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ TREXA™	GR	BR	GG					GR	BR				-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™ TREXA™	GR	BR	GG					GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ TREXA™	GR	BR	GG					GR	BR				-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Essilor Ideal Advanced™ FIT TREXA™	GR	BR	GG										-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ TREXA™	GR	BR	GG					GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT TREXA™	GR	BR	GG					GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport TREXA™	GR	BR	GG					GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	18
Essilor Ideal® Computer TREXA™	GR	BR	GG					GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	13
Essilor Ideal® Driving TREXA™	GR	BR	GG					GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Ideal® TREXA™	GR	BR	GG					GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Ideal® Short TREXA™	GR	BR	GG					GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Natural® Digital TREXA™	GR	BR	GG					GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ TREXA™	GR	BR	GG					GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	15
Essilor Accolade® TREXA™	GR	BR	GG					GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Digital TREXA™	GR	BR	GG					GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital TREXA™	GR	BR	GG					GR	BR				-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Brillenz® TREXA™	GR	BR	GG					GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Brillenz® Short TREXA™	GR	BR	GG					GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Ovation® Digital TREXA™	GR	BR	GG					GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor SmallFit® Digital TREXA™	GR	BR	GG					GR	BR				-10.00 to +7.50	+0.75 to +3.00	-6.00	14

	Transitions Gen#S							Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED					
Progressives													Sphere	Add	Cyl to	MFH
1.59 Polycarbonate																
Varilux® XR Track™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14

Varilux® X Design™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Varilux Comfort DRx™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux® Liberty™ 3.0	GR	BR												-10.00 to +8.00	+0.75 to +1.50	-6.00	14
Varilux Comfort® Airwear®	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short™ Airwear®	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +3.00	-6.00	14
Varilux® Stylistic™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ FIT	GR	BR	GG	S	AME	AMB	E	R						-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Wrap	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Essilor Ideal® FIT Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Sport Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Driving Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor Ideal® Short Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor Natural® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +8.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Ovation® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Ovation® Airwear®	GR	BR												-10.00 to +8.00	+0.75 to +3.50	-4.00	17
Essilor SmallFit® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +8.00	+0.75 to +3.00	-6.00	14

																	
Progressives													Sphere	Add	Cyl to	MFH	
1.60 High Index																	
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14	
Varilux® Physio® W3+ FIT Thin&Lite® 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14	
Varilux® Physio® W3+ Thin&Lite® 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14	
Varilux Comfort® Max Fit Thin&Lite® 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50	-6.00	14	
Varilux Comfort® Max Thin&Lite® 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50	-6.00	14	
Varilux Physio DRx™ Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17	

Essilor Ideal Advanced Wrap™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Sport Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Driving Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Essilor Ideal® Short	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Essilor Natural® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-14.00 to +9.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Brillenz® Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Ovation® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Ovation® Thin&Lite® 1.67	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-4.00	17
Essilor SmallFit® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-14.00 to +9.00	+0.75 to +3.00	-6.00	14

	  										Sphere	Add	Cyl to	MFH	
															
Progressives															
1.74 Ultra High Index															
Varilux® XR Track™	GR	BR										-14.00 to +13.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR										-14.00 to +13.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™ Thin&Lite® 1.74	GR	BR										-14.00 to +13.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ Thin&Lite® 1.74	GR	BR										-14.00 to +13.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ Thin&Lite® 1.74	GR	BR										-14.00 to +13.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.74	GR	BR										-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Thin&Lite® 1.74	GR	BR										-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Thin&Lite® 1.74	GR	BR										-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio DRx™ Thin&Lite® 1.74	GR	BR										-14.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Thin&Lite® 1.74	GR	BR										-14.00 to +13.00	+0.75 to +3.50	-6.00	14
Varilux Comfort DRx™ Thin&Lite® 1.74	GR	BR										-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Thin&Lite® 1.74	GR	BR										-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ FIT Thin&Lite® 1.74	GR	BR										-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Thin&Lite® 1.74	GR	BR										-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.74	GR	BR										-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport Thin&Lite® 1.74	GR	BR										-20.00 to +13.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Computer Thin&Lite® 1.74	GR	BR										-20.00 to +13.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Driving Thin&Lite® 1.74	GR	BR										-20.00 to +13.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Thin&Lite® 1.74	GR	BR										-20.00 to +13.00	+0.75 to +3.50	-6.00	13
Ovation® Digital Thin&Lite® 1.74	GR	BR										-20.00 to +13.00	+0.75 to +3.50	-6.00	14
Essilor Accolade Freedom™ Thin&Lite® 1.74	GR	BR										-20.00 to +13.00	+0.75 to +3.50	-6.00	15

TRANSITIONS® LIGHT INTELLIGENT LENSES™**GR**
Gray**BR**
Brown**GG**
Graphite
Green**S**
Sapphire**AME**
Amethyst**AMB**
Amber**E**
Emerald**R**
RubyAvailable
in style
mirrors

Transitions, *Transitions Signature*, and *XTRActive* are registered trademarks, and *XTRActive Polarized*, the *Transitions* logo and *Transitions Light Intelligent Lenses* are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

Unless indicated otherwise, all registered trademarks and trademarks are the property of Essilor International and/or its subsidiaries in the United States and in other countries.

All other trademarks are the property of their respective owners.