




Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Suppliers](#) / [Signet Armorlite](#)

Signet Armorlite – Entire Product Line

	Transitions Gen#S								Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED			Sphere	Cyl to
	GR	BR	GG	S	AME	AMB	E	R	GR	GR	GR	GR	GR	GR		
Semi-Finished Single Vision																
1.50 Standard Index SFSV																
Signetek Processed SV	GR	BR	GG	S	AME	AMB	E	R	GR						-8.00 to +6.00	-6.00
KODAK Digital SV	GR	BR	GG	S	AME	AMB	E	R	GR						-8.00 to +6.00	-6.00
KODAK PowerUp Lens (.40/.66/1.00D)	GR	BR													-8.00 to +6.00	-6.00
Crossbows Custom SV	GR	BR													-8.00 to +6.00	-6.00
1.53 TRIVEX® Single Vision																
Signetek Processed SV	GR	BR	GG	S	AME	AMB	E	R	GR						-10.00 to +6.00	-6.00
KODAK Digital SV	GR	BR	GG	S	AME	AMB	E	R	GR						-10.00 to +6.00	-6.00
Crossbows Custom SV	GR	BR													-10.00 to +6.00	-6.00
1.59 Polycarbonate Single Vision																
Signetek Processed SV	GR	BR	GG	S	AME	AMB	E	R	GR		GR				-10.00 to +6.50	-6.00
KODAK Digital SV	GR	BR	GG	S	AME	AMB	E	R	GR		GR				-10.00 to +6.50	-6.00
KODAK PowerUp Lens (.40/.66/1.00D)	GR	BR													-10.00 to +6.50	-6.00
Crossbows Custom SV	GR	BR													-10.00 to +6.50	-6.00
1.60 High Index Single Vision																
Signetek Processed SV	GR	BR													-10.00 to +7.50	-6.00
KODAK Digital SV	GR	BR													-10.00 to +7.50	-6.00
Crossbows Custom SV	GR	BR													-10.00 to +7.50	-6.00
1.67 Super High Index Single Vision																
Signetek Processed SV	GR	BR	GG	S	AME	AMB	E	R	GR						-11.00 to +7.50	-6.00
KODAK Digital SV	GR	BR	GG	S	AME	AMB	E	R	GR						-11.00 to +7.50	-6.00
KODAK PowerUp Lens (.40/.66/1.00D)	GR	BR													-11.00 to +6.00	-6.00
Crossbows Custom SV	GR	BR													-11.00 to +7.50	-6.00
1.74 Ultra High Index Single Vision																
Signetek™ Processed SV	GR	BR													-13.00 to +9.00	-6.00
KODAK Digital SV	GR	BR													-13.00 to +9.00	-6.00

	Transitions Gen#S								Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED			Sphere	Add	Cyl to
	GR	BR	GG	S	AME	AMB	E	R	GR	GR	GR	GR	GR				
Bifocal Signature™																	
1.50 Standard Index SFST-28																	
Signetek™ Processed 1.50 FT-28									GR						-10.00 to +6.00	+1.00 to +3.00	-6.00

	Transitions Gen#S								Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED			Sphere	Add	Cyl to	MFH
	GR	BR	GG	S	AME	AMB	E	R	GR	GR	GR	GR	GR					
Progressives																		
1.50 Standard Index																		
KODAK Precise® 1.50 SF	GR														-9.00 to +5.00	+1.00 to +3.00	-4.00	17
KODAK Precise® Digital 1.50	GR														-8.00 to +5.00	+1.00 to +3.00	-6.00	17
KODAK Precise® Plus 1.50	GR	BR	GG	S	AME	AMB	E	R	GR					-8.00 to +4.00	+0.75 to +3.50	-6.00	17	
KODAK Precise® Short Plus 1.50	GR	BR	GG	S	AME	AMB	E	R	GR					-8.00 to +4.00	+0.75 to +3.50	-6.00	13	
KODAK Precise® PB 1.50	GR	BR							GR					-8.00 to +4.00	+0.75 to +3.50	-6.00	17	
KODAK Precise® Short PB 1.50	GR	BR							GR					-8.00 to +4.00	+0.75 to +3.50	-6.00	13	

