


















Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Design](#) / [Essilor Lenses](#) / **Progressives**

Essilor Lenses – Progressives







													Sphere	Add	Cyl to	MFH	
	Progressives																
1.50 Standard Index																	
Varilux® XR Track™	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Varilux® Liberty™ 3.0	GR	BR												-10.00 to +7.50	+0.75 to +1.50	-6.00	14
Varilux Comfort®	GR	BR							GR	BR	GG			-10.00 to +7.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short™	GR	BR							GR	BR	GG			-10.00 to +7.00	+0.75 to +3.00	-6.00	14
Varilux® Stylistic™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ FIT	GR	BR	GG	S	AME	AMB	E	R						-10.00 to +7.50	+0.75 to +3.50	-4.00	17
Essilor Ideal Advanced™ Wrap	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-4.00	17
Essilor Ideal® FIT	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	18
Essilor Ideal® Driving	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	13
Essilor Ideal®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Ideal® Short	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Natural® Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.50	-6.00	15
Essilor Accolade™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Brillenz®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Brillenz® Short	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14

Ovation® Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Ovation®	GR	BR												-10.00 to +6.00	+0.75 to +3.50	-4.00	17
Essilor SmallFit® Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.00	-6.00	14















	Transitions Gen#S								Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED					
Progressives													Sphere	Add	Cyl to	MFH	
1.53 TRIVEX®																	
Varilux® XR Track™	GR	BR											-10.00 to +6.00	+0.75 to +4.00	-6.00	NVB Measurement	
Varilux® XR Design™	GR	BR											-10.00 to +6.00	+0.75 to +4.00	-6.00	14	
Varilux® X 4D™ TREXA™	GR	BR											-10.00 to +6.00	+0.75 to +4.00	-6.00	14	
Varilux® X FIT™ TREXA™	GR	BR											-10.00 to +6.00	+0.75 to +4.00	-6.00	14	
Varilux® X Design™ TREXA™	GR	BR											-10.00 to +6.00	+0.75 to +4.00	-6.00	14	
Varilux® Physio® W3+ Eyecode TREXA™	GR	BR							GR	BR			-10.00 to +6.00	+0.75 to +3.50 to +4.00	-6.00	14	
Varilux® Physio® W3+ FIT TREXA™	GR	BR							GR	BR			-10.00 to +6.00	+0.75 to +3.50 to +4.00	-6.00	14	
Varilux® Physio® W3+ TREXA™	GR	BR							GR	BR			-10.00 to +6.00	+0.75 to +3.50 to +4.00	-6.00	14	
Varilux Comfort® Max Fit TREXA™	GR	BR							GR	BR			-10.00 to +6.00	+0.75 to +3.50	-6.00	14	
Varilux Comfort® Max TREXA™	GR	BR							GR	BR			-10.00 to +6.00	+0.75 to +3.50	-6.00	14	
Varilux Physio DRx™ TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50 to +4.00	-6.00	17	
Varilux Physio Short DRx™ TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.00	-6.00	14	
Varilux Comfort DRx™ TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	17	
Varilux Comfort Short DRx™ TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.00	-6.00	14	
Essilor Ideal Advanced™ FIT TREXA™	GR	BR	GG										-10.00 to +7.50	+0.75 to +3.50	-6.00	17	
Essilor Ideal Advanced™ TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	17	
Essilor Ideal® FIT TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	17	
Essilor Ideal® Sport TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	18	
Essilor Ideal® Computer TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	13	
Essilor Ideal® Driving TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	14	
Essilor Ideal® TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	14	
Essilor Ideal® Short TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	14	
Essilor Natural® Digital TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	17	
Essilor Accolade Freedom™ TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	15	
Essilor Accolade® TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	17	
Adaptar® Digital TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	17	
Adaptar® Short Digital TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.00	-6.00	14	
Brillenz® TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	14	
Brillenz® Short TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	14	
Ovation® Digital TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	17	
Essilor SmallFit® Digital TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.00	-6.00	14	

	Transitions Gen#S								Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED					
Progressives													Sphere	Add	Cyl to	MFH	
1.59 Polycarbonate																	
Varilux® XR Track™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	NVB Measurement	





Varilux® XR Design™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Varilux Comfort DRx™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux® Liberty™ 3.0	GR	BR												-10.00 to +8.00	+0.75 to +1.50	-6.00	14
Varilux Comfort® Airwear®	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short™ Airwear®	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +3.00	-6.00	14
Varilux® Stylistic™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ FIT	GR	BR	GG	S	AME	AMB	E	R						-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Wrap	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Essilor Ideal® FIT Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Sport Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Driving Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor Ideal® Short Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor Natural® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +8.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Ovation® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Ovation® Airwear®	GR	BR												-10.00 to +8.00	+0.75 to +3.50	-4.00	17
Essilor SmallFit® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +8.00	+0.75 to +3.00	-6.00	14

																	
Progressives													Sphere	Add	Cyl to	MFH	
1.60 High Index																	
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14	
Varilux® Physio® W3+ FIT Thin&Lite® 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14	
Varilux® Physio® W3+ Thin&Lite® 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14	

Varilux Comfort® Max Fit Thin&Lite® 1.60	GR	BR	GG														-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max Thin&Lite® 1.60	GR	BR	GG														-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Thin&Lite® 1.60	GR	BR															-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Thin&Lite® 1.60	GR	BR															-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™ Thin&Lite® 1.60	GR	BR															-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Thin&Lite® 1.60	GR	BR															-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux Comfort® Thin&Lite® 1.60	GR																-12.00 to +8.00	+0.75 to +3.50	-4.00	17
Varilux Comfort Short™ Thin&Lite® 1.60	GR																-12.00 to +8.00	+0.75 to +3.00	-4.00	14
Essilor Ideal Advanced™ FIT Thin&Lite® 1.60	GR	BR															-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Thin&Lite® 1.60	GR	BR															-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.60	GR	BR															-12.00 to +8.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Sport Thin&Lite® 1.60	GR	BR															-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Thin&Lite® 1.60	GR	BR															-12.00 to +8.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Driving Thin&Lite® 1.60	GR	BR															-12.00 to +8.00	+0.75 to +3.50	-6.00	13
Essilor Ideal® Thin&Lite® 1.60	GR	BR															-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Short Thin&Lite® 1.60	GR	BR															-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Essilor Natural® Digital Thin&Lite® 1.60	GR	BR															-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Thin&Lite® 1.60	GR	BR															-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Thin&Lite® 1.60	GR	BR															-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Ovation® Digital Thin&Lite® 1.60	GR	BR															-12.00 to +8.00	+0.75 to +3.50	-6.00	17
SmallFit® Digital Thin&Lite® 1.60	GR	BR															-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Thin&Lite® 1.60	GR	BR															-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Thin&Lite® 1.60	GR	BR															-12.00 to +8.00	+0.75 to +3.00	-6.00	14

 essilor	Transitions® Gen#S							Transitions® XTRACTIVE			Transitions® XTRACTIVE POLARIZED										
Progressives														Sphere	Add	Cyl to	MFH				
1.67 Super High Index																					
Varilux® XR Track™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-12.00 to +9.00	+0.75 to +4.00	-6.00				NVB Measurement	
Varilux® XR Design™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-12.00 to +9.00	+0.75 to +4.00	-6.00				14	
Varilux® X 4D™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-12.00 to +9.00	+0.75 to +4.00	-6.00				14	
Varilux® X FIT™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-12.00 to +9.00	+0.75 to +4.00	-6.00				14	
Varilux® X Design™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-12.00 to +9.00	+0.75 to +4.00	-6.00				14	
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00				14	
Varilux® Physio® W3+ FIT Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00				14	
Varilux® Physio® W3+ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00				14	
Varilux Comfort® Max Fit Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-14.00 to +9.00	+0.75 to +3.50	-6.00				14	
Varilux Comfort® Max Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-14.00 to +9.00	+0.75 to +3.50	-6.00				14	
Varilux Physio DRx™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00				17	
Varilux Physio Short DRx™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.00	-6.00				14	
Varilux Comfort DRx™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.50	-6.00				17	
Varilux Comfort Short DRx™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	+0.75 to +3.00	-6.00				14	
Varilux® Liberty™ 3.0	GR	BR												-14.00 to +9.00	+0.75 to +1.50	-6.00				14	
Varilux Comfort® Thin&Lite® 1.67	GR	BR												-14.00 to +9.00	+0.75 to +3.00	-6.00				14	

Varilux Comfort® Short™ Thin&Lite® 1.67	GR	BR																	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Varilux® Stylistic™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR							-5.00 to +5.00	+0.75 to +3.50	-3.00	17
Essilor Ideal Advanced™ FIT Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R											-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal Advanced Wrap™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR							-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR							-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR							-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Sport Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR							-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR							-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Driving Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR							-14.00 to +9.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR							-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Essilor Ideal® Short	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR							-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Essilor Natural® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR							-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG								-14.00 to +9.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG								-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG								-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG								-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Brillenz® Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR							-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR							-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Ovation® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG								-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Ovation® Thin&Lite® 1.67	GR	BR																	-12.00 to +8.00	+0.75 to +3.50	-4.00	17
Essilor SmallFit® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG								-14.00 to +9.00	+0.75 to +3.00	-6.00	14

				Transitions Gen#S		Transitions XTRACTIVE		Transitions XTRACTIVE PGL-ARIZEP														
Progressives										Sphere		Add		Cyl to		MFH						
1.74 Ultra High Index																						
Varilux® XR Track™	GR	BR																	-14.00 to +13.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR																	-14.00 to +13.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™ Thin&Lite® 1.74	GR	BR																	-14.00 to +13.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ Thin&Lite® 1.74	GR	BR																	-14.00 to +13.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ Thin&Lite® 1.74	GR	BR																	-14.00 to +13.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.74	GR	BR																	-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Thin&Lite® 1.74	GR	BR																	-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Thin&Lite® 1.74	GR	BR																	-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio DRx™ Thin&Lite® 1.74	GR	BR																	-14.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Thin&Lite® 1.74	GR	BR																	-14.00 to +13.00	+0.75 to +3.50	-6.00	14
Varilux Comfort DRx™ Thin&Lite® 1.74	GR	BR																	-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Thin&Lite® 1.74	GR	BR																	-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ FIT Thin&Lite® 1.74	GR	BR																	-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Thin&Lite® 1.74	GR	BR																	-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.74	GR	BR																	-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport Thin&Lite® 1.74	GR	BR																	-20.00 to +13.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Computer Thin&Lite® 1.74	GR	BR																	-20.00 to +13.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Driving Thin&Lite® 1.74	GR	BR																	-20.00 to +13.00	+0.75 to +3.50	-6.00	18

Essilor Ideal® Thin&Lite® 1.74	GR	BR																		-20.00 to +13.00	+0.75 to +3.50	-6.00	13
Ovation® Digital Thin&Lite® 1.74	GR	BR																		-20.00 to +13.00	+0.75 to +3.50	-6.00	14
Essilor Accolade Freedom™ Thin&Lite® 1.74	GR	BR																		-20.00 to +13.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Thin&Lite® 1.74	GR	BR																		-14.00 to +9.00	+0.75 to +3.50	-6.00	17

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

Unless indicated otherwise, all registered trademarks and trademarks are the property of Essilor International and/or its subsidiaries in the United States and in other countries.

All other trademarks are the property of their respective owners.