

Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Design](#) / [Progressives](#) / [Hoya Vision Care](#)

Hoya Vision Care — Progressives

HOYA	Transitions® Gen#S								Transitions® XTRACTIVE®			Transitions® XTRACTIVE® POLARIZED®			Sphere	Add	Cyl to	MFH	
	Progressives																		
1.50 Standard Index																			
HOYALUX iD LifeStyle®4 Indoor 1.50	GR!	BR!														-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Outdoor 1.50	GR!	BR!														-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Urban 1.50	GR!	BR!														-8.00 to +6.00	+0.75 to +4.00	-6.00	14
Amp HD3 1.50	GR!	BR!														-8.00 to +6.00	+0.75 to +4.00	-6.00	18
Amp Mini HD3 1.50	GR!	BR!														-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX ARRAY™ 11, 13, 15, 17, VL 1.50	GR	BR							GR	BR		GR				-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL 1.50	GR	BR							GR	BR		GR				-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ 2 11, 13, 15, 17, VL 1.50	GR	BR							GR	BR		GR				-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ 2 Wrap 11, 13, 15, 17, VL 1.50	GR	BR							GR	BR		GR				-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS 1.50	GR	BR							GR	BR		GR				-10.00 to +3.75	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS 1.50	GR	BR							GR	BR		GR				-10.00 to +3.75	+0.75 to +4.00	-6.00	14
Summit cd BKS 1.50	GR	BR							GR	BR		GR				-10.00 to +3.75	+0.75 to +4.00	-6.00	14
Summit ecp BKS 1.50	GR	BR							GR	BR		GR				-10.00 to +3.75	+0.75 to +4.00	-6.00	18
Hoya GP Wide 1.50	GR	BR							GR	BR		GR				-10.00 to +3.75	+0.75 to +4.00	-6.00	18
Tact 40 BKS 1.50	GR	BR							GR	BR		GR				-10.00 to +3.75	+0.75 to +4.00	-6.00	18
Tact 60 BKS 1.50	GR	BR							GR	BR		GR				-10.00 to +3.75	+0.75 to +4.00	-6.00	18
1.53 TRIVEX®																			
HOYALUX iD LifeStyle®4 Indoor Phoenix™	GR!	BR!														-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Outdoor Phoenix™	GR!	BR!														-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Urban Phoenix™	GR!	BR!														-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Indoor Phoenix™	GR!	BR!														-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Outdoor Phoenix™	GR!	BR!														-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Urban Phoenix™	GR!	BR!														-8.00 to +6.00	+0.75 to +4.00	-6.00	14
Amp HD3 Phoenix™	GR!	BR!														-8.00 to +6.00	+0.75 to +4.00	-6.00	18
Amp Mini HD3 Phoenix™	GR!	BR!														-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX ARRAY™ 11, 13, 15, 17, VL Phoenix™	GR	BR							GR	BR		GR				-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL Phoenix™	GR	BR							GR	BR		GR				-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13, 15, 17, VL Phoenix™	GR	BR							GR	BR		GR				-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11, 13, 15, 17, VL Phoenix™	GR	BR							GR	BR		GR				-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS Phoenix™	GR	BR							GR	BR		GR				-10.00 to +6.00	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS Phoenix™	GR	BR							GR	BR		GR				-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Summit cd BKS Phoenix™	GR	BR							GR	BR		GR				-10.00 to +6.00	+0.75 to +4.00	-6.00	18
Summit ecp BKS Phoenix™	GR	BR							GR	BR		GR				-10.00 to +6.00	+0.75 to +4.00	-6.00	18
Hoya GP Wide Phoenix™	GR	BR							GR	BR		GR				-10.00 to +6.00	+0.75 to +4.00	-6.00	18
Tact 40 BKS Phoenix™	GR	BR							GR	BR		GR				-10.00 to +6.00	+0.75 to +4.00	-6.00	18
Tact 60 BKS Phoenix™	GR	BR							GR	BR		GR				-10.00 to +6.00	+0.75 to +4.00	-6.00	18
1.59 Polycarbonate																			
HOYALUX iD LifeStyle®4 Indoor Poly	GR!	BR!														-6.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Outdoor Poly	GR!	BR!														-6.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Urban Poly	GR!	BR!														-6.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Indoor Poly	GR!	BR!														-6.00 to +6.00	+0.75 to +4.00	-6.00	14

HOYALUX iD LifeStyle®3 Outdoor Poly	GR!	BR!										-6.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Urban Poly	GR!	BR!										-6.00 to +6.00	+0.75 to +4.00	-6.00	14
Amp HD3 Poly	GR!	BR!										-6.00 to +6.00	+0.75 to +4.00	-6.00	18
Amp Mini HD3 Poly	GR!	BR!										-6.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX ARRAY™ 11, 13, 15, 17, VL Poly	GR	BR					GR	BR		GR		-13.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL Poly	GR	BR					GR	BR		GR		-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13, 15, 17, VL Poly	GR	BR					GR	BR		GR		-13.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11, 13, 15, 17, VL Poly	GR	BR					GR	BR		GR		-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS Poly	GR	BR					GR	BR		GR		-13.00 to +6.00	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS Poly	GR	BR					GR	BR		GR		-13.00 to +6.00	+0.75 to +4.00	-6.00	14
Summit cd BKS Poly	GR	BR					GR	BR		GR		-13.00 to +6.00	+0.75 to +4.00	-6.00	14
Summit ecp BKS Poly	GR	BR					GR	BR		GR		-13.00 to +6.00	+0.75 to +4.00	-6.00	18
Hoya GP Wide Poly	GR	BR					GR	BR		GR		-13.00 to +6.00	+0.75 to +4.00	-6.00	18
Tact 40 BKS Poly	GR	BR					GR	BR		GR		-13.00 to +6.00	+0.75 to +4.00	-6.00	18
Tact 60 BKS Poly	GR	BR					GR	BR		GR		-13.00 to +6.00	+0.75 to +4.00	-6.00	18
1.60 High Index															
HOYALUX iD LifeStyle®4 Indoor 1.60	GR!	BR!										-11.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Outdoor 1.60	GR!	BR!										-11.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Urban 1.60	GR!	BR!										-11.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Indoor 1.60	GR!	BR!										-11.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Outdoor 1.60	GR!	BR!										-11.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Urban 1.60	GR!	BR!										-11.00 to +8.00	+0.75 to +4.00	-6.00	14
Amp HD3 1.60	GR!	BR!										-11.00 to +8.00	+0.75 to +4.00	-6.00	14
Amp Mini HD3 1.60	GR!	BR!										-11.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX ARRAY 11, 13, 15, 17, VL 1.60	GR	BR										-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY Wrap 11,13,15,17,VL 1.60	GR	BR										-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13, 15, 17, VL 1.60	GR	BR										-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11,13,15,17,VL 1.60	GR	BR										-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS 1.60	GR	BR										-11.00 to +8.00	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS 1.60	GR	BR										-11.00 to +8.00	+0.75 to +4.00	-6.00	14
Summit cd BKS 1.60	GR	BR										-11.00 to +8.00	+0.75 to +4.00	-6.00	14
Summit ecp BKS 1.60	GR	BR										-11.00 to +8.00	+0.75 to +4.00	-6.00	18
Hoya GP Wide 1.60	GR	BR										-11.00 to +8.00	+0.75 to +4.00	-6.00	18
Tact 40 BKS 1.60	GR	BR										-11.00 to +8.00	+0.75 to +4.00	-6.00	18
Tact 60 BKS 1.60	GR	BR										-11.00 to +8.00	+0.75 to +4.00	-6.00	18
1.67 Super High Index															
HOYALUX iD LifeStyle®4 Indoor 1.67	GR!	BR!										-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Outdoor 1.67	GR!	BR!										-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Urban 1.67	GR!	BR!										-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle 3 Indoor 1.67	GR!	BR!										-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle 3 Outdoor 1.67	GR!	BR!										-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle 3 Urban 1.67	GR!	BR!										-13.00 to +8.00	+0.75 to +4.00	-6.00	14
Amp HD3 1.67	GR!	BR!										-13.00 to +8.00	+0.75 to +4.00	-6.00	18
Amp Mini HD3 1.67	GR!	BR!										-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX ARRAY 11, 13,15,17,VL 1.67	GR	BR					GR	BR		GR		-13.00 to +7.50	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY Wrap 11,13,15,17,VL 1.67	GR	BR					GR	BR		GR		-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13,15,17,VL 1.67	GR	BR					GR	BR		GR		-13.00 to +7.50	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11,13,15,17,VL 1.67	GR	BR					GR	BR		GR		-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS 1.67	GR	BR					GR	BR		GR		-13.00 to +7.50	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS 1.67	GR	BR					GR	BR		GR		-13.00 to +7.50	+0.75 to +4.00	-6.00	14
Summit cd BKS 1.67	GR	BR					GR	BR		GR		-13.00 to +7.50	+0.75 to +4.00	-6.00	14
Summit ecp BKS 1.67	GR	BR					GR	BR		GR		-13.00 to +7.50	+0.75 to +4.00	-6.00	18
Hoya GP Wide 1.67	GR	BR					GR	BR		GR		-13.00 to +7.50	+0.75 to +4.00	-6.00	18

Tact 40 BKS 1.67	GR	BR							GR	BR			GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	18
Tact 60 BKS 1.67	GR	BR							GR	BR			GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	18
1.74 ultra High Index																	
HOYALUX iD LifeStyle®4 Indoor 1.74	GR!	BR!												-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX iD LifeStyle®4 Outdoor 1.74	GR!	BR!												-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX iD LifeStyle®4 Urban 1.74	GR!	BR!												-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX iD LifeStyle 3 Indoor 1.74	GR!	BR!												-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX iD LifeStyle 3 Outdoor 1.74	GR!	BR!												-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX iD LifeStyle 3 Urban 1.74	GR!	BR!												-13.00 to +8.00	+0.75 to +3.50	-6.00	14
Amp HD3 1.74	GR!	BR!												-13.00 to +8.00	+0.75 to +3.50	-6.00	14
Amp Mini HD3 1.74	GR!	BR!												-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX ARRAY 11, 13, 15, 17, VL 1.74	GR	BR												-13.00 to +8.00	+0.75 to +3.50	-6.00	14, 17, 19, 21
HOYALUX ARRAY Wrap 11, 13, 15, 17, VL 1.74	GR	BR												-6.00 to +4.00	+0.75 to +3.50	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13, 15, 17, VL 1.74	GR	BR												-13.00 to +8.00	+0.75 to +3.50	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 WRAP 11, 13, 15, 17, VL 1.74	GR	BR												-6.00 to +4.00	+0.75 to +3.50	-6.00	14, 17, 19, 21
Amplitude BKS	GR	BR												-13.00 to +8.00	+0.75 to +3.50	-6.00	19
Amplitude Mini BKS	GR	BR												-13.00 to +8.00	+0.75 to +3.50	-6.00	14
Summit cd BKS 1.74	GR	BR												-13.00 to +8.00	+0.75 to +3.50	-6.00	14
Summit ecp BKS 1.74	GR	BR												-13.00 to +8.00	+0.75 to +3.50	-6.00	19
Hoya GP Wide 1.74	GR	BR												-13.00 to +8.00	+0.75 to +3.50	-6.00	18
Tact 40 BKS 1.74	GR	BR												-13.00 to +8.00	+0.75 to +3.50	-6.00	18
Tact 60 BKS 1.74	GR	BR												-13.00 to +8.00	+0.75 to +3.50	-6.00	18

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.