




Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Design](#) / [Shamir Insights Inc.](#) / **Progressives**

Shamir Insights Inc. — Progressives

												Sphere	Add	Cyl to	MFH						
	Progressives																				
1.50 Standard Index																					
Shamir Autograph Intelligence™	GR	BR										GR					GR	-12.25 to +6.50	+0.75 to +4.00	-6.00	11, 13,15,18
Shamir Autograph® III	GR	BR										GR					GR	-12.25 to +6.00	+0.75 to +4.00	-6.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR										GR					GR	-12.25 to +6.00	+0.75 to +4.00	-6.00	11, 13, 15, 18
Shamir Intouch™	GR	BR										GR					GR	-12.25 to +6.00	+0.75 to +3.50	-6.00	14, 16, 18
Shamir Spectrum+™	GR	BR										GR					GR	-12.25 to +6.00	+0.75 to +3.50	-6.00	14, 16, 18
Shamir Element™	GR	BR										GR					GR	-7.00 to +4.00	+0.75 to +3.50	-6.00	16, 19
Shamir GENESIS HD	GR	BR										GR					GR	-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir Computer™	GR	BR										GR					GR	-7.00 to +4.00	+0.75 to +2.00	-6.00	14 Lower, 10 Upper
Shamir Workspace™	GR	BR										GR					GR	-6.00 to +6.00	+0.75 to +2.50	-6.00	14 Lower, 10 Upper
Shamir Attitude III Sport	GR	BR										GR					GR	-6.50 to +5.25	+0.75 to +3.50	-6.00	18
Shamir Attitude III Fashion	GR	BR										GR					GR	-6.50 to +5.25	+0.75 to +3.50	-6.00	15, 18
1.53 TRIVEX®																					
Shamir Autograph Intelligence™	GR	BR										GR					GR	-13.25 to +7.00	+0.75 to +4.00	-6.50	11, 13, 15, 18
Shamir Autograph® III	GR	BR										GR					GR	-13.25 to +7.00	+0.75 to +4.00	-6.50	11, 13, 15, 18
Shamir Autograph® II+	GR	BR										GR					GR	-13.25 to +6.00	+0.75 to +4.00	-6.50	11, 13, 15, 18
Shamir Intouch™	GR	BR										GR					GR	-13.25 to +6.00	+0.75 to +3.50	-6.50	14, 16, 18
Shamir Golf™	GR	BR										GR					GR	-5.00 to +4.00	+1.00 to +3.50	-6.50	19
Shamir Spectrum+™	GR	BR										GR					GR	-13.25 to +6.00	+0.75 to +3.50	-6.50	14, 16, 18
Shamir Element™	GR	BR										GR					GR	-13.25 to +6.00	+0.75 to +3.50	-6.50	16, 19
Shamir GENESIS HD	GR	BR										GR					GR	-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir Computer™	GR	BR										GR					GR	-6.00 to +8.00	+0.75 to +2.00	-6.50	14 Lower, 10 Upper
Shamir Workspace™	GR	BR										GR					GR	-6.00 to +8.00	+0.75 to +2.50	-6.50	14 Lower, 10 Upper
Shamir Attitude III Sport	GR	BR										GR					GR	-6.50 to +6.00	+0.75 to +3.50	-6.50	18
Shamir Attitude III Fashion	GR	BR										GR					GR	-6.50 to +6.00	+0.75 to +3.50	-6.50	15, 18
1.59 Polycarbonate																					
Shamir Autograph Intelligence™	GR	BR										GR					GR	-14.75 to +7.50	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR										GR					GR	-14.75 to +8.00	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR										GR					GR	-14.75 to +6.50	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Intouch™	GR	BR										GR					GR	-14.75 to +6.00	+0.75 to +3.50	-7.00	15, 18
Shamir Golf™	GR	BR										GR					GR	-6.00 to +5.25	+1.00 to +3.50	-7.00	19
Shamir Spectrum+™	GR	BR										GR					GR	-14.75 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
Shamir Element™	GR	BR										GR					GR	-10.00 to +6.00	+0.75 to +3.50	-7.00	16, 19
Shamir GENESIS HD	GR	BR										GR					GR	-10.00 to +6.00	+0.75 to +3.50	-7.00	16
Shamir Computer™	GR	BR										GR					GR	-6.00 to +8.00	+0.75 to +2.00	-7.00	14 Lower, 10 Upper
Shamir Workspace™	GR	BR										GR					GR	-6.00 to +8.00	+0.75 to +2.50	-7.00	14 Lower, 10 Upper
Shamir Attitude III Sport	GR	BR										GR					GR	-6.50 to +6.00	+0.75 to +3.50	-6.50	18
Shamir Attitude III Fashion	GR	BR										GR					GR	-6.50 to +6.00	+0.75 to +3.50	-6.50	15, 18
1.60 High Index																					
Shamir Autograph Intelligence™	GR	BR										GR						-14.75 to +8.75	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR										GR						-15.00 to +8.00	+0.75 to +3.50	-7.00	15,18
Shamir Autograph® II+	GR	BR										GR						-15.00 to +6.00	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Golf™ SuperLite™ 1.6	GR	BR										GR						-6.00 to +5.25	+1.00 to +3.50	-7.00	19

Shamir Spectrum+™	GR	BR								GR										-15.00 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
Shamir Element™	GR	BR																		-10.00 to +6.00	+0.75 to +3.50	-7.00	16, 19
Shamir GENESIS HD	GR	BR								GR										-6.00 to +6.00	+0.75 to +3.50	-7.00	16
Shamir Computer™	GR	BR								GR										-6.00 to +8.00	+0.75 to +2.00	-7.00	14 Lower, 10 Upper
Shamir Workspace™	GR	BR								GR										-6.00 to +8.00	+0.75 to +2.50	-7.00	14 Lower, 10 Upper
Shamir InTouch™	GR	BR								GR										-15.00 to +6.00	+0.75 to +3.50	-7.00	15, 18
Shamir Attitude III Sport	GR	BR								GR										-6.50 to +6.00	+0.75 to +3.50	-7.00	18
Shamir Attitude III Fashion	GR	BR								GR										-6.50 to +6.00	+0.75 to +3.50	-7.00	15, 18
1.67 Super High Index																							
Shamir Autograph Intelligence™	GR	BR								GR		GR								-16.75 to +10.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR								GR		GR								-16.75 to +9.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR								GR		GR								-16.75 to +8.25	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Spectrum+™	GR	BR								GR		GR								-16.75 to +8.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Element™ SuperLite™ 1.67	GR	BR								GR		GR								-12.50 to +8.00	+0.75 to +3.50	-8.00	16, 19
Shamir GENESIS HD	GR	BR								GR		GR								-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir InTouch™	GR	BR								GR		GR								-16.75 to +8.00	+0.75 to +3.50	-8.00	15, 18
Shamir Attitude III Sport	GR	BR								GR		GR								-6.50 to +6.00	+0.75 to +3.50	-8.00	18
Shamir Attitude III Fashion	GR	BR								GR		GR								-6.50 to +6.00	+0.75 to +3.50	-8.00	15, 18
1.74 Ultra High Index																							
Shamir Autograph Intelligence™	GR	BR																		-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR																		-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR																		-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir GENESIS HD	GR	BR																		-6.00 to +6.00	+0.75 to +4.00	-4.00	16
Shamir Spectrum+™	GR	BR																		-18.75 to +11.50	+0.75 to +3.50	-9.00	14, 16, 18
Shamir InTouch™	GR	BR																		-18.75 to +15.50	+0.75 to +3.50	-9.00	15, 18

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material. All other trademarks are the property of their respective owners.