

Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Design](#) / [Shamir Insights Inc.](#) / **SFSV**

Shamir Insights Inc. – Semi-Finished Single Vision

shamir	Transitions Gen#S								Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED		Sphere	Cyl to	
	GR	BR	GG	S	AME	AMB	E	R	Available in style mirrors	GR	BR	GR				
Semi-Finished Single Vision																
1.50 SFSV																
Shamir Autograph® III Single Vision	GR	BR								GR	BR	GR			-12.00 to +7.25	-6.00
Shamir Relax™	GR	BR								GR	BR	GR			-14.00 to +8.00	-6.00
1.53 SFSV																
Shamir Autograph® III Single Vision	GR	BR								GR	BR	GR			-14.00 to +8.00	-6.50
Shamir Relax™	GR	BR								GR	BR	GR			-14.00 to +8.00	-6.50
1.59 SFSV																
Shamir Autograph® III Single Vision	GR	BR								GR	BR	GR			-14.00 to +8.75	-7.00
Shamir Relax™	GR	BR								GR	BR	GR			-14.75 to +6.00	-7.00
1.60 SFSV																
Shamir Autograph® III Single Vision	GR	BR								GR					-14.00 to +8.75	-7.00
Shamir Relax™	GR	BR								GR					-14.00 to +8.75	-7.00
1.67 SFSV																
Shamir Autograph® III Single Vision	GR	BR								GR	BR	GR			-16.00 to +16.00	-8.00
Shamir Relax™	GR	BR								GR	BR	GR			-16.75 to +16.00	-8.00
1.74 SFSV																
Shamir Autograph® III Single Vision	GR	GR													-20.00 to +13.50	-9.00

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.