





Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Design](#) / **SFSV**





All Suppliers – Semi-Finished Single Vision

 Semi-Finished Single Vision	 Transitions Gen#S							 Transitions XTRACTIVE			 Transitions POLARIZED	Sphere	Cyl to		
	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			GR	
1.50 Standard Index SFSV															
Essilor Spherical Single Vision	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.50 to +10.00	-4.00	
Essilor Stylistic® Single Vision	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	-6.00	
Eyezen® +0	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	-6.00	
Eyezen® +1	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	-6.00	
Eyezen® +2	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	-6.00	
Eyezen® +3	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	-6.00	
Eyezen® +4	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	-6.00	
Eyezen® Start	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +6.00	-6.00	
1.53 TRIVEX® Single Vision															
Essilor Spherical Single Vision TREXA™	GR	BR	GG						GR	BR			-10.00 to +6.00	-4.00	
Essilor Eyecode™ Single Vision TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	-6.00	
Essilor FIT™ Single Vision TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	-6.00	
Essilor Azio™ Single Vision TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	-6.00	
Essilor 360°™ Single Vision TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	-6.00	
Eyezen® +0 TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	-6.00	
Eyezen® +1 TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	-6.00	
Eyezen® +2 TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	-6.00	
Eyezen® +3 TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	-6.00	
Eyezen® +4 TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	-6.00	
Eyezen® Start	GR	BR	GG						GR	BR			-10.00 to +6.00	-6.00	
1.59 Polycarbonate Single Vision															
Essilor Spherical Single Vision Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-4.00	
Essilor Aspheric Single Vision Airwear®	GR	BR	GG						GR	BR	GG		-10.00 to +8.00	-4.00	
Essilor Eyecode™ Aspheric SV Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
Essilor FIT™ Aspheric SV Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
Essilor AZIO™ Aspheric SV Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
Essilor SV 360°™ Aspheric Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
Essilor Stylistic® Single Vision Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
Eyezen® +0 Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
Eyezen® +1 Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
Eyezen® +2 Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
Eyezen® +3 Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
Eyezen® +4 Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
Eyezen® Start Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
Eyezen® Kids	GR	BR	GG	S	AME	AMB	E	R					-10.00 to +6.00	-6.00	
1.60 High Index Single Vision															
Essilor Spherical Single Vision Thin&Lite® 1.60	GR	BR	GG	S	AME	AMB	E	R					-13.00 to +8.00	-4.00	
Eyezen® +0 Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	-6.00	
Eyezen® +1 Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	-6.00	
Eyezen® +2 Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	-6.00	
Eyezen® +3 Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	-6.00	
Eyezen® +4 Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	-6.00	

Eyezen® Start Thin&Lite® 1.60	GR	BR												-12.00 to +8.00	-6.00
1.67 Super High Index Single Vision															
Essilor Spherical Single Vision Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-4.00
Essilor Aspheric Single Vision Thin&Lite® 1.67	GR	BR	GG											-16.00 to +8.00	-4.00
Essilor Eyecode™ SV Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00
Essilor FIT™ SV Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00
Essilor Azio™ SV Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00
Essilor SV 360°™ Spherical Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00
Essilor Stylistic® Single Vision Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00
Eyezen® +0 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00
Eyezen® +1 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00
Eyezen® +2 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00
Eyezen® +3 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00
Eyezen® +4 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00
Eyezen® Start Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-14.00 to +9.00	-6.00
1.74 Ultra High Index Single Vision															
Essilor Spherical Single Vision Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	-4.00
Essilor Eyecode™ SV Thin&Lite® 1.74	GR	BR												-16.00 to +13.00	-6.00
Essilor FIT™ SV Thin&Lite® 1.74	GR	BR												-16.00 to +13.00	-6.00
Essilor 360°™ SV Thin&Lite® 1.74	GR	BR												-16.00 to +13.00	-6.00
Eyezen® +0 Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	-6.00
Eyezen® +1 Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	-6.00
Eyezen® +2 Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	-6.00
Eyezen® +3 Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	-6.00
Eyezen® +4 Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	-6.00
Eyezen® Start Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	-6.00

HOYA		Transitions Gen#S							Transitions XTRACTIVE™			Transitions XTRACTIVE™ POL-WEZE®			
Semi-Finished Single Vision														Sphere	Cyl to
1.50 Standard Index SFSV															
HOYA SFSV 1.50	GR	BR							GR	BR				-8.00 to +9.00	-4.00
HOYA MySV 150	GR	BR							GR	BR		GR		-8.00 to +9.00	-4.00
HOYA SV iQ 1.50	GR	BR							GR	BR		GR		-8.00 to +9.00	-4.00
Sync III	GR	BR							GR	BR		GR		-10.00 to +6.00	-4.00
1.53 TRIVEX® Single Vision															
HOYA SFSV Phoenix™	GR!	BR!							GR	BR				-13.00 to +9.00	-4.00
HOYA MySV Phoenix™	GR!	BR!							GR	BR		GR		-13.00 to +9.00	-4.00
HOYA SV iQ Phoenix™	GR!	BR!							GR	BR				-13.00 to +9.00	-4.00
Sync III	GR!	BR!												-10.00 to +6.00	-4.00
1.59 Polycarbonate Single Vision															
HOYA SFSV Polycarbonate	GR	BR							GR	BR				-6.00 to +6.00	-4.00
HOYA MySV Polycarbonate	GR	BR							GR	BR		GR		-10.00 to +6.00	-4.00
HOYA SV iQ Polycarbonate	GR	BR							GR	BR		GR		-10.00 to +6.00	-4.00
Sync III	GR	BR							GR	BR		GR		-13.00 to +5.25	-4.00
1.60 High Index Single Vision															
HOYA SFSV 1.60	GR!	BR!												-10.00 to +8.00	-4.00
HOYA SV iQ 1.60	GR!	BR!												-10.00 to +8.00	-4.00
Sync III	GR!	BR!												-10.00 to +6.00	-4.00
1.67 Super High Index Single Vision															
HOYA SFSV 1.67	GR	BR							GR	BR				-10.00 to +8.00	-4.00
HOYA MySV 1.67	GR	BR							GR	BR		GR		-10.00 to +8.00	-4.00
HOYA SV iQ 1.67	GR	BR							GR	BR		GR		-10.00 to +8.00	-4.00






	Sync III	GR	BR							GR	BR		GR	-13.00 to +7.50	-4.00	
1.74 ultra High Index Single Vision																
	HOYA SFSV 1.74	GR!	BR!							GR	BR			-13.00 to +10.00	-4.00	
	HOYA SV IQ 1.74	GR!	BR!											-13.00 to +10.00	-4.00	
	Sync III	GR!	BR!											-13.00 to +8.00	-4.00	

																
Semi-Finished Single Vision														Sphere	Cyl to	
1.50 Standard Index SFSV																
	Nikon SeeMax Infinite	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
	Nikon Infinite Single Vision	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
	Nikon Smart Boost	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
1.53 TRIVEX® Single Vision																
	Nikon SeeMax Infinite	GR	BR	GG										-10.00 to +7.00	-6.00	
	Nikon Infinite Single Vision	GR	BR	GG										-10.00 to +7.00	-6.00	
	Nikon Smart Boost	GR	BR	GG										-10.00 to +7.00	-6.00	
1.59 Polycarbonate Single Vision																
	Nikon SeeMax Infinite	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
	Nikon Infinite Single Vision	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
	Nikon Smart Boost	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
1.60 High Index Single Vision																
	Nikon SeeMax Infinite	GR	BR											-12.00 to +8.00	-6.00	
	Nikon Infinite Single Vision	GR	BR											-12.00 to +8.00	-6.00	
	Nikon Smart Boost	GR	BR											-12.00 to +8.00	-6.00	
1.67 Super High Index Single Vision																
	Nikon SeeMax Infinite	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	-6.00	
	Nikon Infinite Single Vision	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	-6.00	
	Nikon Smart Boost	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	-6.00	
1.74 ultra High Index Single Vision																
	Nikon SeeMax Infinite	GR	BR											-16.00 to +10.00	-6.00	
	Nikon Infinite Single Vision	GR	BR											-16.00 to +10.00	-6.00	
	Nikon Smart Boost	GR	BR											-16.00 to +10.00	-6.00	

																
Semi-Finished Single Vision														Sphere	Cyl to	
1.53 TRIVEX® SFSV																
	Seiko PL 1.53/Trivex Transitions®									GR				-10.00 to +7.00	-5.00	
1.67 Super High Index Single Vision																
	Seiko PL 1.67 Transitions®									GR				-12.00 to +10.00	-5.00	

																
Semi-Finished Single Vision														Sphere	Cyl to	
1.50 SFSV																
	Shamir Autograph® III Single Vision	GR	BR							GR	BR		GR	-12.00 to +7.25	-6.00	
	Shamir Relax™	GR	BR							GR	BR		GR	-14.00 to +8.00	-6.00	
1.53 SFSV																
	Shamir Autograph® III Single Vision	GR	BR							GR	BR		GR	-14.00 to +8.00	-6.50	
	Shamir Relax™	GR	BR							GR	BR		GR	-14.00 to +8.00	-6.50	
1.59 SFSV																
	Shamir Autograph® III Single Vision	GR	BR							GR	BR		GR	-14.00 to +8.75	-7.00	
	Shamir Relax™	GR	BR							GR	BR		GR	-14.75 to +6.00	-7.00	







1.60 SFSV												
Shamir Autograph® III Single Vision	GR	BR							GR	-14.00 to +8.75	-7.00	
Shamir Relax™	GR	BR							GR	-14.00 to +8.75	-7.00	
1.67 SFSV												
Shamir Autograph® III Single Vision	GR	BR						GR	BR	GR	-16.00 to +16.00	-8.00
Shamir Relax™	GR	BR						GR	BR	GR	-16.75 to +16.00	-8.00
1.74 SFSV												
Shamir Autograph® III Single Vision	GR	GR									-20.00 to +13.50	-9.00

																				
Semi-Finished Single Vision																			Sphere	Cyl to
1.50 Standard Index SFSV																				
Signetek Processed SV	GR	BR	GG	S	AME	AMB	E	R	GR								-8.00 to +6.00	-6.00		
KODAK Digital SV	GR	BR	GG	S	AME	AMB	E	R	GR								-8.00 to +6.00	-6.00		
KODAK PowerUp Lens (.40/.66/1.00D)	GR	BR															-8.00 to +6.00	-6.00		
Crossbows Custom SV	GR	BR															-8.00 to +6.00	-6.00		
1.53 TRIVEX® Single Vision																				
Signetek Processed SV	GR	BR	GG	S	AME	AMB	E	R	GR								-10.00 to +6.00	-6.00		
KODAK Digital SV	GR	BR	GG	S	AME	AMB	E	R	GR								-10.00 to +6.00	-6.00		
Crossbows Custom SV	GR	BR															-10.00 to +6.00	-6.00		
1.59 Polycarbonate Single Vision																				
Signetek Processed SV	GR	BR	GG	S	AME	AMB	E	R	GR					GR			-10.00 to +6.50	-6.00		
KODAK Digital SV	GR	BR	GG	S	AME	AMB	E	R	GR					GR			-10.00 to +6.50	-6.00		
KODAK PowerUp Lens (.40/.66/1.00D)	GR	BR															-10.00 to +6.50	-6.00		
Crossbows Custom SV	GR	BR															-10.00 to +6.50	-6.00		
1.60 High Index Single Vision																				
Signetek Processed SV	GR	BR															-10.00 to +7.50	-6.00		
KODAK Digital SV	GR	BR															-10.00 to +7.50	-6.00		
Crossbows Custom SV	GR	BR															-10.00 to +7.50	-6.00		
1.67 Super High Index Single Vision																				
Signetek Processed SV	GR	BR	GG	S	AME	AMB	E	R	GR								-11.00 to +7.50	-6.00		
KODAK Digital SV	GR	BR	GG	S	AME	AMB	E	R	GR								-11.00 to +7.50	-6.00		
KODAK PowerUp Lens (.40/.66/1.00D)	GR	BR															-11.00 to +6.00	-6.00		
Crossbows Custom SV	GR	BR															-11.00 to +7.50	-6.00		
1.74 Ultra High Index Single Vision																				
Signetek™ Processed SV	GR	BR															-13.00 to +9.00	-6.00		
KODAK Digital SV	GR	BR															-13.00 to +9.00	-6.00		

																				
Semi-Finished Single Vision																			Sphere	Cyl to
1.50 Standard Index SFSV																				
TRN 1.50 SV (72mm) UC/HC DIA.72									GR	BR							Base 0.50, 2.25, 4.25, 5.25, 6.25, 8.25			
TRN 1.50 DS (76mm) UC/HC DIA.72	GR	BR							GR	BR			GR				Base 0.50, 8.75, 10.50			
TRN 1.50 DS (76mm) UC/HC DIA.76	GR	BR							GR	BR			GR				Base 2.00, 3.50, 4.25, 5.25, 6.25, 7.25			
1.53 TRIVEX® Single Vision																				
TOG 1.53 Trivex (75mm) SV UC/HC DIA.70	GR	BR							GR	BR							Base 1.00			
TOG 1.53 Trivex (75mm) SV UC/HC DIA.75	GR	BR							GR	BR							Base 2.50, 4.00, 5.00, 6.00, 7.00, 8.00			
TRN 1.53 Trivex (78mm) SV UC/HC DIA.74													GR				Base 0.50, 1.00			
TRN 1.53 Trivex (78mm) SV UC/HC DIA.78													GR				Base 2.00, 3.50, 4.25, 5.25, 6.25, 7.25, 8.75			
1.59 Polycarbonate Single Vision																				
TRN 1.59 Polycarbonate (79mm) SV HC DIA.79	GR	BR							GR	BR			GR				Base 1.50, 3.25, 4.25, 4.75, 5.50, 7.00			
TRN 1.59 Polycarbonate (74mm) SV HC DIA.74	GR	BR							GR	BR							Base 7.50, 8.50			

1.60 High Index Single Vision										
TOG 1.60 SV (75mm) UC/HC DIA.65	GR	BR						GR	BR	Base 10.00, 11.00
TOG 1.60 SV (75mm) UC/HC DIA.70	GR	BR						GR	BR	Base 9.00
TOG 1.60 SV (75mm) UC/HC DIA.75	GR	BR						GR	BR	Base 1.00, 2.50, 4.00, 5.00, 6.00, 7.00, 8.00
TOG 1.60 SV (80mm) UC/HC DIA.80								GR	BR	Base 1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.50
TRN 1.60 SV (80mm) UC/HC DIA.76									GR	Base 8.50
TRN 1.60 SV (80mm) UC/HC DIA.80									GR	Base 1.75, 2.75, 3.75, 4.75, 6.00, 6.75
1.67 Super High Index Single Vision										
TOG 1.67 SV (65mm) UC/HC DIA.65	GR	BR						GR	BR	Base 10.00, 11.00
TOG 1.67 SV (70mm) UC/HC DIA.70	GR	BR						GR	BR	Base 9.00
TOG 1.67 SV (75mm) UC/HC DIA.75	GR	BR						GR	BR	Base 1.00, 2.50, 4.00, 5.00, 6.00, 7.00, 8.00
TOG 1.67 SV (80mm) UC/HC DIA.80								GR	BR	Base 1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.50
TRN 1.67 SV (80mm) UC/HC DIA.80									GR	Base 1.75, 2.75, 3.75, 4.75, 6.00
TRN 1.67 SV (80mm) UC/HC DIA.81									GR	Base 6.75, 8.00
1.74 ultra High Index Single Vision										
TOG 1.74 LeafECO SV (75mm) UC/HC DIA.65	GR	BR								Base 11.50, 12.50, 13.50, 14.50, 15.50, 16.50
TOG 1.74 LeafECO SV (75mm) UC/HC DIA.70	GR	BR								Base 0.50, 1.50, 2.25, 6.75, 8.50, 10.25
TOG 1.74 LeafECO SV (75mm) UC/HC DIA.75	GR	BR								Base 3.50, 5.00

											
Semi-Finished Single Vision										Sphere	Cyl to
1.50 Standard Index											
Unity SVx	GR	BR						GR	GR	-12.25 to +6.00	-6.00
Unity SVxtra	GR	BR						GR	GR	-12.25 to +6.00	-6.00
Unity SVxtreme	GR	BR						GR	GR	-6.00 to +6.00	-6.00
1.53 TRIVEX® Single Vision											
Unity SVx	GR	BR						GR	GR	-14.00 to +6.00	-6.00
Unity SVxtra	GR	BR						GR	GR	-14.00 to +6.00	-6.00
Unity SVxtreme	GR	BR						GR	GR	-6.50 to +6.00	-6.00
1.59 Polycarbonate Single Vision											
Unity SVx	GR	BR						GR	GR	-14.00 to +6.00	-6.00
Unity SVxtra	GR	BR						GR	GR	-14.00 to +6.00	-6.00
Unity SVxtreme	GR	BR						GR	GR	-6.50 to +6.00	-6.00
1.60 High Index Single Vision											
Unity SVx	GR	BR								-14.00 to +6.00	-6.00
Unity SVxtra	GR	BR								-14.00 to +6.00	-6.00
Unity SVxtreme	GR	BR								-6.50 to +6.00	-6.00
1.67 Super High Index Single Vision											
Unity SVx	GR	BR						GR	GR	-16.00 to +8.00	-6.00
Unity SVxtra	GR	BR						GR	GR	-16.00 to +8.00	-6.00
Unity SVxtreme	GR	BR						GR	GR	-6.00 to +6.00	-6.00

												
Semi-Finished Single Vision										Sphere	Cyl to	
1.50 Standard Index SFSV												
Younger Optics SFSV (76mm)	GR	BR	GG	SA	AME	AMB	E	R	GR	BR	-9.00 to +6.00	-4.00
Younger Optics SFSV DIGITAL										GR	-10.00 to +8.50	-4.00
1.53 TRIVEX® Single Vision												
Younger Optics Trilogy® Spherical SV	GR	BR							GR	BR	-8.00 to +7.50	-4.00
1.59 Polycarbonate Single Vision												

Younger Optics SV	GR	BR	GG	SA	AME	AMB	E	R	GR	BR			GR	-9.00 to +6.00	-4.00
1.60 High Index Single Vision															
Younger Optics Single Vision 1.60	GR	BR	GG	SA	AME	AMB	E	R						-11.00 to +8.00	-4.00
1.67 Super High Index Single Vision															
Younger Optics Single Vision 1.67 MR-10™	GR	BR	GG	SA	AME	AMB	E	R	GR	BR			GR	-11.00 to +7.00	-4.00
1.74 Ultra High Index Single Vision															
Younger Optics Single Vision 1.74	GR	BR	GG	SA	AME	AMB	E	R	GR	BR				-18.00 to +12.00	-4.00

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

Unless indicated otherwise, all registered trademarks and trademarks are the property of Essilor International and/or its subsidiaries in the United States and in other countries.

All other trademarks are the property of their respective owners.