










# Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Transitions® Signature® GEN 8™](#) / [Shamir Insights Inc.](#)

## Shamir Insights Inc. – Transitions® Signature® GEN 8™

												
Progressives									<b>Sphere</b>	<b>Add</b>	<b>Cyl to</b>	<b>MFH</b>
<b>1.50 Standard Index</b>												
Shamir Autograph® III	GR BR								-12.25 to +6.00	+0.75 to +4.00	-6.00	11, 13, 15, 18
Shamir Autograph® II+	GR BR								-12.25 to +6.00	+0.75 to +4.00	-6.00	11, 13, 15, 18
Shamir Intouch™	GR BR								-12.25 to +6.00	+0.75 to +3.50	-6.00	14, 16, 18
Shamir Spectrum+™	GR BR								-12.25 to +6.00	+0.75 to +3.50	-6.00	14, 16, 18
Shamir Element™	GR BR								-7.00 to +4.00	+0.75 to +3.50	-6.00	16, 19
Shamir GENESIS HD	GR BR								-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir Computer™	GR BR								-7.00 to +4.00	+0.75 to +2.00	-6.00	14 Lower, 10 Upper
Shamir Workspace™	GR BR								-6.00 to +6.00	+0.75 to +2.50	-6.00	14 Lower, 10 Upper
Shamir Attitude III Sport	GR BR								-6.50 to +5.25	+0.75 to +3.50	-6.00	18
Shamir Attitude III Fashion	GR BR								-6.50 to +5.25	+0.75 to +3.50	-6.00	15, 18
<b>1.53 TRIVEX®</b>												
Shamir Autograph Intelligence™	GR BR								-13.25 to +7.00	+0.75 to +4.00	-6.5	11, 13, 15, 18
Shamir Autograph® III	GR BR								-13.25 to +7.00	+0.75 to +4.00	-6.50	11, 13, 15, 18
Shamir Autograph® II+	GR BR								-13.25 to +6.00	+0.75 to +4.00	-6.50	11, 13, 15, 18
Shamir Intouch™	GR BR								-13.25 to +6.00	+0.75 to +3.50	-6.50	14, 16, 18
Shamir Golf™	GR BR								-5.00 to +4.00	+1.00 to +3.50	-6.50	19
Shamir Spectrum+™	GR BR								-13.25 to +6.00	+0.75 to +3.50	-6.50	14, 16, 18
Shamir Element™	GR BR								-13.25 to +6.00	+0.75 to +3.50	-6.50	16, 19
Shamir GENESIS HD	GR BR								-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir Computer™	GR BR								-6.00 to +8.00	+0.75 to +2.00	-6.50	14 Lower, 10 Upper
Shamir Workspace™	GR BR								-6.00 to +8.00	+0.75 to +2.50	-6.50	14 Lower, 10 Upper
Shamir Attitude III Sport	GR BR								-6.50 to +6.00	+0.75 to +3.50	-6.50	18
Shamir Attitude III Fashion	GR BR								-6.50 to +6.00	+0.75 to +3.50	-6.50	15, 18
<b>1.59 Polycarbonate</b>												
Shamir Autograph Intelligence™	GR BR								-14.75 to +7.50	+0.75 to +4.00	-7	11, 13, 15, 18
Shamir Autograph® III	GR BR								-14.75 to +8.00	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® II+	GR BR								-14.75 to +6.50	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Intouch™	GR BR								-14.75 to +6.00	+0.75 to +3.50	-7.00	15, 18
Shamir Golf™	GR BR								-6.00 to +5.25	+1.00 to +3.50	-7.00	19
Shamir Spectrum+™	GR BR								-14.75 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
Shamir Element™	GR BR								-10.00 to +6.00	+0.75 to +3.50	-7.00	16, 19
Shamir GENESIS HD	GR BR								-10.00 to +6.00	+0.75 to +3.50	-7.00	16
Shamir Computer™	GR BR								-6.00 to +8.00	+0.75 to +2.00	-7.00	14 Lower, 10 Upper
Shamir Workspace™	GR BR								-6.00 to +8.00	+0.75 to +2.50	-7.00	14 Lower, 10 Upper
Shamir Attitude III Sport	GR BR								-6.50 to +6.00	+0.75 to +3.50	-6.50	18
Shamir Attitude III Fashion	GR BR								-6.50 to +6.00	+0.75 to +3.50	-6.50	15, 18

**TRANSITIONS® LIGHT INTELLIGENT LENSES™**



*Transitions*, *Transitions Signature*, and *XTRActive* are registered trademarks, and *XTRActive Polarized*, the *Transitions* logo and *Transitions Light Intelligent Lenses* are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.