

Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Sun Lenses/Sun Shields](#)



OUTDOOR LIGHT CHANGES, YOUR SUNWEAR SHOULD TOO

The *Transitions®* line of light intelligent sun lenses and helmet shields automatically lighten and darken, from sunrise to sunset, for an all-day, visually comfortable experience. Our light intelligent sunwear technologies make sports and outdoor activities more convenient by eliminating the need to switch between different pairs of sunwear throughout the day. They help reduce glare, increase contrast, and improve distance and depth perception.



[LEARN MORE](#)
▼




[LEARN MORE](#)
▼



[LEARN MORE](#)
▼



Transitions® Drivewear® Sun Lenses



OLIVE GREEN
To remove glare that would otherwise impact in low light conditions



COPPER
To remove glare for safe driving vision



BROWN
To provide maximum comfort in high light conditions



TRANSITIONS® DRIVEWEAR® ADAPTIVE SUNGLASSES

COLOR: Polarized, Olive Green to Copper to Dark Red-Brown	
FSV	<i>Transitions®</i> Drivewear® Finished Plano HC/Backside AR
Materials: TRIVEX®, Trilogy®	Sphere: Plano to -6.00
Add: N/A	

TRANSITIONS® DRIVEWEAR® ADAPTIVE SUN LENSES

COLOR: Polarized, Olive Green to Copper to Dark Red-Brown		
SFSV	<i>Transitions® Drivewear® SFSV HC</i>	
Material: 1.50, Polycarbonate, Trilogy®	Sphere: -9.00D to +7.50D	
Add: N/A	Cyl to: -4.00	
PROGRESSIVES		
<i>Transitions® Drivewear® Image®</i>		
Material: Polycarbonate	Sphere: -9.00 to +7.00	
Add: +1.00 to +3.00	Cyl to: -4.00	MFH: 18
DIGITAL – Younger Optics		
<i>Camber™</i>		
Material: Polycarbonate	Sphere: -10.00 to +6.25	Cyl to: -4.00



Light Intelligent Sunglasses

Oakley® offers both prescription and plano Oakley® with Transitions® sun lenses and sunglasses, designed to activate to darker, neutral tints to allow athletes to maintain peak performance. Rx sun lenses available through the Oakley lab. Visit Oakley.com or call 800-403-7449 for more info.



OAKLEY® with TRANSITIONS®

Plano	
Color: Clear to Black Iridium®	
Materials: Oakley Plutonite®	Sphere: Plano
Oakley® Flak® 2.0 XL	Oakley® EVZero® Range®
Oakley® Jawbreaker™	Oakley® Radarlock® Edge®
Oakley® Sutro	Oakley® Sutro Lite
Oakley® Radar® EV Advancer	Oakley® Flight Jacket™
Oakley® Radar® Path®	Oakley® Radarlock® Pitch®
Oakley® Radar® EV Path®	Oakley® EVZero® Blades™
Oakley® Radar® EV Pitch®	Oakley® Radar® EV XS
Oakley® Sphaera™ with Transitions®	

OAKLEY® with TRANSITIONS® SUN LENSES

SFSV	
Color: Clear to Black Iridium®	Oakley® with <i>Transitions® SV lenses</i>
Materials: Polycarbonate	Sphere: -3.25 to +2.00
Add: N/A	Cyl to: -3.00
PROGRESSIVES	
Color: Clear to Black Iridium®	Oakley® with <i>Transitions® PAL lenses</i>

Materials: Polycarbonate	Sphere: -3.25 to +2.00
Add: +1.00 to +3.00	Cyl to: -3.00

Light Intelligent Ballistic Grade Sunglasses

The ESS Crossbow® Photochromic™ eyeshield combines ballistic standards with Transitions Optical's light adapting technology to redefine performance and utility for eye protection. ClearZone™ FlowCoat technology minimizes scratches and Tri-Tech™ Fit enhances stability and comfort for any outdoor activity.



ESS Crossbow® Photochromic™ Eyeshield

Color: Clear to Dark Gray	ESS Crossbow One Ballistics Grade Eyeshield
Sphere: Plano	



Light Intelligent Shields

Motorcycle helmet face shields help you take your riding experience to the next level. Unlike clear or tinted shields, these shields automatically adjust their level of darkness to optimize vision in sun and shade, day and night.



TRANSITIONS® MOTORCYCLE FACE SHIELDS

Color: Clear to Dark Gray	Style: Motorcycle Helmet	
SHOEI® Helmets: SHOEI® Transitions® CWR-1 face	KLIM® Helmets: KLIM® Transitions® face shields	Nexx® Transitions® face shield

<p>shield SHOEI® Transitions CWR-F2 face shield</p>		
<p>Availability: Available in CWR-1 Pinlock® face shield for RF-SR, and the CWR-F2 Pinlock® face shield for the RF-1400 and X-Fifteen helmet through SHOEI helmet dealers and Revzilla.com.</p>	<p>Availability: Comes standard with the Krios Pro Helmet and TK1200 Karbon Modular Helmet. Sold as an accessory for the Krios Karbon helmet, Krios Pro Helmet and the TK1200 Karbon Modular Helmet. Available at Klim.com or your local Klim dealer.</p>	<p>Availability: Available for XR3 helmets. Sold as an accessory. Available at your local Nexx dealer.</p>



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.