










Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® XTRActive® Polarized™](#) / [Shamir Insights Inc.](#)

Shamir Insights, Inc. – Transitions® XTRActive® Polarized™

|  | |  | | | |
|---|--|---|----------------|---------------|--|
| Finished Single Vision | |  | Sphere | Cyl to | |
| 1.59 FSV | | | | | |
| Shamir FSV with Glacier PLUS™ AR | | | -6.00 to +4.00 | -2.00 | |

|  | |  | | | |
|---|--|---|------------------|---------------|--|
| Semi-Finished Single Vision | |  | Sphere | Cyl to | |
| 1.50 SFSV | | | | | |
| Shamir Autograph® III Single Vision | | GR | -12.00 to +7.25 | -6.00 | |
| Shamir Relax™ | | GR | -14.00 to +8.00 | -6.00 | |
| 1.53 SFSV | | | | | |
| Shamir Autograph® III Single Vision | | GR | -14.00 to +8.00 | -6.50 | |
| Shamir Relax™ | | GR | -14.00 to +8.00 | -6.50 | |
| 1.59 SFSV | | | | | |
| Shamir Autograph® III Single Vision | | GR | -14.00 to +8.75 | -7.00 | |
| Shamir Relax™ | | GR | -14.75 to +6.00 | -7.00 | |
| 1.67 SFSV | | | | | |
| Shamir Autograph® III Single Vision | | GR | -16.00 to +16.00 | -8.00 | |
| Shamir Relax™ | | GR | -16.75 to +16.00 | -8.00 | |

|  | |  | | | | | |
|---|--|---|-----------------|----------------|---------------|--------------------|--|
| Progressives | |  | Sphere | Add | Cyl to | MFH | |
| 1.50 Standard Index | | | | | | | |
| Shamir Autograph Intelligence™ | | GR | -12.25 to +6.50 | +0.75 to +4.00 | -6.00 | 11, 13,15,18 | |
| Shamir Autograph® III | | GR | -12.25 to +6.00 | +0.75 to +4.00 | -6.00 | 11, 13, 15, 18 | |
| Shamir Autograph® II+ | | GR | -12.25 to +6.00 | +0.75 to +4.00 | -6.00 | 11, 13, 15, 18 | |
| Shamir Intouch™ | | GR | -12.25 to +6.00 | +0.75 to +3.50 | -6.00 | 14, 16, 18 | |
| Shamir Spectrum+™ | | GR | -12.25 to +6.00 | +0.75 to +3.50 | -6.00 | 14, 16, 18 | |
| Shamir Element™ | | GR | -7.00 to +4.00 | +0.75 to +3.50 | -6.00 | 16, 19 | |
| Shamir GENESIS HD | | GR | -6.00 to +6.00 | +0.75 to +3.50 | -4.00 | 16 | |
| Shamir Computer™ | | GR | -7.00 to +4.00 | +0.75 to +2.00 | -6.00 | 14 Lower, 10 Upper | |
| Shamir Workspace™ | | GR | -6.00 to +6.00 | +0.75 to +2.50 | -6.00 | 14 Lower, 10 Upper | |
| Shamir Attitude III Sport | | GR | -6.50 to +5.25 | +0.75 to +3.50 | -6.00 | 18 | |
| Shamir Attitude III Fashion | | GR | -6.50 to +5.25 | +0.75 to +3.50 | -6.00 | 15, 18 | |
| 1.53 TRIVEX® | | | | | | | |
| Shamir Autograph Intelligence™ | | GR | -13.25 to +7.00 | +0.75 to +4.00 | -6.50 | 11, 13, 15, 18 | |
| Shamir Autograph® III | | GR | -13.25 to +7.00 | +0.75 to +4.00 | -6.50 | 11, 13, 15, 18 | |
| Shamir Autograph® II+ | | GR | -13.25 to +6.00 | +0.75 to +4.00 | -6.50 | 11, 13, 15, 18 | |
| Shamir Intouch™ | | GR | -13.25 to +6.00 | +0.75 to +3.50 | -6.50 | 14, 16, 18 | |
| Shamir Golf™ | | GR | -5.00 to +4.00 | +1.00 to +3.50 | -6.50 | 19 | |
| Shamir Spectrum+™ | | GR | -13.25 to +6.00 | +0.75 to +3.50 | -6.50 | 14, 16, 18 | |
| Shamir Element™ | | GR | -13.25 to +6.00 | +0.75 to +3.50 | -6.50 | 16, 19 | |
| Shamir GENESIS HD | | GR | -6.00 to +6.00 | +0.75 to +3.50 | -4.00 | 16 | |

| | | | | | |
|---------------------------------|----|------------------|----------------|-------|--------------------|
| Shamir Computer™ | GR | -6.00 to +8.00 | +0.75 to +2.00 | -6.50 | 14 Lower, 10 Upper |
| Shamir Workspace™ | GR | -6.00 to +8.00 | +0.75 to +2.50 | -6.50 | 14 Lower, 10 Upper |
| Shamir Attitude III Sport | GR | -6.50 to +6.00 | +0.75 to +3.50 | -6.50 | 18 |
| Shamir Attitude III Fashion | GR | -6.50 to +6.00 | +0.75 to +3.50 | -6.50 | 15, 18 |
| 1.59 Polycarbonate | | | | | |
| Shamir Autograph Intelligence™ | GR | -14.75 to +7.50 | +0.75 to +4.00 | -7.00 | 11, 13, 15, 18 |
| Shamir Autograph® III | GR | -14.75 to +8.00 | +0.75 to +4.00 | -7.00 | 11, 13, 15, 18 |
| Shamir Autograph® II+ | GR | -14.75 to +6.50 | +0.75 to +4.00 | -7.00 | 11, 13, 15, 18 |
| Shamir Intouch™ | GR | -14.75 to +6.00 | +0.75 to +3.50 | -7.00 | 15, 18 |
| Shamir Golf™ | GR | -6.00 to +5.25 | +1.00 to +3.50 | -7.00 | 19 |
| Shamir Spectrum+™ | GR | -14.75 to +6.00 | +0.75 to +3.50 | -7.00 | 14, 16, 18 |
| Shamir Element™ | GR | -10.00 to +6.00 | +0.75 to +3.50 | -7.00 | 16, 19 |
| Shamir GENESIS HD | GR | -10.00 to +6.00 | +0.75 to +3.50 | -7.00 | 16 |
| Shamir Computer™ | GR | -6.00 to +8.00 | +0.75 to +2.00 | -7.00 | 14 Lower, 10 Upper |
| Shamir Workspace™ | GR | -6.00 to +8.00 | +0.75 to +2.50 | -7.00 | 14 Lower, 10 Upper |
| Shamir Attitude III Sport | GR | -6.50 to +6.00 | +0.75 to +3.50 | -6.50 | 18 |
| Shamir Attitude III Fashion | GR | -6.50 to +6.00 | +0.75 to +3.50 | -6.50 | 15, 18 |
| 1.67 Super High Index | | | | | |
| Shamir Autograph Intelligence™ | GR | -16.75 to +10.00 | +0.75 to +4.00 | -8.00 | 11, 13, 15, 18 |
| Shamir Autograph® III | GR | -16.75 to +9.00 | +0.75 to +4.00 | -8.00 | 11, 13, 15, 18 |
| Shamir Autograph® II+ | GR | -16.75 to +8.25 | +0.75 to +4.00 | -8.00 | 11, 13, 15, 18 |
| Shamir Spectrum+™ | GR | -16.75 to +8.00 | +0.75 to +4.00 | -8.00 | 11, 13, 15, 18 |
| Shamir Element™ SuperLite™ 1.67 | GR | -12.50 to +8.00 | +0.75 to +3.50 | -8.00 | 16, 19 |
| Shamir GENESIS HD | GR | -6.00 to +6.00 | +0.75 to +3.50 | -4.00 | 16 |
| Shamir InTouch™ | GR | -16.75 to +8.00 | +0.75 to +3.50 | -8.00 | 15, 18 |
| Shamir Attitude III Sport | GR | -6.50 to +6.00 | +0.75 to +3.50 | -8.00 | 18 |
| Shamir Attitude III Fashion | GR | -6.50 to +6.00 | +0.75 to +3.50 | -8.00 | 15, 18 |

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.